

Integra[®] MicroFrance[®]

Laparoscopic Instruments for Urology Procedures*

Designed with world-renowned surgeons, MicroFrance laparoscopic instruments are known for quality and precision. These instruments are designed to perform a variety of laparoscopic procedures.



*Instruments under the medical grouping can be used in other general surgical procedures.

Urology Procedures

Instrument Trays

Reference	Description
3717001	3-Level Aluminum Tray with Plastic Lid; 53cm L x 25.5cm W x 15cm D
7874MF-61	Stainless Steel Tray; Tray, 648 x 187 x 95mm




Modular Instruments

5mm Diameter

Bipolar Forceps, 310mm Working Length

Reference	Description
CEV634-1A	Jaw Insert; Mouiel Fenestrated Bipolar Forceps
CEV6795B	Non-Insulated Tube
CEV669B	Bipolar Angled Handle, Non-Insulated Tube



Monopolar Hook

Reference	Description
CEV2295A	Jaw Insert; Hook Jaws, Size 5, 3-Pack
CEV229-1A	Handle; Round Hook Handle, 350mm Working Length



Modular Instruments

5mm Diameter


Monopolar Forceps, 310mm Working Length – Jaw Inserts

Reference	Description
CEV642-1	Jaw Insert; Delaitre Dissector, Curved, 20mm
CEV532I	Jaw Insert; Short Dissector, Curved, 12mm
CEV625-1	Jaw Insert; Johann Fenestrated Forceps, 20 x 5mm Jaws
CEV611-1	Jaw Insert; Manhès Grasping Forceps
CEV615-1	Jaw Insert; Manhès Alligator Forceps, 15mm Jaws
CEV625E	Jaw Insert; Fenestrated Forceps, 20 x 5mm
CEV605G	Jaw Insert; Insulated Scissors, Very Curved
CEV605GT	Jaw Insert; Insulated Scissors, Very Curved, Titanium
CEV607	Jaw Insert; Insulated Scissors Tip, Very Curved, 6-pack*
CEV607-1	Wire Insert*



Monopolar Forceps, 310mm Working Length – Tubes

Reference	Description
CEV649-5B	Insulated Black
CEV646B5	Insulated Yellow
CEV647B5	Insulated Blue
CEV648B5	Insulated Green
CEV649B5	Non-Insulated*
CEV649C5	Disposable Sheath, Black, 3-pack*



Monopolar Forceps, 310mm Working Length – Handles

Reference	Description
CEV10395R	Rotatable Monopolar Angled Handle
CEV10395D	Rotatable Monopolar Angled Handle, Non-Locking
CEV10195C	Rotatable Monopolar Angled Handle with Internal Ratchet




*Jaw inserts, tubes, and handles in this section are interchangeable. Insulated scissors tip CEV607 must be used with wire insert CEV607-1. Non-insulated tube CEV649B5 must be used with disposable sheath CEV649C5.

Modular Instruments

10mm Diameter

Begin Clamp Applicator, 310mm Working Length

Reference	Description
CEV560-1	Jaw Insert; Begin Clamp Applicator
CEV649-10B	Non-Insulated Tube; Metallic
CEV658-10	Handle; Rotatable Axial Handle



Begin Clamp

Reference	Description
CEV565	Begin Clamp, 30 mm Jaws, 2 Pack
CEV566	Begin Clamp, 40 mm Jaws, 2 Pack
CEV567	Begin Clamp, 50 mm Jaws, 2 Pack
CEV568	Sartori Vascular Clamp, 15mm Jaws, 2 Pack
CEV560	Begin Clamp Applicator, 10mm Diameter, 360mm Working Length



Forceps, 310mm Working Length

Reference	Description
CEV606-1	Jaw insert; Museux Grasping Forceps
CEV649-10B	Non-Insulated tube; Metallic
CEV102910R	Angled Handle with Ratchet for Rotatable Non-Insulated 5mm and 10mm Tubes



Standard Instruments


Clamp, 385mm Working Length, 10mm Diameter

Reference	Description
CEV535	Querleu Vascular Clamp with Axial Ratcheted Handle, 45mm Jaws



Monopolar Forceps, 365mm Working Length, 5mm Diameter

Reference	Description
CEV405	Johann Fenestrated Forceps with Angled Handle and Ratchet, 20mm Jaws



Needle Holder

Reference	Description
CEV700T5	Needle Holder with Carbide Jaws and Ratchet Up, 335mm Length, 5mm Tips



Standard Instruments

10mm Diameter

Retractor

Reference	Description
CEV453	Zimmermann Retractor with 5 Blades, 60mm Jaws, 510mm Length



Trocars

5mm Diameter

Valve Heads

Reference	Description
CEV177A5	Valve Head, Magnetic
CEV127T	Valve Head, Silicone



CEV177A5



CEV127T

Tips

Reference	Description
CEV145A3	Sharp Tip, 134mm Working Length
CEV8145A3	Sharp Tip, 105mm Working Length



CEV145A3



CEV8145A3

Cannula

Reference	Description
CEV145B6	Threaded Cannula, 80mm Long
CEV8145B6	Threaded Cannula, 50mm Long
CEV145B7	Smooth Cannula, 80mm Long
CEV8145B7	Smooth Cannula, 50mm Long



CEV145B6



CEV8145B6



CEV145B7



CEV8145B7

Cables

Monopolar

Reference	Description	Qty
CP391E	For Erbe Generators, 3.5m	1
CP391M	For Martin – Lamidey Generators, 3.5m	1
CP391V	For Valleylab Generators, 3.5m	1

Bipolar

Reference	Description	Qty
CP393-3	For Erbe, Storz, and Wolf Generators, 4.5m	1
CP393-4	For Martin, Berchtold, and Aesculap Generators, 4.5m	1
CP393-2	For Valleylab and Lamidey Generators, 4.5m	1

Miscellaneous

Trocar Accessories


Reference	Description	Qty
CEV179A5	Replacement Valves for Magnetic Head CEV177A5	6
CEV125P6	White Seal for Trocar Valve, 3.5mm, 5mm Diameter	6
CEV105P6	Sealing Caps for 5mm Diameter Trocar Heads, Blue	6
CEV60-120	Insufflation Needle with Luer-Lock, 120mm Length	1

Accessories

5.5mm Diameter


Suction Irrigation, 300mm Working Length, 5mm Diameter

Reference	Description
CEV103-5	Suction Irrigation Tube with Handle
CEV103AH	Handle



Flexible Tension Arm

Reference	Description
CEV984	Operating Table Clamp
CEV980A	Mosnier Flexible Tension Arm, 500mm
CEV980B	Mosnier Flexible Tension Arm, 600mm
CEV980-3	Mosnier Instrument Holder for Arm



Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.
- Consult product labels and inserts for any indication, contraindications, hazards, warnings, precautions, and instructions for use.

For more information or to place an order, please contact:

United States, Canada, Asia, Pacific, Latin America

USA 800-654-2873 ▪ 888-980-7742 fax

International +1 609-936-5400 ▪ +1 609-750-4259 fax

integralife.com/contact

