


# Integra<sup>®</sup> MicroFrance<sup>®</sup> Abitbol Laryngeal Instruments

Developed in conjunction with Jean Abitbol, MD

These refined instruments may help reduce many of the complications associated with laser surgery of the larynx and will add utility to any laryngeal instrument set. The set of grasping forceps is insulated and may help minimize the possibilities of laser-related incidents, has integrated suction for evacuation of laser plume and monopolar jaws to allow hemostasis in the areas of the larynx that are difficult to address.



## Instruments for laser laryngeal surgery

Reference	Description
MCL205	Monopolar Grasping Forceps with Suction, Fenestrated 3mm x 8mm Jaws Overall Length 328mm (12.9") Working Length 200mm (7.9")
MCL219	Monopolar Heart-Shaped Graspers (2mm x 3mm V-Shaped Tip) with Suction, Angled Right Overall Length 320mm (12.6") Working Length 190mm (7.48")
MCL219-1	Monopolar Heart-Shaped Graspers (2mm x 3mm V-Shaped Tip) with Suction, Angled Left Overall Length 320mm (12.6") Working Length 195mm (7.7")
MCL241	Monopolar Alligator Forceps (1mm x 8mm Tip) with Suction, Curved Right Overall Length 328mm (12.9") Working Length 200mm (7.9")
MCL241-1	Monopolar Alligator Forceps (1mm x 8mm Tip) with Suction, Curved Left Overall Length 328mm (12.9") Working Length 200mm (7.9")
MCL210	Suction Tube with Finger Control, Insulated Overall Length 298mm (11.7") Working Length 226mm (8.9")
37-17002	One Level Instrument Tray (not pictured) 25cm wide x 53cm long x 4cm deep
CP391	Cable for Monopolar Coagulating Instruments (not pictured) 185cm




MCL205



MCL219/219-1



MCL241/241-1



MCL210

The MicroFrance® brand is a premier line of ENT instruments, designed with world-renowned ENT surgeons to facilitate new techniques and improve patient outcomes.

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Always refer to the appropriate instructions for use for complete clinical instructions.
- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.
- Please consult product labels and inserts for any indication, contraindications, hazards, warnings, precautions and instructions for use.

**For more information or to place an order, please contact:**

**United States, Canada, Asia, Pacific, Latin America**

USA 800-654-2873 ▪ 888-980-7742 fax

International +1 609-936-5400 ▪ +1 609-750-4259 fax

[integralife.com/contact](http://integralife.com/contact)



Integra MicroFrance SAS  
Le Pavillon, 03160  
Saint-Aubin-le-Monial ▪ France



Integra York PA, Inc.  
589 Davies Drive  
York, PA 17402 ▪ USA