

Integra®

Neuro Spine Kerrison Rongeurs

Limit uncertainty with
durability and quality



INTEGRA®
LIMIT UNCERTAINTY

Integra® Ruggles®-Redmond™ offers a wide selection of kerrison rongeurs. With multiple pattern and coating variations, Integra's Neuro Spine Instruments are sure to meet the needs of the surgeon.

Surface coating options

The Ultra-Coat™ Kerrison Rongeurs have an extremely hard material coating that increases the surface hardness vs. regular stainless steel kerrison rongeurs, which extends the life of the cutting edge.



The Diamond™ Kerrison Rongeurs have a black ceramic coating that increases the surface hardness vs. regular stainless steel kerrison rongeurs, extending the life of the cutting edge and reducing intraoperative glare.



The Stealth™ Kerrison Rongeurs have a surface treatment intended to reduce friction between metal surfaces. The Stealth coating does not increase the hardness of the cutting edge, but reduces intraoperative glare.



All Ultra-Coat™ and Diamond™ Kerrison Rongeurs are covered by a sharpness warranty for a period of two (2) years under normal use. Should sharpening become necessary during the warranty period, Ultra-Coat and Diamond Kerrisons will be resharpened at no charge.

Detach Coated Rongeurs

The instrument disassembles to aid in the cleaning and visual inspection.



Extended Length Kerrison Rongeurs

Additional lengths to accommodate increasing patient BMI's.
Available in Ultra-Coat™, Diamond™ and Stealth™ coatings.



Ultra-Coat Detach Kerrison Rongeurs - 9mm jaw opening

Catalog Number	Description
RB5760DT	Ultra-Coat, 7", 40 degree up, 1mm
RB5762DT	Ultra-Coat, 7", 40 degree up, 2mm
RB5763DT	Ultra-Coat, 7", 40 degree up, 3mm
RB5764DT	Ultra-Coat, 7", 40 degree up, 4mm
RB5765DT	Ultra-Coat, 7", 40 degree up, 5mm
RB5766DT	Ultra-Coat, 7", 40 degree up, 6mm
RB5860DT	Ultra-Coat, 8", 40 degree up, 1mm
RB5862DT	Ultra-Coat, 8", 40 degree up, 2mm
RB5863DT	Ultra-Coat, 8", 40 degree up, 3mm
RB5864DT	Ultra-Coat, 8", 40 degree up, 4mm
RB5865DT	Ultra-Coat, 8", 40 degree up, 5mm
RB5866DT	Ultra-Coat, 8", 40 degree up, 6mm
RB5960DT	Ultra-Coat, 9", 40 degree up, 1mm
RB5962DT	Ultra-Coat, 9", 40 degree up, 2mm
RB5963DT	Ultra-Coat, 9", 40 degree up, 3mm
RB5964DT	Ultra-Coat, 9", 40 degree up, 4mm
RB5965DT	Ultra-Coat, 9", 40 degree up, 5mm
RB5966DT	Ultra-Coat, 9", 40 degree up, 6mm
RB5160DT	Ultra-Coat, 11", 40 degree up, 1mm
RB5162DT	Ultra-Coat, 11", 40 degree up, 2mm
RB5163DT	Ultra-Coat, 11", 40 degree up, 3mm
RB5164DT	Ultra-Coat, 11", 40 degree up, 4mm
RB5165DT	Ultra-Coat, 11", 40 degree up, 5mm
RB5166DT	Ultra-Coat, 11", 40 degree up, 6mm

Diamond Detach Kerrison Rongeurs - 9mm jaw opening

Catalog Number	Description
RD4771DT	Diamond, 7", 40 degree up, 1mm
RD4752DT	Diamond, 7", 40 degree up, 2mm
RD4753DT	Diamond, 7", 40 degree up, 3mm
RD4754DT	Diamond, 7", 40 degree up, 4mm
RD4755DT	Diamond, 7", 40 degree up, 5mm
RD4756DT	Diamond, 7", 40 degree up, 6mm
RD4891DT	Diamond, 8", 40 degree up, 1mm
RD4882DT	Diamond, 8", 40 degree up, 2mm
RD4853DT	Diamond, 8", 40 degree up, 3mm
RD4854DT	Diamond, 8", 40 degree up, 4mm
RD4855DT	Diamond, 8", 40 degree up, 5mm
RD4856DT	Diamond, 8", 40 degree up, 6mm
RD4991DT	Diamond, 9", 40 degree up, 1mm
RD4982DT	Diamond, 9", 40 degree up, 2mm
RD4953DT	Diamond, 9", 40 degree up, 3mm
RD4954DT	Diamond, 9", 40 degree up, 4mm
RD4955DT	Diamond, 9", 40 degree up, 5mm
RD4956DT	Diamond, 9", 40 degree up, 6mm
RD41191DT	Diamond, 11", 40 degree up, 1mm
RD41192DT	Diamond, 11", 40 degree up, 2mm
RD41193DT	Diamond, 11", 40 degree up, 3mm
RD41194DT	Diamond, 11", 40 degree up, 4mm
RD41195DT	Diamond, 11", 40 degree up, 5mm
RD41196DT	Diamond, 11", 40 degree up, 6mm

For a complete listing of kerrison rongeurs refer to the catalog

For more information or to place an order, please contact:
 Integra ■ 311 Enterprise Drive, Plainsboro, NJ 08536
 USA and Canada: 877.468.5572 ■ 914.592.9050 (Outside USA)
 800.275.7192 (Fax)
integralife.com

Stealth Detach Kerrison Rongeurs - 9mm jaw opening

Catalog Number	Description
RB4891DT	Stealth, 8", 40 degree up, 1mm
RB4882DT	Stealth, 8", 40 degree up, 2mm
RB4853DT	Stealth, 8", 40 degree up, 3mm
RB4854DT	Stealth, 8", 40 degree up, 4mm
RB4855DT	Stealth, 8", 40 degree up, 5mm
RB4856DT	Stealth, 8", 40 degree up, 6mm

Detach Medium Kerrison Rongeurs - 15mm jaw opening

Catalog Number	Description
RB5860DTM	Ultra-Coat, 8", 40 degree up, 1mm
RB5862DTM	Ultra-Coat, 8", 40 degree up, 2mm
RB5863DTM	Ultra-Coat, 8", 40 degree up, 3mm
RB5864DTM	Ultra-Coat, 8", 40 degree up, 4mm
RB5865DTM	Ultra-Coat, 8", 40 degree up, 5mm
RB5866DTM	Ultra-Coat, 8", 40 degree up, 6mm
RD4821DTM	Diamond, 8", 40 degree up, 1mm
RD4822DTM	Diamond, 8", 40 degree up, 2mm
RD4823DTM	Diamond, 8", 40 degree up, 3mm
RD4824DTM	Diamond, 8", 40 degree up, 4mm
RD4825DTM	Diamond, 8", 40 degree up, 5mm
RD4826DTM	Diamond, 8", 40 degree up, 6mm
RB4821DTM	Stealth, 8", 40 degree up, 1mm
RB4822DTM	Stealth, 8", 40 degree up, 2mm
RB4823DTM	Stealth, 8", 40 degree up, 3mm
RB4824DTM	Stealth, 8", 40 degree up, 4mm
RB4825DTM	Stealth, 8", 40 degree up, 5mm
RB4826DTM	Stealth, 8", 40 degree up, 6mm

Extended Length Kerrison Rongeurs - 9mm jaw opening

Catalog Number	Description
RB5162	Ultra-Coat, 11", 40 degree up, 2mm
RB5163	Ultra-Coat, 11", 40 degree up, 3mm
RB5164	Ultra-Coat, 11", 40 degree up, 4mm
RB5165	Ultra-Coat, 11", 40 degree up, 5mm
RB5166	Ultra-Coat, 11", 40 degree up, 6mm
RD41192	Diamond, 11", 40 degree up, 2mm
RD41193	Diamond, 11", 40 degree up, 3mm
RD41194	Diamond, 11", 40 degree up, 4mm
RD41195	Diamond, 11", 40 degree up, 5mm
RD41196	Diamond, 11", 40 degree up, 6mm
RB5362	Ultra-Coat, 13", 40 degree up, 2mm
RB5363	Ultra-Coat, 13", 40 degree up, 3mm
RB5364	Ultra-Coat, 13", 40 degree up, 4mm
RB5365	Ultra-Coat, 13", 40 degree up, 5mm
RB5366	Ultra-Coat, 13", 40 degree up, 6mm

Extended Length Bayonet Kerrison Rongeurs - 9mm jaw opening

Catalog Number	Description
RB5662	Ultra-Coat, Bayonet 11", 40 degree up, 2mm
RB5663	Ultra-Coat, Bayonet 11", 40 degree up, 3mm
RB5664	Ultra-Coat, Bayonet 11", 40 degree up, 4mm
RB5665	Ultra-Coat, Bayonet 11", 40 degree up, 5mm
RB5666	Ultra-Coat, Bayonet 11", 40 degree up, 6mm