


MAXIMUM SPEED AND DURABILITY* FOR ENT SURGERY

Average Resecting Efficiency (g/min)




XPS burs operating at 12K RPM and 30K Burs at 30K RPM

Data on file

30K Burs

For the StraightShot™ M5
Microdebrider



Medtronic

MAXIMUM SPEED AND DURABILITY* FOR ENT SURGERY.

Finding the right bur for the job has never been easier. With the largest selection of burs for frontal sinus surgery, you can choose a bur for the specific anatomic needs of each patient or procedure.

30K Burs have multiple technology innovations, offering:

- **Speeds up to 30,000 rpm**
- **Improved durability and stability***
- **Up to 10 times more effective at removing material***
- **Best selection available for frontal sinus drilling**
- **Integrated irrigation**
- **Distal suction**

*Compared to Medtronic XPS™ Burs

Rx only. Refer to product instruction manual/package insert for instructions, warnings, precautions, and contraindications.

For further information, please call Medtronic ENT at 800.874.5797 or consult Medtronic's website at www.medtronicent.com or www.medtronic.com/30Kburs.

Product #	Description	Length	Qty
1884015RTD 	30K Tapered Diamond Choanal Atresia Bur, 4.0 mm, 15°	13 cm	1
1883040BLD 	30K Bullet Diamond Bur, 3.0 mm, 40°	13 cm	1
1883040BLC 	30K Bullet Cutting Bur, 3.0 mm, 40°	13 cm	1
1884040RTC 	30K Reverse Tapered Cutting Bur, 4.0 mm, 40°	13 cm	1
1883040BRC 	30K Frontal Finesse Barrel Cutting Bur, 3.0 mm, 40°	13 cm	1
1883655BRC 	30K Barrel Cutting Bur, 3.6 mm, 55°	13 cm	1
1884070RTD 	30K Reverse Tapered Diamond Bur, 4.0 mm, 70°	13 cm	1
1884070RTC 	30K Reverse Tapered Cutting Bur, 4.0 mm, 70°	13 cm	1
1883070BLD 	30K Bullet Diamond Bur, 3.0 mm, 70°	13 cm	1
1883070BLC 	30K Bullet Cutting Bur, 3.0 mm, 70°	13 cm	1

Medtronic

Medtronic plc ENT

6743 Southpoint Drive N
Jacksonville, FL 32216
USA
Toll free: (800) 874-5797
Telephone: (904) 296-9600
Fax: (800) 678-3995

International Telephone Numbers

Europe HQ Switzerland 41-21-802-7000
Latin America HQ 305-500-9328
Adriatic Region 385-1-488-1120
Argentina 54-11-4898-5700
Australia 1-800-668-670
Baltic Region 37-1-67560226
Belgium 32-2456-09-09
Brazil 55-11-2182-9200
Canada 1-800-268-5346
Chile 56-2-2655-5110
China 86-10-5869-8989
Colombia 57-1-742-7300
Czech Republic 420-233-059-111
France 33-470-679-800
Germany 49-2159-8149-353
Greece 30-210-67-79-099
Hong Kong 852-2919-1300
Hungary 36-188-90600

India 91-22-33074700
Israel 972-9-972-4400
Italy 39-02-24137-324
Japan 81-3-6774-4611
Korea 82-2-3404-3600
Lebanon 961-1-370-670
Luxembourg 32-2456-09-09
Malaysia 60-37-953-4800
Mexico 52-55-11-02-90-30
Netherlands 31-45-566-8800
Poland 48-22-4656900
Russian Fed. 7-495-580-73-77
Singapore 65-6436-5000
South Africa 27-11-260-9300
Spain 34-91-625-05-40
Taiwan 886-2-2183-6000
Thailand 662-232-7400
UK 44-1923-205-166

Part Numbers Defined

All start with 188 Diameter (4.0 mm) Angle (15°)

1884015RTD

BL = Bullet D = Diamond
RT = Reverse Taper C = Cutter or Fluted (synonymous)
BR = Barrel

© 2017 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. UC201500639d EN 03.2017