




Malleable SuctionTM Instruments

BENDABLE SUCTION WITH BUILT-IN NAVIGATION




Re-shape a suction instrument during image-guided surgery without having to re-verify it.

Malleable Suction™ Instruments

New Malleable Suction™ Instruments feature built-in AxiEM™ navigation technology within the instrument tip. Simply bend the instrument as necessary, and continue operating. There's no need to re-verify.

- **Plug-and-play with Medtronic ENT electro-magnetic (EM) navigation**
- **Multiple lengths, diameters, and tips for precise localization and access**
- **Standard suction design with built-in EM tracker**



Standard tip



Angle tip



Ball tip

Part Number	Description	Length
9735015	Small standard tip	14.5 cm
9735016	Medium standard tip	16.5 cm
9735017	Medium angle tip	16.5 cm
9735018	Medium ball tip	16.5 cm
9735019	Large standard tip	20 cm

For further information, please call Medtronic ENT at 800.874.5797 or 904.296.9600. You may also consult our website at www.MedtronicENT.com.

Medtronic ENT

6743 Southpoint Drive N
Jacksonville, FL 32216
USA
www.MedtronicENT.com
Toll free: (800) 874-5797
Fax: (800) 678-3995

International Telephone Numbers

Adriatic Region 385-1-488-1120	Italy 39-02-24137-324
Australia 1-800-668-670	Japan 81-3-6430-2017
Baltic Region 37-1-67560226	Korea 82-2-3404-3600
Belgium 32-2456-09-09	Lebanon 961-1-370-670
Canada 1-800-217-1617	Luxembourg 32-2456-09-09
China 86-21-50800998	Netherlands 31-45-566-8800
Czech Republic 420-2-9657-9580	Poland 48-22-4656900
France 33-470-679-800	Russian Fed. 7-495-580-73-77
Germany 49-2159-8149-209	Singapore 65-6776-6255
Greece 30-210-67-79-099	South Africa 27-11-466-1820
Hong Kong 852-2919-1312	Spain 34-91-625-05-40
Hungary 36-30-5052987	Taiwan 886-2-2183-6000
India 91-22-26836733	UK 44-1923-205-166
Israel 972-9-972-4400	USA 1-904-296-9600

Plug-and-play with
the Fusion® ENT
Navigation System

