

## Integra<sup>®</sup> MicroFrance<sup>®</sup>

### Articulated Frontal Sinus Punch Forceps

The Articulated Frontal Sinus Punch Forceps offer a method of accessing the frontal sinus for those who prefer to use manual instruments.

The spring-loaded handle provides leverage for cutting bone while the distal tip remains stable. Additionally, the handle includes a locking mechanism.



# Integra® MicroFrance®

## Articulated Frontal Sinus Punch Forceps

### Ordering Information

Reference	Description	Quantity
3712118	Articulating Frontal Sinus Punch Forceps, 2.0mm, Angled 70° Up-Bite	1
3712119	Articulating Frontal Sinus Punch Forceps, 3.0mm, Angled 70° Up-Bite	1
3712124	Articulating Frontal Sinus Punch Forceps, 4.0mm, Angled 70° Up-Bite	1

The MicroFrance® brand is a premier line of ENT instruments, designed with world-renowned ENT surgeons to facilitate new techniques and improve patient outcomes.

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.
- Consult product labels and inserts for any indication, contraindications, hazards, warnings, precautions, and instructions for use.

**For more information or to place an order, please contact:**

**United States, Canada, Asia, Pacific, Latin America**

USA 800-654-2873 ▪ 888-980-7742 fax

International +1 609-936-5400 ▪ +1 609-750-4259 fax

[integralife.com/contact](http://integralife.com/contact)



Integra York, PA Inc.  
589 Davies Drive  
York, PA 17402 ▪ USA