


Integra[®] MicroFrance[®]


Developed in conjunction with Professor Ludovic De Gabory

Limit uncertainty with
nasal surgery instruments.




De Gabory Fine Sinus Curette


- Interchangeable inserts
- Stylet handle




Jaws



Interchangeable inserts



Stylet handle



Cleaning adaptor MCEN39L (pack of 3)




Integra® MicroFrance®

Instruments developed by MicroFrance® in collaboration with Professor Ludovic De Gabory

Control handle for Luer-lock suction tubes

- Stop or regulate suction as needed



Micro curette with “K” handle

- Handle without closed loops to help free up the hand

Jaws


MCEN4DG



MCEN5DG



MCEN91DG



Integra® MicroFrance®

Instruments developed by MicroFrance® in collaboration with Professor Ludovic De Gabory

Reference	Descriptions (WL: Working Length)
MCEN103LDG	De Gabory Suction Tube Handle 8mm, for Suction Irrigation
MCEN4DG	De Gabory Laryngeal Cup Forceps w/Cleaning Port; WL 110mm (4 ³ / ₈ "), 3.5mm x 9mm, Straight
MCEN5DG	De Gabory Laryngeal Cup Forceps w/Cleaning Port; WL 110mm (4 ³ / ₈ "), 3.5mm x 9.4mm, Angled Up 45°
MCEN68ADG	De Gabory Sinus Curette, 3.4 x 6.4mm Cup Insert, Angled Up, WL 122mm (4 ³ / ₄ ") for MCEN728DG Handle
MCEN68BDG	De Gabory Sinus Curette, 3.5mm x 5.8mm Cup Insert, Angled Up, WL 122 mm (4 ³ / ₄ ") for MCEN728DG Handle
MCEN68CDG	De Gabory Sinus Curette, 3.5mm x 5.8mm Cup Insert, Straight, WL 122 mm (4 ³ / ₄ ") for MCEN728DG Handle
MCEN728DG	De Gabory Modular Handle for Cup Inserts
MCEN91DG	De Gabory Laryngeal Cup Forceps w/Cleaning Port; WL 110mm (4 ³ / ₈ "), 2.8mm x 11mm Cup Jaws, Angled Up 45°

The MicroFrance® brand is a premier line of ENT instruments, designed with world-renowned ENT surgeons to facilitate new techniques and improve patient outcomes.

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.
- Consult product labels and inserts for any indication, contraindications, hazards, warnings, precautions, and instructions for use.

For more information or to place an order, please contact:

United States, Canada, Asia, Pacific, Latin America

USA 800-654-2873 ▪ 888-980-7742 fax

International +1 609-936-5400 ▪ +1 609-750-4259 fax

integralife.com/contact



Integra MicroFrance SAS

Le Pavillon, 03160

Saint-Aubin-le-Monial ▪ France