


## Integra<sup>®</sup> MicroFrance<sup>®</sup>

### Kennedy Frontal Sinus Surgery Instrument Set

Developed in conjunction with David W. Kennedy, MD

#### Developed for Minimally Invasive Frontal Sinus Surgery


This set of malleable instruments has been designed to aid in preservation of mucosa within the frontal recess and frontal sinus.



## Frontal Sinus Surgery Instrument Set

Reference	Description	Quantity
MCSKo**	Kennedy Frontal Sinus Surgery Instrument Set	1 set
MCSK1	Malleable Blunt Probe, Round Tip, Overall Length 152 mm (6")	1
MCSK2	Malleable Sharp Knife, Overall Length 152 mm (6")	1
MCSK3	Malleable Hook, Vertical, Overall Length 152 mm (6")	1
MCSK4	Malleable Hook, Back, Overall Length 152 mm (6")	1
MCSK5	Malleable Suction Tube, 6 French, 2 mm Diameter, Overall Length 123mm (4 7/8")	1
MCSK6	Malleable Curette, 2.0mm Cup, Overall Length 152 mm (6")	1
MCSKoSC	Instrument Tray	1

\*\*Not for sale outside US



MCSK1

MCSK2

MCSK3

MCSK4

The MicroFrance® brand is a premier line of ENT instruments, designed with world-renowned ENT surgeons to facilitate new techniques and improve patient outcomes.

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Always refer to the appropriate instructions for use for complete clinical instructions.
- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.
- Please consult product labels and inserts for any indication, contraindications, hazards, warnings, precautions and instructions for use.

**For more information or to place an order, please contact:**

**United States, Canada, Asia, Pacific, Latin America**

USA 800-654-2873 ■ 888-980-7742 fax

International +1 609-936-5400 ■ +1 609-750-4259 fax

[integralife.com/contact](http://integralife.com/contact)



Integra MicroFrance SAS  
Le Pavillon, 03160  
Saint-Aubin-le-Monial ■ France