

Integra[®] MicroFrance[®]

Lusk MicroBite[™] Forceps

Developed in conjunction with Rodney P. Lusk, MD


The Lusk MicroBite[™] Through-Cutting Sinus Forceps are designed to assist in sharp resection of the anterior and medial wall of the ethmoid bulla while helping to preserve delicate mucosa overlaying the lamina papyracea.



Ordering Information

Reference	Description	Quantity
MCEN76*	Lusk MicroBite Forceps Set	1
Set Includes:		
MCEN73	1.4 x 5.0mm, Straight, 130mm (5 1/8")	1
MCEN74	1.4 x 5.0mm, Angled Up 45°, 126mm (5")	1
MCEN75	1.4 x 5.0mm, Angled Up 90°, 126mm (5")	1

*This set is only available in the US. If you are not in the US, please order the above items individually.



The MicroFrance® brand is a premier line of ENT instruments, designed with world-renowned ENT surgeons to facilitate new techniques and improve patient outcomes.

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.
- Consult product labels and inserts for any indication, contraindications, hazards, warnings, precautions, and instructions for use.

For more information or to place an order, please contact:

United States, Canada, Asia, Pacific, Latin America

USA 800-654-2873 ▪ 888-980-7742 fax

International +1 609-936-5400 ▪ +1 609-750-4259 fax

integralife.com/contact

