


## Integra<sup>®</sup> MicroFrance<sup>®</sup>

### Mini-Struycken Through-Cutting Forceps

These through-cutting forceps with a slim profile assist in accessing the frontal recess and frontal sinus to help remove cell walls and bony partitions. The forceps are designed in both side-to-side and front-to-back biting configurations with 90 and 45 degree bends. They are also designed with 3mm wide rectangular jaws and 1.5mm wide bite. Five instruments are offered to assist the surgeon in working in all quadrants of the frontal ostium and frontal sinus.



# Integra® MicroFrance®

## Mini-Struycken Through-Cutting Forceps

### Ordering Information

Reference	Description	Quantity
MCEN65	1.5 x 7mm, Straight, 116mm (4 <sup>5</sup> / <sub>8</sub> " )	1
MCEN66	1.5 x 6mm, Angled Up 45°, 108mm (4 <sup>1</sup> / <sub>4</sub> " )	1
MCEN67	1.5 x 6mm, Angled Up 90°, 93mm (3 <sup>5</sup> / <sub>8</sub> " )	1
MCEN67-1	1.5 x 6mm, Angled Right 90°, 91mm (3 <sup>5</sup> / <sub>8</sub> " )	1
MCEN67-2	1.5 x 6mm, Angled Left 90°, 91mm (3 <sup>5</sup> / <sub>8</sub> " )	1

The MicroFrance® brand is a premier line of ENT instruments, designed with world-renowned ENT surgeons to facilitate new techniques and improve patient outcomes.

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.
- Consult product labels and inserts for any indication, contraindications, hazards, warnings, precautions, and instructions for use.

**For more information or to place an order, please contact:**

**United States, Canada, Asia, Pacific, Latin America**

USA 800-654-2873 ▪ 888-980-7742 fax

International +1 609-936-5400 ▪ +1 609-750-4259 fax

[integralife.com/contact](http://integralife.com/contact)

**Manufacturer:**



Integra MicroFrance SAS

Le Pavillon, 03160

Saint-Aubin-le-Monial ▪ France