

Integra[®] MicroFrance[®]

Toriumi Rhinoplasty Instrument Set

Designed in conjunction with Dean Toriumi, MD

This 50-piece set comprises a combination of ENT instruments, including the Cartilage Template and Converse Retractor.

The guard and tip of the Toriumi Osteotome is specially shaped to glide along the lateral surface of the bone during lateral osteotomy. Toriumi Curved Rasp's small cutting surface helps with targeted reduction of the dorsal hump and radix. The Suction Elevators help provide visibility during septoplasty.



Integra® MicroFrance®

Toriumi Rhinoplasty Instrument Set

Ordering Information

Reference	Description	Qty
3717008**	Toriumi Rhinoplasty Instrument Set	1 set

Set Includes:

Reference	Description	Qty
CP1100	Curved Rasp, 100mm (3 ⁷ / ₈ "	1
CP1102	Suction Elevator, Sharp, 82mm (3 ¹ / ₄ "	1
CP1103	Suction Elevator, Blunt, 82mm (3 ¹ / ₄ "	1
CP1109	Cartilage Template, Black Delrin, 80 x 90mm, 1 x 1cm Gridlines, Templates, 6.87, 7.65 and 10.3mm	1
CP604BR	Adson Bipolar Forceps, 120mm (4 ³ / ₄ "	1
CP393	Bifurcated Banana Plug Bipolar Cable, not shown	1
3714038	Freer Septum Elevator, Double-Ended, 4.5, 180mm (7"	1
3714125	Joseph Periosteal Elevator, 5mm, 170mm (6 ³ / ₄ "	1
3714127	Cottle Septum Elevator, Graduated, Double-Ended, 217mm (8 ¹ / ₂ "	1
3711002	Lucae Ear Forceps, Delicate, Bayonet, 146mm (5 ³ / ₄ "	1
3711021	Takahashi Nasal Forceps, 4 x 9mm, 115mm (4 ¹ / ₂ "	1
3740039	Brown-Adson Tissue Forceps, 7 x 7 Teeth, 121mm (4 ³ / ₄ "	2
3740028	Adson Tissue Forceps, Delicate, 2 x 3 Teeth, 120mm (4 ³ / ₄ "	1
3740051	Bishop-Harmon Tissue Forceps, 0.3mm, 1 x 2 Teeth, 86mm (3 ³ / ₈ "	1
3740053	Bishop-Harmon Tissue Forceps, Curved, 1 x 2 Teeth, 86mm (3 ³ / ₈ "	1
3755101	Halstead Mosquito Forceps, Straight, 120mm (4 ³ / ₄ "	1
3755102	Halstead Mosquito Forceps, Curved, 125mm (4 ⁷ / ₈ "	1
3714045	Double Fixation Hook, 3mm, 122mm (4 ³ / ₄ "	1
3714021	Joseph Skin Hook, 2-Prong, 5mm, 159mm (6 ¹ / ₄ "	1
3714023	Joseph Skin Hook, 2-Prong, 10mm, 159mm (6 ¹ / ₄ "	1
3714110	Cottle-Joseph Double Hook, 2-Prong, 12mm, 137mm (5 ³ / ₈ "	1
3714020	Joseph Skin Hook, 2-Prong, 2mm, 159mm (6 ¹ / ₄ "	1
3714179	Converse Nasal Retractor, 110mm (4 ³ / ₈ "	1
3714681	Neivert Knife Guide and Retractor, Double-Ended 200mm (7 ³ / ₄ "	1
3744012	Webster Needle Holder, Smooth, 130mm (5 ¹ / ₈ "	1
3714059	Hexagonal Handle Osteotome, 2mm, 159mm (6 ¹ / ₄ "	1


Reference	Description	Qty
3714076	Rubin Osteotome, Finger Grip, 14mm	1
3714077	Rubin Osteotome, Finger Grip, 16mm, not shown	1
3714060	Hexagonal Handle Osteotome, 3mm, 159mm (6 ¹ / ₄ "	1
3740086	Ruler 6", Stainless Steel, Millimeters and Inches	1
3752008	Castroviejo Needle Holder, Delicate, Straight w/Lock, 10mm, 140mm (5 ¹ / ₂ "	1
3714148	Fomon Rasp, Fine, Double-Ended, 210mm (8 ¹ / ₄ "	1
3714149	Fomon Rasp, Coarse, Double-Ended, 210mm (8 ¹ / ₄ "	1
3714063	Hexagonal Handle Osteotome, 5mm, 159mm (6 ¹ / ₄ "	1
3713005	Iris Scissors, Curved, Lower Blade Serrated, 105mm (4 ¹ / ₈ "	1
3713013	Stevens Tenotomy Scissors, Standard, Straight, 112mm (4 ³ / ₈ "	1
3713014	Stevens Tenotomy Scissors, Standard, Curved, 110mm (4 ³ / ₈ "	1
3741033	Stitch Scissors, Ribbon Type, Curved, 98mm (3 ⁷ / ₈ "	1
3741076	Metzenbaum Scissors, Delicate, Curved, 133mm (5 ¹ / ₄ "	1
3741164	Wilmer Conjunctival Scissors, Angled, Sharp, 93mm (3 ⁵ / ₈ "	1
3714201	Vienna Nasal Speculum, Size 3, 30mm Blades	1
3714472	Cottle Septum Speculum, 50mm Blades, 135mm (5 ³ / ₈ "	1
3714206	Killian Septum Speculum, 65mm Blades	1
3714004	Frazier Suction Tube, 10FR. Angled, 100mm (3 ⁷ / ₈ "	1
3714353	Cottle Mallet, 7 OZ., 185mm (7 ¹ / ₄ "	1
3740000	Knife Handle #3, Overall length 125mm (4 ⁷ / ₈ "	2
3750012	Cartilage Crusher, 65mm (2 ¹ / ₂ "	1
3717011	Rhinoplasty Tray, 3-Level, 52.8cm L x 25cm W x 12.3cm D (20 ³ / ₄ " L x 9 ³ / ₄ " W x 4 ⁷ / ₈ " D), not shown	1
3714001	Frazier Suction Tube, 7FR. Angled, 100mm (3 ⁷ / ₈ "	1
CP313-1	Becker Osteotome, Guarded, 3mm, 200mm (7 ³ / ₄ "	1



CP1100
Curved Rasp



CP1102
Suction Elevator, Sharp



CP1103
Suction Elevator, Blunt

**Not for sale outside US





CP1100



CP1102



CP1103



CP1109



3714038



3714125



3714127



3711002



3711021



3740028



3740051



3755101



3714045



3714110



3714020



3714179




3714681



3744012



3714059



3714076



3740086



3752008



3714149



3713013



3714201



3714472



3714206



3714004



3714353



3740000



3750012



CP313-1

Integra® MicroFrance®

Toriumi Rhinoplasty Instrument Set



The MicroFrance® brand is a premier line of ENT instruments, designed with world-renowned ENT surgeons to facilitate new techniques and improve patient outcomes.

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.
- Consult product labels and inserts for any indication, contraindications, hazards, warnings, precautions, and instructions for use.

For more information or to place an order, please contact:

United States, Canada, Asia, Pacific, Latin America

USA 800-654-2873 ▪ 888-980-7742 fax

International +1 609-936-5400 ▪ +1 609-750-4259 fax

integralife.com/contact

