

## Integra<sup>®</sup> MicroFrance<sup>®</sup>

### Wormald Bipolar Forceps with Suction

Developed in conjunction with P.J. Wormald, MD, FRACS, FRCS

Available in two lengths, the MicroFrance<sup>®</sup> Wormald Bipolar Forceps with Suction help achieve hemostasis during endoscopic sinus and anterior skull base surgery. The longer forceps are able to reach the pituitary and are available in straight and up-curved versions. These forceps feature a slim profile and fine tips to help provide access within difficult-to-reach areas. The integrated suction can be advanced independently of the cautery tips to help provide visibility of the operative site, and then retracted from the line of sight during cauterization. An integrated flush port facilitates cleaning and debris removal.



# Integra® MicroFrance®

## Wormald Bipolar Forceps with Suction

Reference	Description	Qty
<b>MCEN64K</b>	<b>Wormald Bipolar Forceps with Suction Set, 105 mm</b>	<b>1 Set</b>
Set includes:		
MCEN64	Wormald Bipolar Forceps with Suction, 105 mm	1
MCEN640	Instrument Tray	1
CP393	Bifurcated (Banana Plug) Bipolar Cable	1

MCEN641



Reference	Description	Qty
<b>MCEN641K</b>	<b>Wormald Bipolar Forceps with Suction Set, 170 mm</b>	<b>1 Set</b>
Set includes:		
MCEN641	Wormald Bipolar Forceps with Suction, 170 mm	1
MCEN640	Instrument Tray	1
CP393	Bifurcated (Banana Plug) Bipolar Cable	1

MCEN643 open tip



Reference	Description	Qty
<b>MCEN643K</b>	<b>Wormald Bipolar Forceps with Suction Set, Up-curved, 170 mm</b>	<b>1 Set</b>
Set includes:		
MCEN643	Wormald Bipolar Forceps, Up-curved, 170 mm	1
MCEN644	Disposable Suction Tip	5
MCEN640	Instrument Tray	1
CP393	Bifurcated (Banana Plug) Bipolar Cable	1

MCEN643 closed tip



The MicroFrance® brand is a premier line of ENT instruments, designed with world-renowned ENT surgeons to facilitate new techniques and improve patient outcomes.

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Always refer to the appropriate instructions for use for complete clinical instructions.
- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.

**For more information or to place an order, please contact:**

**United States, Canada, Asia, Pacific, Latin America**

USA 800-654-2873 ▪ 888-980-7742 fax

International +1 609-936-5400 ▪ +1 609-750-4259 fax

[integralife.com/contact](http://integralife.com/contact)

**Manufacturer:**



Integra MicroFrance SAS

Le Pavillon, 03160

Saint-Aubin-le-Monial ▪ France