

## Integra® MicroFrance®

### Wormald Frontal Sinus Surgery Instrument Set

Developed in conjunction with P.J. Wormald, MD, FRACS, FRCS

This set of instruments is designed to aid in preservation of mucosa within the frontal recess and frontal sinus.

The set includes five instruments and an instrument tray. Each instrument has a working length of 11.0 cm with an outer diameter ranging from 2.0 mm to 2.5 mm. The last 4.0 cm of each instrument is malleable for access across a wide range of anatomies.

A variable suction design on the malleable suction probe, freer suction elevator, and sharp-cupped suction curette allows the surgeon to control the amount of suction through a thumb-controlled aperture.



# Integra® MicroFrance®

## Wormald Frontal Sinus Surgery Instrument Set

### Ordering Information

Product	Description	Quantity
MCEN136**	Wormald Frontal Sinus Surgery Instrument Set	1
Set includes:		
MCEN130	Malleable Sinus Probe, Double-Ended, 260mm (10¼")	1
MCEN131	Malleable Suction Probe, 115mm (4½")	1
MCEN133	Malleable Suction Elevator, Blunt, 115mm (4½")	1
MCEN134	Malleable Suction Curette, Sharp Cup, 115mm (4½")	1
MCEN135	Knot Pusher, 3.0mm, 92mm (3⅝")	1
MCEN149W	Bending Tool, 150mm (5⅞")	1
MCEN136SC	285mm L x 175mm W x 50mm D (11¼" L x 7" W x 2" D)	1

\*\*Not for sale outside US

The MicroFrance® brand is a premier line of ENT instruments, designed with world-renowned ENT surgeons to facilitate new techniques and improve patient outcomes.

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Always refer to the appropriate instructions for use for complete clinical instructions.
- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.
- Please consult product labels and inserts for any indication, contraindications, hazards, warnings, precautions and instructions for use.

For more information or to place an order, please contact:

United States, Canada, Asia, Pacific, Latin America

USA 800-654-2873 ▪ 888-980-7742 fax

International +1 609-936-5400 ▪ +1 609-750-4259 fax

[integralife.com/contact](http://integralife.com/contact)



Integra MicroFrance SAS  
Le Pavillon, 03160  
Saint-Aubin-le-Monial ▪ France



Integra York PA, Inc.  
589 Davies Drive  
York, PA 17402 ▪ USA