


Integra[®] MicroFrance[®]

Wormald Hemorrhage Control Set

Developed in conjunction with P.J. Wormald, MD, FRACS, FRCS

Surgeons can use the vascular clamps to control a large vessel hemorrhage, enabling visualization and sustainable blood flow.



MicroFrance Wormald Hemorrhage Control Set



MCEN190



MCEN191



MCEN192



MCEN193



MCEN194



MCEN195

Wormald Clamps

- Multiple tips and orientations for use with vessel configurations
- Debaquey-style jaws designed to help grasp blood vessels atraumatically
- Ratchet lock on handle helps maintain jaw closure
- Integrated flush port facilitates cleaning and helps debris removal
- Rotatable shaft allows for orientation of jaw or handle

Design helps enable visibility, access and control.



MCEN211* and MCEN211S**




MCEN212* and MCEN212S**

MCEN213* and MCEN213S**

Wormald Needle Holder


- Tungsten-carbide jaw inserts
- Rotatable shaft helps enable orientation
- Pistol-grip ring handle
- Ratchet lock on handle helps maintain jaw closure



MCEN218



MCEN218



MCOHoSC

Wormald Clip-Aplying Forceps

- Rotation design helps allow 360-degree adjustment of the tip and shaft
- Available in straight, 45° jaw and 45° shaft versions

Instrument Tray

54.6 cm L x 25.9 cm W x 12.0 cm D

- Both upper insert and base include silicone mats
- Two 3717022 pick trays will fit inside base
- Includes 6 silicone straps, 12 snap-in posts, and a post removal tool to help stabilize ring-handled instruments

*Compatible with Ethicon® Medium Clips LT200

**Compatible with Ethicon Small Clips LT100

Integra® MicroFrance®

Wormald Hemorrhage Control Set

Ordering Information

Reference	Description	Quantity
MCEN189	Wormald Hemorrhage Control Set*	1
Set Includes:		
MCEN190	Clamp, Straight, Rotatable, 173mm	1
MCEN191	Clamo, Curved, Small, 173mm	1
MCEN192	Clamp, Curved, Long, 167mm	1
MCEN193	Clamp, 45°, Straight, 163mm	1
MCEN194	Clamp, 45°, Curved, Small, 162mm	1
MCEN195	Clamp, 45°, Curved, Long, 162mm	1
MCEN211	Clip-Appling Forceps, Rotatable, 45° Jaw**	1
MCEN211S	Clip-Appling Forceps, Rotatable, 45° Jaw***	1
MCEN212	Clip-Appling Forceps, Rotatable, Straight**	1
MCEN212S	Clip-Appling Forceps, Rotatable, Straight***	1
MCEN213	Clip-Appling Forceps, Rotatable, 45° Shaft**	1
MCEN213S	Clip-Appling Forceps, Rotatable, 45° Shaft***	1
MCEN218	Needle Holder, Rotatable, 170mm	1
MCOHoSC	Instrument Tray, 54.6 cm L x 25.9 cm W x 12.0 cm D	1

*Note: This set is only available in the US. If you are not in the US, please order the above items individually.

** Compatible with Ethicon Medium Clips LT-200

***Compatible with Ethicon Small Clips LT100

The MicroFrance® brand is a premier line of ENT instruments, designed with world-renowned ENT surgeons to facilitate new techniques and improve patient outcomes.

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.
- Consult product labels and inserts for any indication, contraindications, hazards, warnings, precautions, and instructions for use.

For more information or to place an order, please contact:

United States, Canada, Asia, Pacific, Latin America

USA 800-654-2873 ▪ 888-980-7742 fax

International +1 609-936-5400 ▪ +1 609-750-4259 fax

integralife.com/contact



MicroFrance, Integra and the Integra logo are registered trademarks of Integra LifeSciences Corporation or its subsidiaries in the United States and/or other countries. All other trademarks and trade names are the property of their respective owners.

©2017 Integra LifeSciences Corporation. All rights reserved. Printed in USA. oM

0582632-1-EN