

## Integra<sup>®</sup> MicroFrance<sup>®</sup>

### California Ear Institute Ossiculoplasty Sizer Set

A Sizer Set for Ossicular Reconstruction Surgery

Developed in conjunction with Joseph B. Roberson, MD,  
of the California Ear Institute at Stanford

This ossiculoplasty sizer set includes measuring instruments to use as a guide for cutting prostheses that may help ensure the selection of the correct pre-sized prostheses. The instrument set includes an instrument tray and 13 malleable sizers from 2.0 mm to 8.0 mm in 0.5 mm increments.



# Integra® MicroFrance®

## California Ear Institute Ossiculoplasty Sizer Set

### Ordering Information

Product	Description	Quantity
MCOR15*	California Ear Institute Ossiculoplasty Sizer Set	1 set
<b>Set Includes:</b>		
MCOR20	Ossiculoplasty Sizer, 2.0 mm	1
MCOR25	Ossiculoplasty Sizer, 2.5mm	1
MCOR30	Ossiculoplasty Sizer, 3.0 mm	1
MCOR35	Ossiculoplasty Sizer, 3.5 mm	1
MCOR40	Ossiculoplasty Sizer, 4.0 mm	1
MCOR45	Ossiculoplasty Sizer, 4.5 mm	1
MCOR50	Ossiculoplasty Sizer, 5.0 mm	1
MCOR55	Ossiculoplasty Sizer, 5.5 mm	1
MCOR60	Ossiculoplasty Sizer, 6.0 mm	1
MCOR65	Ossiculoplasty Sizer, 6.5 mm	1
MCOR70	Ossiculoplasty Sizer, 7.0 mm	1
MCOR75	Ossiculoplasty Sizer, 7.5 mm	1
MCOR80	Ossiculoplasty Sizer, 8.0 mm	1
MCOR15SC	Instrument Tray	1

\*This set is only available in the US. If you are not in the US, please order the above items individually.

The MicroFrance® brand is a premier line of ENT instruments, designed with world-renowned ENT surgeons to facilitate new techniques and improve patient outcomes.

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.
- Consult product labels and inserts for any indication, contraindications, hazards, warnings, precautions, and instructions for use.

**For more information or to place an order, please contact:**

**United States, Canada, Asia, Pacific, Latin America**

USA 800-654-2873 ▪ 888-980-7742 fax

International +1 609-936-5400 ▪ +1 609-750-4259 fax

[integralife.com/contact](http://integralife.com/contact)



Integra MicroFrance SAS  
Le Pavillon, 03160  
Saint-Aubin-le-Monial ▪ France



Integra York, PA Inc.  
589 Davies Drive  
York, PA 17402 ▪ USA