

Integra[®] MicroFrance[®] Josephson Suction Round Knife

Developed in conjunction with Gary D. Josephson, MD

The Josephson Suction Round Knife is used for otologic surgery. It has suction on the knife plane which may help the surgeon's line of sight stay clear of fluid and may help prevent the surgeon from having to hold a separate suction tube in the other hand. This may thus assist the surgeon in utilizing more instruments simultaneously during otologic procedures.



Integra® MicroFrance®

Josephson Suction Round Knife

Ordering Information

Product	Description	Quantity
MCO5D	Josephson Round Knife with Suction, 2 mm, Overall Length 180mm (7")	1
MCO5E	Josephson Round Knife with Suction, 3 mm, Overall Length 180mm (7")	1

The MicroFrance® brand is a premier line of ENT instruments, designed with world-renowned ENT surgeons to facilitate new techniques and improve patient outcomes.

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.
- Consult product labels and inserts for any indication, contraindications, hazards, warnings, precautions, and instructions for use.

For more information or to place an order, please contact:

United States, Canada, Asia, Pacific, Latin America

USA 800-654-2873 ▪ 888-980-7742 fax

International +1 609-936-5400 ▪ +1 609-750-4259 fax

integralife.com/contact



Integra MicroFrance SAS

Le Pavillon, 03160

Saint-Aubin-le-Monial ▪ France