

# Integra® MicroFrance®

## Kraus TM MicroCut® Forceps

Developed in conjunction with Eric M. Kraus, MD, MS, FACS

The Kraus TM MicroCut Forceps were designed to help remove the epithelial edge of a tympanic membrane perforation without over-enlarging the existing perforation.

- Useful during tympanic membrane patching, myringoplasty, and tympanoplasty procedures
- Designed for precision edge cutting
- Low profile rotation knob allows for universal 360° cutting ability and promotes visualization
- Thumb-activated reverse-action handle allows for controlled movement and helps reduce distal tip motion
- Distal tip basket captures resected tissue to help prevent tissue from entering the middle ear



# Integra® MicroFrance®

Kraus TM MicroCut® Forceps

## Ordering Information

Reference	Description	Quantity
MCO270	Kraus TM MicroCut Forceps, 90° Up-Biting	1
MCO271	Kraus TM MicroCut Forceps, Forward-Biting	1
MCO272	Kraus TM MicroCut Forceps, Backward-Biting	1

MCO270



MCO271



MCO272



The MicroFrance® brand is a premier line of ENT instruments, designed with world-renowned ENT surgeons to facilitate new techniques and improve patient outcomes.

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.
- Consult product labels and inserts for any indication, contraindications, hazards, warnings, precautions, and instructions for use.

For more information or to place an order, please contact:

United States, Canada, Asia, Pacific, Latin America

USA 800-654-2873 ▪ 888-980-7742 fax

International +1 609-936-5400 ▪ +1 609-750-4259 fax

[integralife.com/contact](http://integralife.com/contact)

**INTEGRA**   
LIMIT UNCERTAINTY