


Integra[®] MicroFrance[®]

Teixido Cartilage Cutter

Developed in conjunction with Michael Teixido, MD

The Teixido Cartilage Cutter provides a simple solution to preparing cartilage grafts for tympanoplasty and ossicular chain reconstruction procedures. This device slices cartilage in a standard thickness of 0.3 mm in one swift motion. Working space for trimming and measuring the cartilage is also incorporated to facilitate graft preparation.

- Minimum number of parts to clean
- Angled, replaceable blade
- 0.3mm cartilage slice allows for 2-3 slices per tragus



Integra® MicroFrance®

Teixido Cartilage Cutter

Teixido Cartilage Cutter

Reference	Description	Qty
MCO255	Teixido Cartilage Cutter with One Blade	1
MCO251	Replacement Blades *	10

*Not for sale in the US



The MicroFrance® brand is a premier line of ENT instruments, designed with world-renowned ENT surgeons to facilitate new techniques and improve patient outcomes.

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.
- Consult product labels and inserts for any indication, contraindications, hazards, warnings, precautions, and instructions for use.

For more information or to place an order, please contact:

United States, Canada, Asia, Pacific, Latin America

USA 800-654-2873 ■ 888-980-7742 fax

International +1 609-936-5400 ■ +1 609-750-4259 fax

integralife.com/contact



Integra MicroFrance SAS

Le Pavillion, 03160

Saint-Aubin-Le-Monial ■ FRANCE