

Integra[®] MicroFrance[®]

Benninger Microlaryngeal Instruments

Developed in conjunction with Michael Benninger, MD



A New Level of Versatility

Laryngeal instruments that combine quality, precision, and stability

The Benninger Microlaryngeal Instruments are designed to help provide stabilization of the surgeon's hand. The reverse action handles are thumb-activated to help reduce the amount of distal tip movement and wrist movement required by the surgeon for atraumatic dissection. This set includes forceps, knives, retractors, injection instruments, suction cannulas, and spatulas for microlaryngeal manipulation and dissection.

Ordering Information

Product	Description	Qty
MCLB2**	MicroFrance® Benninger Microlaryngeal Instrument Set	1 set
Set Includes:		
MCL-S21R	Grasper, Right, Heart-Shaped	1
MCL-S22R	Grasper, Left, Heart-Shaped	1
MCL-S23R	Alligator Forceps, Right	1
MCL-S24R	Alligator Forceps, Left	1
MCLS18R	Cup Forceps, Downbiting, 2.3 mm	1
MCLB94	Laryngeal Needle Driver, Angled 70° Up	1
MCLB30	Laryngeal Suture Knot Pusher	1
MCLB16	Laryngeal Micro-Scissors, Curved Right, 1.5mm to 8.5mm Blades	1
MCLB17	Laryngeal Micro-Scissors, Curved Left, 1.5mm to 8.5mm Blades	1
MCLB31	Laryngeal Micro-Scissors, Angled 70° Up, 1.5 mm Blades	1
Sataloff Microlaryngeal Instruments Included in this Set:		
MCL-S0	Suction/Irrigation Adapter	1
MCL-S1	Elevator, 1mm x 9mm, Curved, 252mm (10")	1
MCL-S7	Spatula, Straight, Elevator, 1mm x 10mm, Straight, 252mm (10")	1
MCL-S4	Injection Needle, Straight	6
MCL-S25	Sickle Knife <i>For use with MCL-S19</i>	6
MCL-S19	Universal Knife Handle, Ring Handle with Finger Loop	1
MCL-S16	Suction, 3 French	1
3734109	Suction, 5 French	1
MCLB2SC	Benninger Instrument Tray	1

**Not for sale outside US

The MicroFrance® brand is a premier line of ENT instruments, designed with world-renowned ENT surgeons to facilitate new techniques and improve patient outcomes.

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.
- Consult product labels and inserts for any indication, contraindications, hazards, warnings, precautions, and instructions for use.
- Please consult product labels and inserts for any indication, contraindications, hazards, warnings, precautions and instructions for use.

For more information or to place an order, please contact:

United States, Canada, Asia, Pacific, Latin America

USA 800-654-2873 ■ 888-980-7742 fax

International +1 609-936-5400 ■ +1 609-750-4259 fax


integralife.com/contact



Integra MicroFrance SAS
Le Pavillon, 03160
Saint-Aubin-le-Monial ■ France



Integra York PA, Inc.
589 Davies Drive
York, PA 17402 ■ USA



- Reverse action handle
- Helps reduce distal tip movement