


## Integra<sup>®</sup> MicroFrance<sup>®</sup>

### Laryngeal Bipolar Forceps

The Laryngeal Bipolar Forceps are designed for laryngology and phonosurgery procedures.

The coagulating tip is 5mm in length, and the design helps ensure that the bipolar energy is delivered in a controlled manner. The 29cm length and the 3mm outer diameter make these forceps useful for delicate work. A specially designed stainless steel instrument tray is also included.



# Integra® MicroFrance®

## Laryngeal Bipolar Forceps

### Ordering Information

Product	Description	Quantity
MCL34	Bipolar Forceps with Box & Tray, 290mm (11½")	1
CP393	Bifurcated (Banana Plug) Bipolar Cable, 4.5m (14'–¾")	1

The MicroFrance® brand is a premier line of ENT instruments, designed with world-renowned ENT surgeons to facilitate new techniques and improve patient outcomes.

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.
- Consult product labels and inserts for any indication, contraindications, hazards, warnings, precautions, and instructions for use.

**For more information or to place an order, please contact:**

**United States, Canada, Asia, Pacific, Latin America**

USA 800-654-2873 ▪ 888-980-7742 fax

International +1 609-936-5400 ▪ +1 609-750-4259 fax

[integralife.com/contact](http://integralife.com/contact)



Integra MicroFrance SAS

Le Pavillon, 03160

Saint-Aubin-le-Monial ▪ France