


# Integra<sup>®</sup> ENT

## Mustard Table

Helps enhance Surgeon Control for Suspension Procedures





**Developed in conjunction with  
David S. Parsons, MD, FAAP, FACS**

The Mustard Table is an over-the-chest operating table that helps improve suspension during surgical procedures, such as: direct laryngoscopy, tonsillectomy, adenoidectomy, cleft palate repair, and other oral procedures.

When the patient is suspended, the Mustard Table allows you to move the OR table and patient in multiple directions.

- Easily attaches directly to OR table sideframe.
- Helps enhance surgeon control during surgical procedures such as laryngoscopy
- Edge slots enhance ease of adjustment and disconnection as needed

**Ordering Information**

Reference	Description	Qty.
3734500	Mustard Table	1
3734510	Mustard Table adapter	1

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.
- Consult product labels and inserts for any indication, contraindications, hazards, warnings, precautions, and instructions for use.

**For more information or to place an order, please contact:**

**United States, Canada, Asia, Pacific, Latin America**

USA 800-654-2873 ▪ 888-980-7742 fax  
International +1 609-936-5400 ▪ +1 609-750-4259 fax  
[integralife.com/contact](http://integralife.com/contact)

 Integra York PA, Inc.  
589 Davies Drive  
York, PA 17402 ▪ USA