

Integra[®] MicroFrance[®]

Woo Laryngeal Scissors

The MicroFrance[®] Woo Laryngeal scissors have a double-action lever mechanism, which allows for effective cutting of fibrous tissue. This design provides stability and limited rotation during the cutting process.

The scissors are also short enough to be balanced on the edge of the laryngoscope. The fine triangular tip may be used for fine dissection.



Integra® MicroFrance®

Woo Laryngeal Scissors

Reference	Description	Qty
MCLW10	MicroFrance Woo Laryngeal Scissors, Straight	1
MCLW11	MicroFrance Woo Laryngeal Scissors, 30° Up	1

MCLW11



MCLW10



The MicroFrance® brand is a premier line of ENT instruments, designed with world-renowned ENT surgeons to facilitate new techniques and improve patient outcomes.

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Always refer to the appropriate instructions for use for complete clinical instructions.
- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.

For more information or to place an order, please contact:

United States, Canada, Asia, Pacific, Latin America

USA 800-654-2873 ▪ 888-980-7742 fax

International +1 609-936-5400 ▪ +1 609-750-4259 fax

integralife.com/contact

Manufacturer:



Integra MicroFrance SAS

Le Pavillon, 03160

Saint-Aubin-le-Monial ▪ France