

Integra[®] MicroFrance[®]

Wormald Anterior Skull Base Clip-Appling Forceps

Developed in conjunction with P.J. Wormald, MD, FRACS, FRCS

The Wormald Clip-Appling Forceps help enable access to the anterior skull base region. The rotation design allows you to adjust the orientation of the tip as well as the handle 360° for multiple instrument access. The forceps have a working length of 17 centimeters and are available in straight 45° jaw, and 45° shaft versions with small and medium clip options.



Integra® MicroFrance®

Wormald Anterior Skull Base Clip-Aplying Forceps

Ordering Information

Compatible with Ethicon® Medium Clips LT200

Working length 182mm (7¼"), W/Cleaning Port

Reference	Description
MCEN211	Clip-Aplying Forceps, Rotatable, 45° Jaw
MCEN212	Clip-Aplying Forceps, Rotatable, Straight
MCEN213	Clip-Aplying Forceps, Rotatable, 45° Shaft



MCEN211 and MCEN211S

Compatible with Ethicon Small Clips LT100

Working length 182mm (7¼"), W/Cleaning Port

Reference	Description
MCEN211S	Clip-Aplying Forceps, Rotatable, 45° Jaw
MCEN212S	Clip-Aplying Forceps, Rotatable, Straight
MCEN213S	Clip-Aplying Forceps, Rotatable, 45° Shaft



MCEN212 and MCEN212S



MCEN213 and MCEN213S

The MicroFrance® brand is a premier line of ENT instruments, designed with world-renowned ENT surgeons to facilitate new techniques and improve patient outcomes.

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.
- Consult product labels and inserts for any indication, contraindications, hazards, warnings, precautions, and instructions for use.

For more information or to place an order, please contact:

United States, Canada, Asia, Pacific, Latin America

USA 800-654-2873 ■ 888-980-7742 fax

International +1 609-936-5400 ■ +1 609-750-4259 fax

integralife.com/contact

