

# Integra<sup>®</sup> MicroFrance<sup>®</sup>

## Wormald Anterior Skull Base and Pituitary Instrument Set


Developed in conjunction with P.J. Wormald, MD, FRACS, FRCS

This instrument set is intended for both otolaryngologists and neurosurgeons in a combined approach to endoscopic anterior skull base and pituitary procedures.


This set includes elongated malleable instruments with suction. The set also features a variety of reverse-action scissors of various sizes and orientations and forceps.

### Key Features

- Scissors have long, slim, profiles to aid in visibility within the intracranial space
- Blakesley-style forceps to aid in tissue removal
- Smaller diameter instruments help enable access to the anterior skull region
- Pituitary forceps aid in tissue removal
- Malleable instruments
- Large on/off suction aperture helps prevent unintentional suction



## Scissors & Forceps



Scissors, Up,  
5.0 mm Blade  
**MCLW11**



Scissors, Straight,  
5.0 mm Blade  
**MCLW10**



Scissors, Up,  
8.0 mm Blade  
**MCEN217U**



Scissors, Straight,  
8.0 mm Blade  
**MCEN217S**



Scissors, Right,  
8.0 mm Blade  
**MCEN217R**



Scissors, Left,  
8.0 mm Blade  
**MCEN217L**



Scissors, Right,  
6.0 mm Blade  
**MCEN216R**



Scissors, Left,  
6.0 mm Blade  
**MCEN216L**



Forceps, 45°  
**MCEN215**




Forceps, Straight  
**MCEN214**

## Suction



Malleable Suction Probe  
MCEN138



Malleable Probe, 90°  
MCEN139



Malleable Probe, Straight  
MCEN139A



Malleable Suction  
Curette, 45°, Small  
MCEN140



Malleable Suction  
Curette, 45°, Large  
MCEN140A



Malleable Suction  
Curette, 90°, Small  
MCEN141



Malleable Suction  
Curette, 90°, Large  
MCEN141A



Malleable Suction Basket  
MCEN142



Malleable Suction Elevator  
MCEN147

### Ordering Information

Product	Description	Quantity
MCEN138	Malleable Suction Probe, 3.0mm, 176mm (7")	1
MCEN139	Malleable Probe, 2.5mm, 90°, 170mm (6¾")	1
MCEN139A	Malleable Probe, 1.0mm, Straight, 170mm (6¾")	1
MCEN140	Malleable Suction Curette, 45°, Small, 175mm (7")	1
MCEN140A	Malleable Suction Curette, 45°, Large, 175mm (7")	1
MCEN141	Malleable Suction Curette, 90°, Small, 175mm (7")	1
MCEN141A	Malleable Suction Curette, 90°, Large, 175mm (7")	1
MCEN142	Malleable Suction Basket, 177mm (7")	1
MCEN143	Instrument Tray, 52cm L x 24.5cm W x 12cm D, (20½" L x 9¾" W x 4¾" D)	1
MCEN147	Malleable Suction Elevator, 3.5mm, 176mm (7")	1
MCEN214	Pituitary Forceps, 1.5 x 6.0mm, Straight, 168mm (6½")	1
MCEN215	Pituitary Forceps, 1.5 x 4.5mm, 45° Up, 166mm (6½")	1
MCEN216L	Pituitary Scissors, 6.0mm, Left, 168mm (6½")	1
MCEN216R	Pituitary Scissors, 6.0mm, Right, 168mm (6½")	1
MCLW10	Pituitary Scissors, 5.0mm, Straight, 202mm (8")	1
MCLW11	Pituitary Scissors, 5.0mm, Up, 202mm (8")	1
MCEN217L	Pituitary Scissors, 8.0mm, Left, 170mm (6¾")	1
MCEN217R	Pituitary Scissors, 8.0mm, Right 170mm (6¾")	1
MCEN217S	Pituitary Scissors, 8.0mm, Straight, 170mm (6¾")	1
MCEN217U	Pituitary Scissors, 8.0mm, Up, 170mm (6¾")	1
MCEN149W	Bending Tool, 150mm (5⅞")	1

The MicroFrance® brand is a premier line of ENT instruments, designed with world-renowned ENT surgeons to facilitate new techniques and improve patient outcomes.

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Always refer to the appropriate instructions for use for complete clinical instructions.
- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.
- Please consult product labels and inserts for any indication, contraindications, hazards, warnings, precautions and instructions for use.

### For more information or to place an order, please contact:

#### United States, Canada, Asia, Pacific, Latin America

USA 800-654-2873 ▪ 888-980-7742 fax

International +1 609-936-5400 ▪ +1 609-750-4259 fax

[integralife.com/contact](http://integralife.com/contact)



Integra MicroFrance SAS  
Le Pavillon, 03160  
Saint-Aubin-le-Monial ▪ France



Integra York PA, Inc.  
589 Davies Drive  
York, PA 17402 ▪ USA