


Integra[®] MicroFrance[®]

Stamm Skull Base Instrument Set

Developed in conjunction with Aldo Stamm, MD, PhD

The Stamm Skull Base Instrument Set is commonly used by otolaryngologists and neurosurgeons in a combined approach to transnasal endoscopic skull base surgery and pituitary procedures. This set consists of 16 instruments, including elevators, scissors, suctions, cup forceps, a hook, probes, and bipolar forceps.

- Pencil-grip squeeze handles help reduce distal tip movement in delicate intracranial dissection
- Elongated working lengths (140 mm – 160 mm) to access challenging skull base anatomy
- Tips facilitate fine, complex intradural and extradural dissection
- Keyhole aperture on suction instruments allows for variable suction capability
- Flush ports on each lumen aid in cleaning and debris removal



Integra® MicroFrance®

Stamm Skull Base Instrument Set

Ordering Information

Product	Description	Working Length	Qty
MCEN320*	Stamm Skull Base Instrument Set	-	1
Set Includes:			
MCEN328S	Cup Forceps, Straight	155 mm	1
MCEN328U	Cup Forceps, Up	155 mm	1
MCEN329	Hook, 90-Deg, Sharp	160 mm	1
MCEN330	Ball Seeker, 90-Deg	160 mm	1
MCEN344S	Scissors, Bayonet, Straight, 5 mm Blades	155 mm	1
MCEN344U	Scissors, Bayonet, 45-Deg, 5 mm Blades	155 mm	1
MCEN344V	Scissors, Bayonet, Curved Left, 5 mm Blades	155 mm	1
MCEN345	Elevator, 2 mm Width	160 mm	1
MCEN346	Elevator, 4 mm Width	160 mm	1
MCEN347	Malleable Suction, Straight, Large	150 mm	1
MCEN348	Malleable Suction, Straight, Small	150 mm	1
MCEN349	Malleable Suction, Curved	150 mm	1
MCEN350	Suction Tube, Angulated, Curved Down	150 mm	1
MCEN351	Suction Tube, Angulated, Curved Up	150 mm	1
MCEN364	Bipolar Forceps, 45-Deg	140 mm	1
MCEN365	Bipolar Forceps, Bayonet, Straight	140 mm	1
MCOS10SC	Instrument Tray, Two-Level 38.5 cm L x 25.6 cm W x 12.8 cm D	-	1
CP393	Bifurcated (Banana Plug) Bipolar Cable	-	1

*This set is only available in the US. If you are not in the US, please order the above items individually.




MCOS10SC

The MicroFrance® brand is a premier line of ENT instruments, designed with world-renowned ENT surgeons to facilitate new techniques and improve patient outcomes.

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Always refer to the appropriate instructions for use for complete clinical instructions.
- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.
- Please consult product labels and inserts for any indication, contraindications, hazards, warnings, precautions and instructions for use.



For more information or to place an order, please contact:

United States, Canada, Asia, Pacific, Latin America

USA 800-654-2873 ▪ 888-980-7742 fax

International +1 609-936-5400 ▪ +1 609-750-4259 fax

integralife.com/contact



Integra MicroFrance SAS
Le Pavillon, 03160
Saint-Aubin-le-Monial ▪ France



Integra York PA, Inc.
589 Davies Drive
York, PA 17402 ▪ USA