


Integra[®] MicroFrance[®] Robier Retractor

The MicroFrance Robier Retractor has features designed for access to the middle fossa. Adjustable components allow customization to your patient's anatomy to help improve visualization and stability.

- Malleable blade with integrated suction is adjustable in length, height and rotation
- Hinged arms help reduce profile to allow more working space at the surgical site



Ordering Information

Reference	Description	Qty
MCO357A	Retractor with Malleable Integrated Suction Blades, Overall length 181mm (7¼")	1



The MicroFrance® brand is a premier line of ENT instruments, designed with world-renowned ENT surgeons to facilitate new techniques and improve patient outcomes.

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.
- Consult product labels and inserts for any indication, contraindications, hazards, warnings, precautions, and instructions for use.

For more information or to place an order, please contact:

United States, Canada, Asia, Pacific, Latin America
USA 800-654-2873 ▪ 888-980-7742 fax
International +1 609-936-5400 ▪ +1 609-750-4259 fax
integralife.com/contact

