



Integra[®] MicroFrance[®]

Integra Laparoscopy catalogue

Table of contents	ii	Scissors ø 10 mm	29
Ordering and warranty informations	iv	Complete non modular instrument	29
Indications – Contraindications	v	Modular instruments	29
MicroFrance® modular Instruments	vi	Vascular clamp 5 mm, ø 12 mm	30
Fenestrated forceps ø 5 mm	2	Modular instruments	30
Complete non modular instrument	3, 5	Vascular clamp ø 5 mm	31
Modular instruments	2, 4, 6	Vascular clamp ø 12 mm	31
Graduated fenestrated forceps ø 5 mm	7	Vascular clamp ø 10 mm	32
Fenestrated forceps ø 10 mm	8	Complete non modular instrument	33
Complete non modular instrument	9	Modular instrument	32
Modular instruments	8	Intestinal clamp	34
Babcock forceps	10	Complete non modular instrument	35
Complete non modular instrument	11	Modular instruments	34
Modular instruments	10	Aortic clamp ø 6mm	36
Alligator forceps ø 5 mm	12	Temporary Clamps, Clamp applicator	38
Complete non modular instrument	13	Clamp applicator	38
Modular instruments	12	Clamp applicator insert	38
Teeth forceps	14	Temporary Clamps	38
Complete non modular instrument	15	Clips applier, Circular suturing holder insert	39
Modular instruments	14	Circular suturing holder insert	39
Dissector ø 5 mm	17	Clips applier	39
Complete non modular instrument	19	Rotatable angled handles, ø 5 mm	40
Modular instruments	17	Tubes, ø 5 mm for rotatable angled handles	41
Dissector ø 10 mm	20	Insulated tubes	41
Complete non modular instrument	21	Metal tube and Pack of 3 black sheaths	41
Modular instruments	20	Monopolar cables	41
Needle-Holder 5 mm	22	Rotatable axial handles, ø 5 mm	42
Complete non modular instrument	23	Tubes, ø 5 mm for rotatable axial handles	43
Modular instruments	22	Insulated tubes	43
Needle-Holder 10 mm	24	Metal tube and Pack of 3 black sheaths	43
Knot guide scissors, Knot Pusher	25	Monopolar cables	43
Knot guide scissors	25	Rotatable ergonomic handles, ø 5 mm	44
Knot guide scissors insert	25	Tubes, ø 5 mm for rotatable ergonomic handles	45
Knot Pusher	25	Insulated tubes	45
Scissors ø 5 mm	26	Metal tube and Pack of 3 black sheaths	45
Black sheath for metal tube	27	Monopolar cables	45
Complete non modular instrument	28	Rotatable handles, for 5 and ø 10 mm non-insulated tube	46
Insulated scissors inserts	27	Angled handle	46
Metal tube	27	Axial handle	46
Modular instruments	26	Ergonomic handle	46
Wire insert	27		

Non-insulated tubes	47	Electric irrigation suction pump	60	3,5mm diameter instrumentation	72
Non rotatable axial handles, Ø 5 mm and insulated tubes	48	Cleaning sleeve for telescope	60	Adapter cleaning luer-lock	72
Angled handle	48	Electric handle	60	Bipolar forceps	74
Axial handle	48	Electric irrigation suction pump	60	Complete non modular instrument	73
Black insulated tube for non-rotatable instrument	48	Sterile suction irrigation tubes	60	Forceps without port	75
Non rotatable ergonomic handles and non insulated tubes	49	Suction irrigation tubes	60	Insert with insulated tube	72
Ergonomic handle	49	Suction irrigation tube	61	Knot pusher	75
Non-insulated tube for non-rotatable instrument	49	Handles for tube	61	Monopolar hook	75
Modular rotatable bipolar forceps	50	Modular suction irrigation tubes	61	Palpation probe	75
Bipolar ergonomic handle	50	Suction irrigation tubes	61	Suction irrigation tubes	75
Bipolar inserts	50	Needles - Wire-Pass	62	Trocar	76
Protection part	50	Articulated wire-pass	62	Pediatric instrumentation, 3,5mm diameter	77
Tubes, Ø 5 mm for rotatable bipolar forceps	51	Insufflation needle	62	Adapter cleaning luer-lock	78
Bipolar cables	51	Mintz Set for puncture	62	Bipolar forceps	80
Insulated tubes	51	Puncture needle	62	Complete non modular instrument	79
Metal tube and Pack of 3 black sheaths	51	Retractable wire pass	62	Discision forceps	77
Modular non rotatable bipolar forceps, Ø 5 mm	52	Cholangiography - Choledoscopy	63	Insert with insulated tube	78
Complete non modular instrument	53	Cholangiography	63	Monopolar hook	77
Handles	52	Choledoscopy	63	Monopolar instruments	81
Inserts	52	Modular trocars - 5,5 mm	64	Needle-holder	81
Tubes	52	Modular trocars - 10,5 mm	65	Retractable knife	77
Monopolar hook	56	Modular trocars - 12/13 mm	66	Suction irrigation tubes	77
Complete hook	56	Trumpet valve trocars	67	Henry set, 3,0 mm diameter	82
Hook handle	56	Complete trocar	67	Fenestrated forceps	82
Modular hook	56	Sealing cap	67	Forceps, alligator jaws	82
Knives, Monopolar electrode, Palpation probe	57	Automatic valve trocars	68	Monopolar hook	82
Knife blades	57	Trocar accessories	69	Round probe	82
Monopolar cables	57	Double retractor	69	Scissors	82
Monopolar probe	57	Reducer sleeve for trocar	69	Trocar	82
Retractable knife	57	Sealing cap for trocar	69	Trocar holder ring	82
Retractable monopolar electrode	57	Thoracic instrumentation	70	Terris set	83
Retractors - Palpation probes	58	Fenestrated forceps	70	Double retractor	83
Articulated retractor	58	Forceps, curette jaws	70	Elevator	83
Flat palpation probe	58	Insert, curette jaws	70	Nerve suction with finger control	83
Flexible retractor	58	Trocar	70	Peanut holder	83
Round palpation probe	58	Irrigation bipolar forceps	71	Suction tube elevator	83
Flexible tension arm	59	3,5mm diameter instrumentation	72	Scheye rectal forceps	84
Flexible tension arm	59	Adapter cleaning luer-lock	72	Axial handle	84
Instruments holder	59	Bipolar forceps	74	Insert, biopsy jaws and tube 5mm diam	84
Operating table clamp	59	Complete non modular instrument	73	Sterilization case, other accessories	85
		Forceps without port	75	Cap	85
		Insert with insulated tube	72	Plastic sterilization case	85
		Knot pusher	75	Stainless steel sterilization case	85
		Monopolar hook	75	Stainless steel sterilization cases	85
		Palpation probe	75	Surgical instrument oil	85
		Suction irrigation tubes	75	Tray	85
		Trocar	76	Numerical INDEX	86

Pricing :

Prices stated in specific written quotations are firm for thirty days from the date given, and are otherwise subject to change without prior notice. Pricing terms stated on written agreements are governed by such agreements.

Minimum Order Requirements :

The minimum ordering amount is €250 (Integra reserves the right to apply €20 handling fees and transport fees for any order below €250). For partial unit orders, please contact the customer service to obtain a quote.

Administrative fees :

Administrative fees may be charged for any order.

Ordering Procedure :

A written purchase order is required for the purchase of Integra MicroFrance® products.

Acceptance of Orders :

Orders are accepted upon approval by Integra Customer Service.

Return Policy :

- Authorization from Customer Service must be obtained prior to returning product.
- Sterile product must be returned in unopened, undamaged cartons, packed to prevent damage.
- Non-sterile product must be returned in unused saleable condition in original package.
- A credit note will be issued for products returned less than ninety days from shipment date with a 20% restocking fee, provided that the product being returned is not damaged and it is possible to verify that the product has not been used or opened.

Integra Limited Warranty

Please refer to the warranty included in the packaging of each product to know the applicable warranty.

Contact Customer service :

France / International

+33 (0) 470 67 98 00
+33 (0) 470 67 98 20 (Fax)
microfrance@integralife.com

Benelux

+32 (0) 2 257 4130
+32 (0) 2 253 2466 (Fax)
custsvcbenelux@integralife.com

United Kingdom

+44 (0) 1264 312 725
+44 (0) 1264 312 821 (Fax)
custsvcuk@integralife.com

Germany

+49 (0) 2102 5535 6200
+49 (0) 2102 5536 636 (Fax)
custsvcgermany@integralife.com

MicroFrance® Laparoscopy Instrumentation

INDICATIONS

- All the medical devices showed in this catalogue have indications for diagnosis and/or laparoscopy surgery.
- Sizes are not contractual and are given for information purposes only.

CONTRAINDICATIONS

- Do not use in indications other than the indications mentioned in the DFU or Labeling.

How to read this catalogue

Title and/or product family Symbol
"complete instrument"

"Modular" symbol

Instrument to a 1/1 scale except where mentioned otherwise

Mors Echelle (s)		Longueur utile			N° Pages
		Court	Standard	Long	
Démontable					
Johann 200/30mm					
		CEV625D 345 mm	CEV625 345 mm	CEV625+ 465 mm	
					Rotatif 40 -> 47
					Non rotatif 48 - 49
		CEV625D 345 mm	CEV625 345 mm	CEV625D 465 mm	
		X	CEV625+ 345 mm	X	
		CEV625C 345 mm	CEV625C 345 mm	CEV625C 465 mm	

Mors Echelle (s)		Longueur utile			
		Court	Standard	Long	
NON Démontable					
Instruments complets					
Johann 200/30mm					
		Monoplaire à ciseaux	Montsupet CEV625M 345 mm	CEV625 345 mm	X
		A ressort	X	Gayet CEV625G 345 mm	X
		A ciseaux	Montsupet CEV625 345 mm	X	X
		A ciseaux	CEV625C 345 mm	CEV625 345 mm	X
		Monoplaire à ciseaux	X	CEV625M 345 mm	X
		A ciseaux longue	X	CEV625 345 mm	X
		X	X	X	X
		Monoplaire à ciseaux	X	CEV625 345 mm	X
		A ciseaux	X	CEV625C 345 mm	X

DOCUMENT DESTINE A L'EUROPE, AU MOYEN-ORIENT ET A L'AFRIQUE.

Except contrary indications, pages of this catalogue are presented as below :

- left page : modular instruments
- right page : complete non modular instrument

In the majority of cases, modular instruments are linked to complete instruments. If we don't have any correspondence, you will have an «X» on the line.


Example : CEV625C doesn't have complete version.

In the modular inserts table, handles to be used are precised with specific pages. Tubing compatible with handles are mentioned on pages close to handles' pages.

Please refer to modular instruments composition on following page.










MicroFrance® laparoscopic dismantled instruments are composed of three pieces : Jaws, tube and handle..

The black portion of the left tab indicates the section featured in the current page :




In the majority of cases, modular instruments are linked to complete instruments

Handles
































Jaws 















Modular scissors insert : insulated scissors insert blades and wire insert for scissors insert

Tubes 












Non rotatable instrument	Rotatable instrument	Insulated tube	Non insulated tube	Tube with independant sheath	Used for bipolar instrument
					
					
					
					
					














Modular				
Inserts				Handle
Jaws Scale 1/1	Working length			Page #
	Short	Standard	Large	
 Johann 20x5mm	CEV8625-1 245 mm	CEV625-1 345 mm	CEV9625-1 445 mm	   Rotatable 40 -> 47
 Johann 20x5mm		CEV625P1 345 mm Atraumatic grip		Non rotatable 48 - 49  
 Double-action jaws, 20x5mm	CEV8625D 245 mm	CEV625D 345 mm	CEV9625D 445 mm	
 Concave 20x5mm	X	CEV670-1 345 mm	X	
 Curved jaws, 20x5mm	CEV8625C 245 mm	CEV625C 345 mm	CEV9625C 445 mm	

NON modular










Complete instruments










Jaws Scale 1/1	Handle		Working length		
			Short	Standard	Large
 Johann 20x5mm	 Monopolar with ratchet	Montupet CEV815M 240 mm	CEV405 365 mm	X	
	 With spring	X	Gayet CEV405BG 365 mm	X	
	 With ratchet	Montupet CEV815 240 mm	X	X	
	 With ratchet	CEV915G 240 mm	CEV525 345 mm	X	
	 Monopolar with ratchet	X	CEV525M 345 mm	X	
	 With long ratchet	X	CEV725 345 mm	X	
X	X	X	X	X	X
 Concave 20x5mm	 Monopolar with ratchet	X	CEV470 365 mm	X	
 Curved jaws, 20x5mm	 With ratchet	X	CEV525C 345 mm	X	

Modular				
Inserts				Handle
Jaws Scale 1/1	Working length			Page #
	Short	Standard	Large	
 15x5mm	X	CEV4061 340 mm	X	
 Demaizières 20x3mm	CEV84071 245 mm	CEV4071 345 mm	CEV94071 445mm	
 40x5mm	CEV8660-1 270 mm	CEV660-1 370 mm	CEV9660-1 470 mm	 Rotatable 40 -> 47
 Concave 40x5mm	X	CEV671-1 370 mm	X	Non rotatable 48 - 49 
 Crovella 40x5mm smooth jaws	X	CEV6611 370 mm	X	 

NON modular






Complete instruments






Jaws Scale 1/1	Handle		Working length		
			Short	Standard	Large
 15x5mm		Monopolar with ratchet		CEV406 360mm	X
 Demaizières 20x3mm		Monopolar with ratchet	X	Demaizières CEV407 365mm	X
 40x5mm		Monopolar with ratchet	X	Desplantez CEV460 385mm	X
		Axial with long ratchet	X	CEV760 365mm	X
 Concave 40x5mm		Monopolar with ratchet	X	CEV470-1 385mm	X

Modular				
Inserts				Handle
Jaws Scale 1/1	Working length			Page #
	Short	Standard	Large	
 Round double-action jaws 15x5mm	X	CEV433D 350 mm	X	   Rotatable 40 -> 47 Non rotatable 48 - 49
 Smooth jaws 20x3mm	CEV8625E 245 mm	CEV625E 345 mm	X	
 Porpiglia one fenestrated jaws 20x5mm	X	CEV625P 345 mm	X	
 Concave double action jaws 15x5mm	CEV8670D 240 mm	CEV670D 340 mm	CEV9670D 440 mm	 











Modular






Inserts		Handle
Jaws Scale 1/1		Page #
 <p>De Manzini CEV7625I 435 mm - 20x5mm</p>	 <p>Scopinaro CEV7660I 435 mm - 40x5mm</p>	
Tubes Scale 1/3		
CEV649M5	CEV649G5	
		 <p>Rotatable 46 - 47</p>
Graduation 50 / 100 / 150 mm	Graduation 100 / 150 / 200 mm	

Modular			
Inserts		Handle	
Jaws Scale 1/1	Working length	Page #	
 Round jaws, double action jaws	CEV433-1 370 mm	 Rotatable 46 -> 47	
 35x10mm	CEV610-1 375 mm		
X	X		
X	X		
 Hüscher Round jaws, double action jaws	CEV433H 380 mm		
 Round jaws, double action jaws with tube CEV864910B	CEV8434I 295 mm Short		

NON modular

Complete instruments








Jaws Scale 1/1	Handle		Working length
 Round jaws, double action jaws		Axial with ratchet	CEV533-2 370 mm
 35x10mm		Angled, with long ratchet	CEV410 360 mm
		Axial with ratchet	CEV530 370 mm
 Duval		Angled, with long ratchet	CEV432 365 mm
 Round jaws		Angled, with long ratchet	CEV433 365 mm
X	X	X	X
X	X	X	X










Modular				
Inserts			Handle	
Jaws Scale 1/1	Working length			Page #
 ø 10mm				 Delaitre CEV602-1 370 mm Rotatable 46 - 47
 Delaitre ø 5mm	Short	Standard	Large	 Rotatable 40 -> 47 Non rotatable 48 - 49 
	CEV8621-1 250 mm	CEV621-1 350 mm	X	

NON modular










Complete instruments














Jaws Scale 1/1	Handle		Working length		
 ø 10mm	 Angled, with long ratchet		Delaitre CEV420 360 mm		
	 Axial with ratchet		Delaitre CEV502 370 mm		
 ø 5mm	 Angled, with long ratchet		Short	Standard	Large
			Montupet CEV821G 245 mm	Delaitre CEV421 370 mm	X
	 Axial with ratchet		X	Delaitre CEV521 350 mm	X
	 Axial, long ratchet		X	CEV721 350 mm	X

Modular				
Inserts				Handle
Jaws Scale 1/1	Working length			Page #
	Short	Standard	Large	
 Manhès 10mm	X	CEV614-1 335 mm	CEV9614-1 435 mm	   Rotatable 40 -> 47 Non rotatable 48 - 49  
 Manhès 15mm	CEV8615-1 240 mm	CEV615-1 340 mm	CEV9615-1 440 mm	
 Manhès 20mm	X	CEV620-1 345 mm	CEV9620-1 445 mm	
 Dauphin 20x3mm	X	CEV617l 340 mm	CEV9617l 440 mm	














NON modular
Complete instruments





Jaws Scale 1/1	Handle		Working length		
			Short	Standard	Large
 10mm		Monopolar with ratchet	X	Manhès CM110 330 mm	X
 15mm		Monopolar with ratchet	X	Manhès CM115 360 mm	X
		Axial with ratchet	X	CEV515 340 mm	X
		Axial, monopolar with ratchet	X	CEV515M 340 mm	X
 20mm		Monopolar, with spring	X	Manhès CEV220 365 mm	X
		Monopolar with ratchet	X	Manhès CM120 365 mm	X
X	X	X	X	X	X

Modular				
Inserts			Handle	
Jaws Scale 1/1	Working length			Page #
 ø 10mm Grasping insert	Museux CEV606-1 370 mm			 Rotatable 46 -> 47
 Manhès Grasping insert ø 5mm	Short	Standard	Large	   Rotatable 40 -> 47
	X	Manhès CEV611-1 335 mm	CEV9611-1 435 mm	
 Manhès Cobra ø 5mm	X	Manhès CEV612-1 340 mm	X	Non rotatable 48 - 49  
 Manhès Tigre ø 5mm	X	Manhès CEV613-1 340 mm	X	









NON modular








Complete instruments

Jaws Scale 1/1	Handle		Working length		
			Short	Standard	Large
 ø 10mm	 Angled with spring	X	Museux CEV321 360 mm	X	
	 Axial with ratchet	X	Museux CEV506 370 mm	X	
 Manhès Grasping forceps with teeth ø 5mm	 Monopolar, with spring	X	CEV211 355 mm	X	
	 Monopolar with ratchet	X	Manhès CM111 355 mm	X	
	 Axial with ratchet	X	Manhès CEV511 335 mm	X	
	 Monopolar with ratchet	X	Manhès CEV511M 335 mm	X	
 Manhès Cobra ø 5mm	 Monopolar with ratchet	X	Manhès CM112 355 mm	X	
	 Axial with ratchet	X	Manhès CEV512 335 mm	X	
 Manhès Tigre ø 5mm	 Monopolar with ratchet	X	Manhès CM113 360 mm	X	

Modular				
Inserts				Handle
Jaws Scale 1/1	Working length			Page #
	Short	Standard	Large	
 Cobra Double action jaws ø 5mm	CEV8612D 240 mm	CEV612D 340 mm	CEV9612D 440 mm	 Rotatable 40 -> 47
 Vallancien Fenestrated jaws insert with teeth ø 5mm	CEV8616I 245 mm	CEV616I 345 mm	X	Non rotatable 48 - 49 









Modular				
Inserts				Handle
Jaws Scale 1/1	Working length			Page #
	Short	Standard	Large	
 Delaitre Double action jaws, 22mm ø 5mm	X	CEV641D 355 mm	X	  Rotatable 40 -> 47 Non rotatable 48 - 49  
 Double action jaws, 20mm ø 5mm	X	CEV641F 350 mm	X	
 Allal Double action jaws 15mm ø 5mm	X	CEV642TF 345 mm	X	
 Double action jaws, angled 20mm jaws 90deg ø 5mm	X	CEV643I 345 mm	X	







Modular				
Inserts				Handle
Jaws Scale 1/1	Working length			Page #
	Short	Standard	Large	
 Short curved Double action jaws 12mm	X	CEV532I 335 mm	X	  Rotatable 40 -> 47 Non rotatable 48 - 49  
 Medium curved Double action jaws 15mm	CEV8641-1 245 mm	CEV641-1 345 mm	CEV9641-1 445 mm	
 Long curved Double action jaws 18mm	CEV8642-1 245 mm	CEV642-1 345 mm	X	

NON modular

Complete instruments


















Jaws Scale 1/1	Handle		Working length		
			Short	Standard	Large
 Short curved Double action jaws 12mm		Axial with ratchet	X	CEV532 340 mm	X
 Medium curved Double action jaws 15mm		Axial with ratchet	X	CEV531-2 340 mm	X
 Long curved Double action jaws 18mm		Axial, monopolar with ratchet	X	CEV531-3 350 mm	X

	Modular		
	Inserts		Handle
	Jaws Scale 1/1	Working length	Page #
 <p>Delaitre Medium curved 50mm</p>	<p>CEV601-1 380 mm</p>		
 <p>Delaitre Medium curved - Double action jaws 50mm</p>	<p>CEV651-1 380 mm</p>		
 <p>Pignata Short curved 90°, 45mm</p>	<p>CEV601P 380 mm</p>	<p>Rotatable 46 -> 47</p>	
 <p>Puglisi Long curved Double action jaws 60mm</p>	<p>CEV652I 380 mm</p>		

NON modular

















Complete instruments







	NON modular			
	Complete instruments			
Jaws Scale 1/1	Handle		Working length	
 Medium curved 50mm		Angled	Delaitre CEV430 370 mm	
		Axial with ratchet	Nouaille CEV501 380 mm	
 Delaitre Medium curved - Double action jaws 50mm		Axial with ratchet	CEV501-2 380 mm	
X	X	X	X	
X	X	X	X	

Modular				
Inserts				Handle
Jaws Scale 1/1	Working length			Page #
	Short	Standard	Large	
 straight tungsten jaws	X	CEV500T 335 mm	CEV9500T 435 mm	   Rotatable 40 -> 47 Non rotatable 48 - 49  
 tungsten curved left jaws	CEV8538T 240 mm	CEV538T 340 mm	X	
 tungsten right curved jaws	CEV8544I 235 mm	CEV544I 335 mm	CEV9544I 435 mm	
 titanium automatic jaws	CEV8650T 235 mm	CEV650T 335 mm	CEV9650T 435 mm	

NON modular





Complete instruments

Jaws Scale 1/1	Handle		Working length		
			Short	Standard	Large
 straight tungsten jaws	 Axial with ratchet	CEV965 230 mm	CEV500T5 335 mm	X	
	 Large ergonomic handle with ratchet up	X	CEV700T5 335 mm	X	
	 Large ergonomic handle with ratchet down	X	CEV700T5D 335 mm	CEV9700T5D 450 mm	
 tungsten curved left jaws	 Axial with ratchet	X	CEV538T5 335 mm	X	
	 Large ergonomic handle with ratchet down	X	CEV738T5 335 mm	CEV9738T5 455 mm	
	 Fine ergonomic handle	X	CEV738-5 350 mm	X	
 tungsten right curved jaws	 Axial with ratchet	X	CEV544 335 mm	X	
 titanium automatic jaws	 Axial with ratchet	X	Mouiel CEV540-5 335 mm	X	
	 Large ergonomic handle with ratchet down		CEV650T5D 335 mm	X	
	 Fine ergonomic handle	X	CEV740-5 350 mm	X	
 tungsten right curved jaws	 Large ergonomic handle with ratchet down		CEV744T5 335 mm	CEV9744T5 455 mm	








NON modular			
Complete instruments			
Jaws Scale 1/1	Handle		Working length
 Nouaille straight jaws		Axial with ratchet	Nouaille CEV500-A10 345 mm
 left curved jaws		Axial with ratchet	CEV538-10 345 mm
 Mouiel straight automatic jaws		Axial with ratchet	Mouiel CEV540-10 345 mm













Knot Pusher - 5 mm

NON modular	
Complete instruments	
 <p>Delaitre CEV442 340 mm</p>	 <p style="text-align: right;">Scale 1/3</p>
<p>CEV442E 365 mm</p>	
<p>CEV442Z 350 mm</p>	












Knot guide scissors - 5 mm

NON modular					
Complete instruments					
Jaws Scale 1/1	Handle		Working length		
			Short	Standard	Large
 <p>Drouard Complete non modular instrument</p>		<p>Angled</p>	<p>X</p>	<p>CEV445-5 355 mm</p>	<p>X</p>
Modular					
 <p>Drouard Instrument démontable</p>	  <p>Rotatable 40 -> 47</p>	  <p>Non rotatable 48 - 49</p>	<p>CEV8445I 235 mm</p>	<p>CEV445I 335 mm</p>	<p>CEV9445I 435 mm</p>

Modular				
Inserts				Handle
Jaws Scale 1/1	Working length			Page #
	Short	Standard	Large	
 Straight, serrated blades 10mm	CEV8104I 240 mm	CEV104I 340 mm	CEV9104I 440 mm	
 Straight, Double action jaws 10mm	X	CEV623-1 340 mm	X	
 Curved, Double action jaws 15mm	CEV8605-1 245 mm	CEV605-1 345 mm Pack 1 CEV605-3 345 mm Pack 3	CEV9605-1 445 mm	 Rotatable 40 -> 47
 Curved, insulated, double action jaws 15mm	CEV8605G 250 mm Pack 3	CEV605G 350 mm Pack 3	CEV9605G 450 mm Pack 3	Non rotatable 48 - 49 
 Curved, insulated, double action titanium jaws 15mm	CEV8605GT 250 mm	CEV605GT 350 mm	CEV9605GT 450 mm	
 Curved, Double action jaws 10mm	CEV8604-1 240 mm	CEV604-1 340 mm Pack 1 CEV604-3 340 mm Pack 3	X	













Modular

Inserts				Handle
Jaws Scale 1/1	Wire insert			Page #
	Working length			
	Short	Standard	Large	
	CEV8607-1 245 mm	CEV607-1 345 mm	CEV9607-1 445 mm	
 CEV607 Curved insulated scissors insert blades - Pack 6 Double action jaws 15 mm for wire insert				  Rotatable 40 -> 47 Non rotatable 48 - 49
 CEV607C Curved insulated scissors insert blades - Pack 6 Double action jaws 10 mm for wire insert				
 CEV607D Straight insulated scissors insert blades - Pack 6 Double action jaws 10 mm for wire insert				
 CEV605BGP (Pack 1) CEV605GP (Pack 3)				 
				
Metal tube for insert CEV605BGP / CEV605GP and black sheath	CEV8649GP Pack 1 CEV8649BGP Pack 3	CEV649GP Pack 1 CEV649BGP Pack 3	CEV9649GP Pack 1 CEV9649BGP Pack 3	
				
Black sheath for metal tube - Pack 3	CEV8649C5 210mm	CEV649C5 310mm	CEV9649C5 410mm	




Monopolar cables - 4,5 m

		
CP391E Erbe Storz Wolf	CP391M Aesculap Berchtold Martin Erbe T Serie Lamidey	CP391V Bowa Emc Lamidey Valleylab



NON modular					
Complete instruments					
Jaws Scale 1/1	Handle		Working length		
			Short	Standard	Large
  Straight scissors serrated blades 10mm	 Angled, monopolar	X	CEV104M 355 mm	X	
	 Angled	Montupet CEV804 235 mm	Manhès CM104 355 mm	X	
	 Axial	CEV904 235 mm	X	X	
 Straight scissors Delicate serrated blades 10mm	 Angled	X	Manhès CM106 355 mm	X	
 Curved scissors serrated blades 10mm	 Angled monopolaire	X	CEV114M 355 mm	X	
	 Axial, monopolar	X	CEV514M 335 mm	X	





NON modular

Complete instruments

Jaws Scale 1/1	Handle		Working length
 <p>Curved scissors serrated blades 25 mm</p>		Angled	CEV404 340 mm
		Axial	CEV503 350 mm



Modular

 <p>Curved scissors serrated blades Double action jaws 25 mm</p>	 <p>Rotatable 46 -> 47</p>	Working length



Modular				
Inserts				Handle
Jaws Scale 1/1	Working length			Page #
	Short	Standard	Large	
 Straight clamp insert 25mm jaws ϕ 5mm	CEV8626I 250 mm	CEV626-1 350 mm	X	 Rotable 40 -> 47
 Curved clamp insert 25mm jaws Double action jaws ϕ 5mm	CEV8627I 255 mm	CEV627I 355 mm	CEV9627I 455 mm	Non rotatable 48 - 49 






Vascular clamp - ø 5 mm

NON modular					
Complete instruments					
Jaws Scale 1/1	Handle		Working length		
			Short	Standard	Large
 <p>Debakey Straight clamp 20mm jaws ø 5mm</p>		Axial with ratchet	X	Debakey CEV526 350 mm	X

Vascular clamp - ø 12 mm










NON modular			
Complete instruments			
Jaws Scale 1/1	Handle		Working length
 <p>Dulucq Left curved 60mm ø 12mm</p>		Axial with ratchet	Dulucq CEV570 400 mm




	Modular			
	Inserts			Handle
	Jaws Scale 1/1	Working length		Page #
	Short	Standard		
X	X	X		
 Delaitre Left curved 45mm \varnothing 10mm	X	CEV635-1 385 mm	 Rotatable 46 -> 47	
X	X	X		
X	X	X		

NON modular

Complete instruments
















Jaws Scale 1/1	Handle		Working length
 <p>Delaitre Left curved 45mm \varnothing 10mm</p>		Axial with ratchet with safety lock	Delaitre CEV505 385 mm
 <p>Querleu Left curved 45mm \varnothing 10mm</p>		Angled with long ratchet	Querleu CEV435-1 375 mm
		Axial with ratchet with safety lock	Querleu CEV535 385 mm
 <p>Querleu Right curved 45mm \varnothing 10mm</p>		Angled with long ratchet	Querleu CEV435-2 375 mm
 <p>Curved 75mm \varnothing 10mm</p>		Axial with ratchet with safety lock	CEVA112 320 mm

Modular		
Inserts		Handle
Jaws Scale 1/1	Working length	Page #
 <p style="text-align: center;">Straight intestinal clamp insert 65x7mm jaws - ø 10mm</p>	<p>CEV6o8-1 395 mm</p>	 <p style="text-align: center;">Rotatable 46 - 47</p>
 <p style="text-align: center;">Annibale Straight intestinal clamp insert Double action jaws 65x7mm jaws - ø 10mm</p>	<p>Annibale CEV6o8D 395 mm</p>	

NON modular

Complete instruments









	Jaws Scale 1/1	Handle		Working length
	 <p data-bbox="443 952 654 1003">Straight intestinal clamp 65x7mm jaws - ø 10mm</p>		Angled with long ratchet	Delaitre CEV431-65 390 mm
<p data-bbox="539 1550 558 1576">X</p>		Axial with ratchet with safety lock	CEV508 400 mm	
	<p data-bbox="981 1550 1000 1576">X</p>	<p data-bbox="1257 1550 1276 1576">X</p>	<p data-bbox="1409 1550 1428 1576">X</p>	

NON modular			
Complete instruments			
Jaws Scale 1/1		Handle	Working length
 Satinsky Angled clamp curved tube 110mm jaws		 Axial with ratchet with safety lock	CEVA101 260 mm
 Angled clamp curved tube 90mm jaws			CEVA102 250 mm
 Satinsky Angled clamp curved tube 80mm jaws			CEVA103 240 mm
 Curved clamp curved tube 110mm jaws			CEVA104 290 mm








NON modular


Complete instruments




	Jaws Scale 1/1		Handle	Working length
	 <p>Demaizières Curved clamp 85mm jaws curved tube 85mm jaws</p>			<p>CEVA110 270 mm</p>
 <p>Demaizières Curved clamp curved tube 90mm jaws</p>		 <p>Axial with ratchet with safety lock</p>	<p>CEVA528 260 mm</p>	
 <p>Demaizières Curved clamp curved tube 60mm jaws</p>			<p>CEVA529 220 mm</p>	

Temporary Clamps



 <p>CEV565 + CEV560</p>		
<p>Jaws Scale 1/1</p>		
	<p>Straight vascular clamp 15x3mm Pack 2</p>	<p>Sartori CEV568</p>
	<p>Straight vascular clamp 25x3mm Pack 2</p>	<p>Begin CEV565</p>
	<p>Choledoque clamp 35x4mm Pack 2</p>	<p>Begin CEV566</p>
	<p>Intestinal clamp 45x5mm Pack 2</p>	<p>Begin CEV567</p>

Clamp applicator - 10 mm diameter






NON modular			
	Jaws Scale 1/1	Handle	Working length
	 <p>Clamp applicator ø 10mm 360mm</p>	 <p>Axial with ratchet</p>	<p>Begin CEV560 360 mm</p>

Modular		
	 <p>Clamp applicator insert for tube ø 10mm 310mm</p>	 <p>Rotatable 46 -> 47</p>
		<p>CEV560-1 355 mm</p>

Clips applier









NON modular			
Complete instruments			
Jaws Scale 1/1	Handle		Working length
 <p>Clip Applier for T1300 et LC300 clips Ø 10mm</p>		<p>Angled with spring</p>	<p>CEV310 335 mm</p>

Circular suturing holder insert





Modular			
 <p>Morino Circular suturing holder insert Double action jaws Ø 10mm</p>	 <p>Rotatable 46 -> 47</p>	<p>Morino CEV603I 375 mm</p>	
 <p>Massard Circular suturing holder insert Ø 5mm</p>	 <p>Rotatable 40 -> 47 Non rotatable 48 - 49</p> 	<p>Massard CEV622I 335 mm</p>	

Modular



Rotatable instruments

Handle Scale 1/3	Reference	Handle Scale 1/3	Reference
 <p>Monopolar With ratchet Angled handle With black rotatable ring</p>	CEV1019-5B	 <p>Monopolar With ratchet Angled handle With black rotatable ring</p>	CEV10195C
 <p>Monopolar With ratchet Angled handle With black rotatable ring</p>	CEV10195D	 <p>Monopolar With ratchet Angled handle With black rotatable ring</p>	CEV10195R
 <p>Monopolar Without ratchet Angled handle With black rotatable ring</p>	CEV1039-5-B	 <p>Monopolar Without ratchet Angled handle With black rotatable ring</p>	CEV10395D
 <p>Monopolar Without ratchet Angled handle With black rotatable ring</p>	CEV10395R	 <p>Monopolar Without ratchet Angled handle With black rotatable ring</p>	CEV1039G5

Insulated tubes - 5 mm

Modular			
Rotatable instruments			
Working length			
Short	Standard	Large	
			
CEV8649-5-B 210mm	CEV649-5B 310mm	CEV9649-5-B 410mm	
			
CEV8646B5 210mm	CEV646B5 310mm	CEV9646B5 410mm	
			
CEV8647B5 210mm	CEV647B5 310mm	CEV9647B5 410mm	
			
CEV8648B5 210mm	CEV648B5 310mm	CEV9648B5 410mm	

Metal tube and Pack of 3 black sheaths

			
CEV8649B5 210mm	CEV649B5 310mm	CEV9649B5 410mm	
			
CEV8649C5 210mm	CEV649C5 310mm	CEV9649C5 410mm	

Monopolar cables 4,5 m








CP391E
Erbe Storz Wolf







CP391M
Aesculap Berchtold Martin
Erbe T Serie Lamidey





CP391V
Bowa Emc Lamidey Valleylab

Modular			
Rotatable instruments			
Handle Scale 1/3	Reference	Handle Scale 1/3	Reference
 <p>3 teeth</p> <p>Monopolar With ratchet With black rotatable ring</p>	CEV6185R	 <p>With ratchet With black rotatable ring</p>	CEV6285R
 <p>Monopolar Without ratchet With black rotatable ring</p>	CEV6385R	 <p>With ratchet With black rotatable ring</p>	CEV728T5
 <p>5 teeth</p> <p>With ratchet With black rotatable ring</p>	CEV628P5		

Insulated tubes, 5 mm diameter

Modular		
Rotatable instruments		
Working length		
Short	Standard	Large
		
Black insulated tube		
CEV8649-5-B 210mm	CEV649-5B 310mm	CEV9649-5-B 410mm
		
Yellow insulated tube		
CEV8646B5 210mm	CEV646B5 310mm	CEV9646B5 410mm
		
Blue insulated tube		
CEV8647B5 210mm	CEV647B5 310mm	CEV9647B5 410mm
		
Green insulated tube		
CEV8648B5 210mm	CEV648B5 310mm	CEV9648B5 410mm

Metal tube and Pack of 3 black sheaths

		
Metal tube 5mm diam. WL 310mm for black sheath		
CEV8649B5 210mm	CEV649B5 310mm	CEV9649B5 410mm
		
Black sheath 5mm diam. WL 310mm for metal tube Pack 3		
CEV8649C5 210mm	CEV649C5 310mm	CEV9649C5 410mm

Monopolar cables 4,5 m







CP391E
Erbe Storz Wolf







CP391M
Aesculap Berchtold Martin
Erbe T Serie Lamidey





CP391V
Bowa Emc Lamidey Valleylab

Modular			
Rotatable instruments			
Handle Scale 1/3	Reference	Handle Scale 1/3	Reference
 Large ergonomic handle Monopolar With ratchet With black rotatable ring	CEV720BR	 Fine ergonomic handle Monopolar With ratchet With black rotatable ring	CEV720R
 Elhadad Fine ergonomic handle Monopolar With ratchet With black rotatable ring	CEV7285R	 Large ergonomic handle Monopolar With ratchet With black rotatable ring	CEV728BR5

Insulated tubes - 5 mm

Modular			
Rotatable instruments			
Working length			
Short	Standard	Large	
			
Black insulated tube			
CEV8649-5-B 210mm	CEV649-5B 310mm	CEV9649-5-B 410mm	
			
Yellow insulated tube			
CEV8646B5 210mm	CEV646B5 310mm	CEV9646B5 410mm	
			
Blue insulated tube			
CEV8647B5 210mm	CEV647B5 310mm	CEV9647B5 410mm	
			
Green insulated tube			
CEV8648B5 210mm	CEV648B5 310mm	CEV9648B5 410mm	

Metal tube and Pack of 3 black sheaths

			
Metal tube 5mm diam. WL 310mm for black sheath			
CEV8649B5 210mm	CEV649B5 310mm	CEV9649B5 410mm	
			
Black sheath 5mm diam. WL 310mm for metal tube Pack 3			
CEV8649C5 210mm	CEV649C5 310mm	CEV9649C5 410mm	

Monopolar cables 4,5 m










CP391E
Erbe Storz Wolf



CP391M
Aesculap Berchtold Martin
Erbe T Serie Lamidey



CP391V
Bowa Emc Lamidey Valleylab

Modular			
Rotatable instruments			
Handle Scale 1/3	Reference	Handle Scale 1/3	Reference
 Angled handle for 5 and 10mm diam. non-insulated tube With long ratchet with grey rotatable ring	CEV102910R	 Angled handle for 5 and 10mm diam. non-insulated tube Without ratchet with grey rotatable ring	CEV1059-10B
 Axial handle for 5 and 10mm diam. non-insulated tube With ratchet with grey rotatable ring	CEV628-10	 Axial handle for 5 and 10mm diam. non-insulated tube Without ratchet with grey rotatable ring	CEV658-10
 Fine ergonomic handle for 5 and 10mm diam. non-insulated tube With ratchet with grey rotatable ring	CEV72810R	 Large ergonomic handle for 5 and 10mm diam. non-insulated tube With ratchet with grey rotatable ring	CEV728BR
 Axial handle for 5 and 10mm diam. non-insulated tube With long ratchet with grey rotatable ring	CEV728T10	Axial handle for 5 and 10mm diam. non-insulated tube With long ratchet with grey rotatable ring	





Modular

Rotatable instruments

Working length

Short

Standard

Large



Non-insulated tube for rotatable instrument \varnothing 10mm

CEV864910B
215mm

CEV649-10B
305mm

X









Non-insulated tube for rotatable instrument \varnothing 5mm

CEV8649X5
210mm


CEV649X5
310mm

CEV9649X5
410mm

Non Rotatable handles - 5 mm

Modular			
Non rotatable instruments			
Handle Scale 1/3	Reference	Handle Scale 1/3	Reference
 Angled handle for non-rotatable insulated tube Monopolar - With long ratchet	CEV1019-5N	 Axial handle for non-rotatable insulated tube Monopolar - With long ratchet	CEV618-5N
 Angled handle for non-rotatable insulated tube Monopolar - With long ratchet	CEV10195NA	 Axial handle for non-rotatable insulated tube With ratchet	CEV628-5N
 Angled handle for non-rotatable insulated tube Monopolar - Without ratchet	CEV1039-5N	 Axial handle for non-rotatable insulated tube Monopolar - Without ratchet	CEV638-5N

Black insulated tube for non-rotatable instrument ø 5mm



 Black insulated tube for non-rotatable instrument ø 5mm		
Working length		
Short	Standard	Large
CEV8649-5N 210mm	CEV649-5N 310mm	CEV9649-5N 410mm



Non Rotatable handles - 5 mm











Modular

Non rotatable instruments

Handle Scale 1/3	Reference	Handle Scale 1/3	Reference
 <p>Poignée ergonomique fine pour tube métallique ø 5mm With long ratchet</p>	CEV720	 <p>Poignée ergonomique grip large pour tube métallique ø 5mm With long ratchet</p>	CEV720B

Non-insulated tube for non-rotatable instrument ø 5mm





Working length		
Short	Standard	Large
CEV8649P5 210mm	CEV649P5 310mm	CEV9649P5 410mm

Modular					
Inserts					Handle
Jaws Scale 1/1		Working length			
		Short	Standard	Large	
	Manhès Grasping forceps	CEV8611H 235 mm	CEV611H 335 mm	CEV9611H 435 mm	 CEV669E Large ergonomic handle bipolar Without ratchet with blue rotatable ring for insulated tube ø 5mm
	Manhès Cobra	CEV8612LH 245 mm	CEV612LH 345 mm	CEV9612LH 445 mm	
	Alligator	CEV8615H 240 mm	CEV615H 340 mm	CEV9615H 440 mm	
	Johann 20x5mm	CEV8625H 245 mm	CEV625H 345 mm	CEV9625H 445 mm	
	Botella Triangle fenestrated jaws	CEV8633H 245 mm	CEV633H 345 mm	CEV9633H 445 mm	
	Dissector curved left	CEV8642H 245 mm	CEV642H 345 mm	CEV9642H 445 mm	
	Dissector curved right	CEV8642HD 245 mm	CEV642HD 345 mm	CEV9642HD 445 mm	





CEV600H
Protection part
for tip jaws Bipolar
Delivered with bipolar insert
For cleaning and sterilizing

Insulated tubes - 5 mm

Modular		
Rotatable instruments		
Working length		
Short	Standard	Large
		
Black insulated tube		
CEV8649-5-B 210mm	CEV649-5B 310mm	CEV9649-5-B 410mm
		
Yellow insulated tube		
CEV8646B5 210mm	CEV646B5 310mm	CEV9646B5 410mm
		
Blue insulated tube		
CEV8647B5 210mm	CEV647B5 310mm	CEV9647B5 410mm
		
Green insulated tube		
CEV8648B5 210mm	CEV648B5 310mm	CEV9648B5 410mm

Metal tube and Pack of 3 black sheaths

		
Metal tube 5mm diam. WL 310mm for black sheath		
CEV8649B5 210mm	CEV649B5 310mm	CEV9649B5 410mm
		
Black sheath 5mm diam. WL 310mm for metal tube Pack 3		
CEV8649C5 210mm	CEV649C5 310mm	CEV9649C5 410mm

Bipolar cables 4,5 m











CP393-2
Bowa Emc Lamidey Valleylab



CP393-3
Erbe Storz Wolf



CP393-4
Aesculap Berchtold Martin









Modular								Handle Non rotatable
Inserts & Tubes								
Jaws Scale 1/1	Working length							
	Short		Standard		Large			
	Mors	Tube 1/3	Mors	Tube 1/3	Mors	Tube 1/3		
	 Mouiel 15x5mm	CEV86341A 245 mm		CEV634-1A 340 mm		CEV9634-1A 430 mm		 CEV669B Angled
	 with hook 15x5mm	X		CEV632-1A 340 mm		CEV96321A 430 mm		 CEV669BG Gayet Ergonomic
	 Gayet 20x5mm	CEV86361A 250 mm		CEV636-1A 345 mm		CEV96361A 435 mm		 MCEN54P Angled
	 Vallancien with teeth 20x5mm	X		CEV637I 345 mm		X		
			CEV86795B 210 mm		CEV6795B 310 mm		CEV96795B 395 mm	











NON modular

Complete instruments












Jaws Scale 1/1	Handle	Working length		
		Short	Standard	Large
 15x5mm	 Angled	CEV834 240 mm	Mouiel CEV134 360 mm	X
	 Ergonomic	X	Gayet CEV134BG 360 mm	Gayet CEV9134BG 450 mm
	 Axial	X	CEV580-2 365 mm	X
 15x5mm	 Angled	X	CEV132 360 mm	X
 Gayet 20x5mm	 Angled	X	Gayet CEV136 365 mm	X
X	X	X	X	X










Modular								
Inserts & Tubes							Handle Non rotatable	
Jaws Scale 1/1	Working length							
	Short		Standard		Large			
	Mors	Tube 1/3	Mors	Tube 1/3	Mors	Tube 1/3		
 Manhès Spatula jaws 3mm	X		CEV630-1A 340 mm		CEV9630-1A 425 mm		 CEV669B Angled	
 Botella Triangle fenestrated jaws	X		CEV633-1A 340 mm		CEV96331A 430 mm			 Gayet CEV669BG Ergonomic
 Spatula angled jaws 2mm	X		CEV630C 340 mm		X			 MCEN54P Angled
	X	 CEV86795B 210 mm	X	 CEV6795B 310 mm	X		X	
	X		X		X		X	

NON modular





Complete instruments

Jaws Scale 1/2	Handle		Working length		
			Short	Standard	Large
 Spatula jaws 3mm		Angled	Montupet CEV830 240 mm	Manhès CEV130 360 mm	X
 Botella Triangle fenestrated jaws		Angled	X	Botella CEV133 360 mm	X
		Axial	X	Botella CEV580-3 365 mm	X
X	X	X	X	X	X
 Delicate triangular jaws 15x5mm		Angled	X	Botella CEV133-1 360 mm	X
 Manhès 2mm		Angled	X	Manhès CM131 360 mm	X



Modular			
Hook jaws			
Scale 1/1			
			
CEV229-2A	CEV229-3A	CEV229-4A	CEV2295A
Hook # 2	Hook # 3	Hook # 4	Hook # 5
Pack 3	Pack 3	Pack 3	Pack 3
Hook handles			
Scale 1/3			
CEV225A	Working length 430mm		
CEV82291	Working length 230mm		
CEV229-1A	Working length 350mm		
CEV9229-1	Working length 450mm		
CEV229B	Working length 380mm		

Monopolar hook non modular			
Scale 1/3			
	CEV225	Working length 430mm	
	Périssat CEV230-1	Working length 330mm	
	Montupet CEV825	Working length 220mm	


Knife, monopolar electrode

NON modular				
Complete instruments				
1/3				
		CEV135-5	Retractable knife ø 5mm 300mm	
1/3				
1/1		CM109 Manhès	Retractable monopolar electrode ø 5mm 320mm	

Monopolar probe

1/3				
		CEV443 Delaître	Monopolar palpation probe ø 5mm 350mm	

Monopolar cables - 4,5 m










Knife blades










^(D) Distributed by Integra MicroFrance SAS (France)

Palpation probe

	1/3			
	1/1		CEV204 320 mm - Ø 5mm	Round probe
	1/1		CEV204G45 450 mm - Ø 5mm	Round probe, graduated
	1/1		CEV205 340 mm - Ø 5mm	Flap probe
	1/3			
	1/1		CEV207 410 mm - Ø 5mm	Double ended probe round end flap graduated

Articulated retractor

		CEV416 380mm - Ø 5mm	Flexible retractor	1/3
		CEV413 370mm - Ø 5mm	Retractor with 3 blades 35mm	1/3
		CEV425 320mm - Ø 5mm	Nouaille Retractor losenge design	1/3
		CEV423 400mm - Ø 5mm	Retractor with 3 blades 65mm	1/3
		CEV453 400mm Ø 10mm	Retractor with 5 blades 65mm	1/3
		CEV453A 330mm - Ø 10mm	Articulated retractor with 5 blades 65mm	1/3

NON modular

Complete instruments

Flexible tension arm delivered with instruments holder CEV980-3

Accessories for flexible tension arm



Mosnier
500mm : CEV980A
600mm : CEV980B

CEV980-3
ø 5 et 10mm
Mosnier



CEV984
Operating
table clamp
with screw



CEV985
Operating
table clamp



Electric irrigation suction pump and accessories

Electric irrigation suction pump

Suction irrigation tubes
Sterile - Pack 10

Electric handle

CM201

CM217S^(D)

CM201-9



Eyes suction irrigation tubes for handle CM201-9

1/1

«L» Working length

1/3



CM20132H10

ø 10mm L 320mm



CM20128H-10

ø 10mm L 280mm



CM20132H-5

ø 5mm L 320mm



CM20128H-5

ø 5mm L 280mm



CM20140H5

ø 5mm L 400mm



CM20132L-5

ø 5mm L 320mm



Cleaning sleeve for telescope




Cleaning sleeve for telescope 0°
ø 10mm 300mm
Scale 1/3












CEV81

^(D) Distributed by Integra MicroFrance SAS (France)






Suction irrigation tubes




Modular suction irrigation tubes

Eyes suction irrigation tube ø 5mm		Handles for tube ø 5mm
<p>1/1</p> <p>CEV8103D5 ø 5mm 195mm</p> 		<p>1/3</p>  <p>CEV103AH</p>
<p>CEV103-5 ø 5mm 295mm</p> 		
<p>CEV9103D5 ø 5mm 395mm</p> 		
<p>1/1</p> <p>CEV8103G5 ø 5mm 195mm</p> 		<p>1/3</p>  <p>CEV103GH</p>
<p>CEV103G5 ø 5mm 295mm</p> 		
<p>CEV9103G5 ø 5mm 395mm</p> 		
Eyes suction irrigation tube ø 10mm		Handle for tube ø 10mm
<p>1/1</p> <p>CEV103-10 295mm</p> 	1/3	<p>1/3</p>  <p>CEV103AH</p>
<p>1/1</p> <p>CEV103P Prete 410mm</p> 		


Wire-pass

	<p>Articulated wire-pass</p> <p>CEV438 Bertrand</p> <p>460mm - ø 10mm</p>		1/3
	<p>Articulated wire-pass</p> <p>CEV438B</p> <p>LU 335mm - ø 5mm</p>		1/3
	<p>Retractable wire pass</p> <p>CEV438D5 Morino</p> <p>430mm - ø 5mm</p>		1/3
	<p>Retractable wire pass</p> <p>CEV456</p> <p>170mm - ø 2,5mm</p>		1/3

Mintz Set for punction


<p>Sleeve - ø 2,5mm - 240mm</p> <p>CM170-1</p> <p>Sharp trocar</p> <p>CM170-2</p> <p>Smooth trocar</p> <p>CM170-3</p>		1/3
--	--	-----

Needles

<p>Insufflation needle with luer-lock</p> <p>120mm</p> <p>CEV60-120</p> <p>150mm</p> <p>CEV60-150</p>		1/3
<p>Punction needle with luer-lock</p> <p>ø 5mm - 300mm</p> <p>CEV65</p> <p>ø 1,2mm - 40mm</p> <p>CEV24F</p> <p>ø 2,0mm - 40mm</p> <p>CEV24G</p>		1/3

Cholangiography

1/3



1/1



Cholangiography forceps

CEV448

350mm - ø 2,5mm

1/3



Cholangiography needle

CEV440

230mm - int ø 2,5mm

Joint

CEV440P6

Pack 6



1/3



Guide sleeve for cholangiography

CEV455

170mm - ø 4,0mm

Falchero

Joint

CEV103P6

Pack 6



Choledoscopy

1/3



Mauras


Guide sleeve for choledoscopy delivered with one guide CEV449G

CEV449

190mm - ø 8mm

1/3

Spare parts



White guide sleeve for CEV449

CEV449G

400mm

Sealing cap













CEV105P6

Pack 6

















Scale 1/1

				Working length			
				50 mm	80 mm	100 mm	120 mm
Trocart tip		blunt		CEV8145A5	CEV145A5	CEV7145A5	CEV9145A5
		sharp		CEV8145A3	CEV145A3	CEV7145A3	CEV9145A3
Sealing cap	blue		CEV105P6 Pack 6				
Head of trocar	magnetic valve		CEV177A5				
	silicone valve		CEV127T				
Valve for head of trocar	magnetic		CEV179A5 Pack 6				
	silicone		CEV125P6 Pack 6				
Sleeve	with stopcock	Screw		CEV8145A6	CEV145A6	CEV7145A6	CEV9145A6
		Smooth		CEV8145A7	CEV145A7	CEV7145A7	CEV9145A7
	without stopcock	Screw		CEV8145B6	CEV145B6	CEV7145B6	CEV9145B6
		Smooth		CEV8145B7	CEV145B7	CEV7145B7	CEV9145B7





















Scale 1/1

				Working length	
				90 mm	130 mm
Trocart tip		blunt		CEV150-4	CEV9150A4
		sharp		CEV150-3	CEV9150A3
Sealing cap		red		CEV110P6 Pack 6	
Head of trocar		magnetic valve		CEV177A10	
		silicone valve		CEV127T10	
Valve for head of trocar		magnetic		CEV179A10 Pack 6	
		silicone		CEV125Q6 Pack 6	
Sleeve	with stopcock	Screw		CEV150A6	CEV9150A6
		Smooth		CEV150A7	CEV9150A7
	without stopcock	Screw		CEV150B6	CEV9150B6
		Smooth		CEV150B7	CEV9150B7



Scale 1/1

				Working length	
				90 mm	130 mm
Trocart tip		blunt		CEV172A4	CEV9172A4
		sharp		CEV172A3	CEV9172A3
Sealing cap		Orange		CEV112P6 Pack 6	
Head of trocar		magnetic valve		CEV177A12	
		silicone valve		CEV127T12	
Valve for head of trocar		magnetic		CEV179A12 Pack 6	
		silicone		CEV125Q6 Pack 6	
Sleeve	with stopcock	Screw		CEV172A6	CEV9172A6
		Smooth		CEV172A7	CEV9172A7
	without stopcock	Screw		CEV172B6	CEV9172B6
		Smooth		CEV172B7	CEV9172B7

1/3	CEV105	CEV105I	CEV110I	CEV122I	CEV112I
					
	∅ 5,5mm 80mm	∅ 5,5mm 80mm	∅ 10,5mm 80mm	∅ 12,5mm 80mm	∅ 10,5mm 120mm
	without stopcock	with stopcock	with stopcock	with stopcock	with stopcock
1/3	CEV805I				
					
	Montupet				
	∅ 5,5mm 50mm				
	with stopcock				

Sealing cap for trumpet valve trocar



CEV105P6
∅ 5,5mm
blue
Pack 6



CEV110P6
∅ 10,5mm
red
Pack 6



CEV112P6
∅ 12mm
orange
Pack 6



				Diameter		
				ø5,5mm.	ø 10mm	ø12,5mm
Trocart tip		blunt		CEV145-4	CEV150-4	-
		sharp		CEV145-3	CEV150-3	CEV172-3
Sealing cap				 CEV105P6 Pack 6	 CEV110P6 Pack 6	 CEV112P6 Pack 6
Head of trocar	with stopcock			 CEV178-5	 CEV278-10	 CEV278-12
	without stopcock			 CEV177-5	-	 CEV177-12
Sleeve	without stopcock	Screw	 CEV145-6	 CEV150-6	 CEV172-1	
		Smooth	 CEV145-7	 CEV150-7	-	

Sealing cap



CEV103P6

ø 3,5mm

Pack 6



CEV105P6

ø 5,5mm

Pack 6



CEV110P6

ø 10,5mm

Pack 6



CEV112P6

ø 12mm

Pack 6

Reducer sleeve for Trocar

Silicone ring for
attaching 15mm
reducer sleeve

CEV155B

Pack 6

1/1



CEV155A

from ø 10,5 to 5,5mm



CEV156A

from ø 12/13 to 5,5mm



CEV181A

from ø 12/13 to 10,5mm

1/3



CEV153

135 mm

from ø 5,5 to
3,5mm

CEV155

135 mm

from ø 10,5 to
5,5mm

CEV155-1

65 mm

from ø 10,5 to
5,5mm

CEV155P

135 mm

from ø 12/13 to
5,5mmInternal
tube
ø 5mm

CEV156

135 mm

from ø 12/13 to
10,5mm

CEV181

135 mm

from ø 12/13 to
10,5mm

CEV181-1

65 mm

from ø 12/13 to
10,5mm

Double retractor

1/3



Descottes

CEV417

200 mm

Paire

Double retractor blunt



CEV418

200 mm

Paire

Double retractor blunt



Mosnier










CEV419

160 mm

Paire







Double retractor blunt

Trocars 85mm



Scale 1/3	Trocart tip Blunt	Trocart tip Sharp	Flexible threaded cannula
ø 5,5mm	CET105-1 	CET105-2 	CET105-3 
	CET110-1 	CET110-2 	CET110-3 
ø 12,5mm	CET112-1 	CET112T2 	CET112-3 

Sealing cap : page 69

Forceps, non modular instruments

Jaws Scale 1/1	Handle	Working length		
		Short	Standard	Large
 Manhès Curette forceps 12x5mm ø 5mm	 Angled, monopolar	X	CM107 355mm	X
 Curette forceps 12x5mm, double action jaws ø 5mm	 Angled	X	CEV444 355mm	X
 Fenestrated forceps 20x5mm, curved tube ø 5mm	 Axial with ratchet Curved tube	X	CET525 340mm	X

Modular forceps

 Manhès Curette insert, 12x5mm ø 5mm	 Rotatable Page 40 - 47 Non rotatable Page 48 - 49	CEV8107I 245mm	CEV107I 345mm	CEV9107I 445mm
---	--	--------------------------	-------------------------	--------------------------













1/2



Belghiti	Soubrane		Dousset	Gayet	Partenski
215mm	215mm	215mm	290mm	220mm	190mm
CEV392D20	CEV392S20	CEV392P20	CEV392D30	CEV392G24	CEV392R15
Straight	Straight	Straight	Straight	Straight	Bayonet
Jaws 2mm	Jaws 2mm double irrigation	Serrated jaws 2mm	Jaws 3mm double irrigation	Jaws 5x15mm	Jaws 4x10mm












Semi-Modular










			Handle
Inserts with tubes			Page #
	Jaws with insulated tube Scale 1/1	Working length	
	 Babcock insert	CEV391 335 mm	 Rotatable 40 -> 44
	 Curved scissors insert serrated blades 8mm	CEV393 330 mm	
	 Straight scissors insert serrated blades 8mm	CEV394 330 mm	
	 Fenestrated insert 12x3.5mm jaws	CEV395 335 mm	
	 Alligator jaws insert 6mm	CEV396 330 mm	
	 Dissector insert medium curved 10mm jaws	CEV398 340 mm	
	-	-	
	-	-	
		CEV39L Adapter cleaning luer-lock for tube insert	

NON modular

Complete instruments

Jaws Scale 1/1	Handle		Working length
 Babcock forceps		Angled with long ratchet, monopolar	CEV391B 355 mm
 Straight scissors serrated blades 8mm		Angled monopolar	CEV394B 350 mm
 Fenestrated forceps 12x3mm jaws		Angled with long ratchet, monopolar	CEV395B 355 mm
 Needle-holder straight tungsten jaws		Axial	CEV992B 340 mm
 Needle-holder left curved tungsten jaws			CEV992C 340 mm



NON modular		
Complete instruments : bipolar forceps		
Jaws Scale 1/1	Handle	Working length
 Gayet Forceps fenestrated jaws graduated tube 100mm	 Ergonomic bipolar	CEV397BG 305mm
 1/3		
 Forceps 2mm spatula jaws	 Ergonomic bipolar	CEV397C 305mm
 1/3		
 Forceps 2mm serrated jaws	 Angled bipolar	CEV397 305mm
 1/3		

Cadière forceps

Forceps without port
Fenestrated forceps with sharp jaws

CEV350
1/1 ø 3,0mm 195mm

Cadière



1/3

Palpation probe

Round probe ø 3,5mm 330mm

CEV203
1/1




1/3

Monopolar hook

Monopolar hook ø3,5mm 330mm

CEV390B Round tip
1/1

CEV390C 90°
1/1




1/3

Suction irrigation tubes

Eyes suction irrigation tube, graduated 50/100/150mm

CEV103G3 Eyes suction irrigation tube, blunt graduated 50/100/150mm for handle CEV103GH
ø 3,5mm 295mm

CEV103GH Handle for suction tube




1/3

Knot pusher


CEV342 lanelli
440mm L 340mm ø 3,5mm

CEV342E lanelli
360mm L 340mm ø 3,5mm













1/3


Trocars




Scale 1/1

				Working length	
				50 mm	80 mm
Trocart tip		blunt		CEV8143A4	CEV143A4
		sharp		CEV8143A3	CEV143A3
Sealing cap		green		CEV103P6 Pack 6	
Head of trocar		magnetic valve		CEV177A3	
		silicone valve		CEV127T3	
Valve for head of trocar		magnetic		CEV179A5 Pack 6	
		silicone		CEV125P6 Pack 6	
Sleeve		with stopcock		CEV8143A7	CEV143A7
		without stopcock		CEV8143B7	CEV143B7


Monopolar hook



Retractable knife













Suction irrigation tubes









Discision forceps



Semi-Modular			
Inserts with tubes			Handle
	Jaws with insulated tube Scale 1/1	Working length	Page #
	 Babcock insert	CEV891I 235 mm	 Rotatable 40 -> 44
	 Curved scissors insert serrated blades 8mm	CEV893I 230 mm	
	 Straight scissors insert serrated blades 8mm	CEV894I 230 mm	
	 Fenestrated insert 12x3.5mm jaws	CEV895I 235 mm	
	 Alligator jaws insert 6mm	CEV896I 235 mm	
	 Dissector insert medium curved 10mm jaws	CEV898I 230 mm	
		CEV39L Adapter cleaning luer-lock for tube insert	

NON modular

Complete instruments







Jaws Scale 1/1	Handle		Working length
 Montupet Babcock forceps		Angled with long ratchet, monopolar	CEV891M 200 mm
 Montupet Straight scissors serrated blades 8mm		Angled monopolar	CEV894M 220 mm
 Montupet Fenestrated forceps 12mm jaws 12x3,5mm		Angled with long ratchet, monopolar	CEV895M 195 mm





NON modular

Complete instruments : bipolar forceps

Jaws Scale 1/1	Handle		Working length
 Bipolar forceps 1.5mm serrated jaws		Angled bipolar	CEV897 250 mm
 1/3			
 Bipolar forceps spatula serrated jaws		Ergonomic bipolar	CEV897C 250 mm
 1/3			




















NON modular

Complete instruments

Monopolar instruments
Needle-holder

Jaws Scale 1/1	Handle	Working length
 De Lagausie Fenestrated forceps 12mm jaws		CEV995M 220 mm
 De Lagausie Straight scissors serrated blades 8mm		CEV996M 220 mm
 De Lagausie Alligator forceps 6mm serrated jaws		CEV997M 210 mm
 De Lagausie Dissector medium curved 10mm		CEV998M 220 mm
 Page 43		
 Montupet Needle-holder straight tungsten jaws		CEV892 220 mm
 Reinberg Needle-holder straight tungsten jaws		CEV992 220 mm
 Needle-holder left curved tungsten jaws		CEV992D 220 mm

Fenestrated forceps



MCLP20

180mm

Monopolar angled handle

mors 8x3mm

Fenestrated forceps

MCLP25

L 100mm

Monopolar axial handle

mors 8x3mm



Scale 1/2
Working length «L»

MCLP20

Forceps, alligator jaws



MCLP30

L 95mm

Monopolar axial handle

mors 8x3mm



MCLP30

Curved micro scissors



MCLP40

L 95mm

Monopolar axial handle

blades 4mm



MCLP40

Monopolar hook



MCLP50

L 105mm



MCLP50

Blunt tip trocar 55mm



MCLP10

Trocar holder ring
ø 10mm



MCLP70

Scale 1/2

Double blunt retractor
180mm

MCLT10

small size

MCLT11

Large size



MCLT10

MCLT11

Blunt curved elevator 190mm

MCLT20

2mm

MCLT21

4mm



MCLT20

MCLT21

Peanut holder



200mm - \varnothing 3mm

MCLT30

Axial handle



MCLT30

Suction tube elevator

200mm - \varnothing 3mm

MCLT40

malleable



MCLT40

Nerve eye suction tube
with finger control

200mm - \varnothing 3mm

MCLT50

malleable




MCLT50



1/3



1/3



1/3


CEV504PS




the handle is resterilizable

CEV504PS


Axial handle for inserts CEV504S3 and CEV504S5






CEV100-3

1/3



7874MF-34
7874MF-40



7874MF-61

Stainless steel sterilization cases

7874MF-34	340x200x80mm
7874MF-40	400x200x90mm
7874MF-61	610x180x80mm
CEV100-3	3 perforated trays for 7874MF-61



3717001

Case 530x255x150mm with 3 levels



MCOS1oSC

Selesnick sterilization case 385x256x128mm with 2 levels



MCO731^(D)

Surgical instrument oil



CEV56

Stainless steel cap
for luer-lock cleaning port



CEV57

silicone cap
for luer-lock cleaning port
Pack 6

^(D) Distributed by Integra MicroFrance SAS (France)

3		
3717001.....		85
7		
7874MF-34.....		85
7874MF-40.....		85
7874MF-61.....		85
C		
CET105-1.....		70
CET105-2.....		70
CET105-3.....		70
CET110-1.....		70
CET110-2.....		70
CET110-3.....		70
CET112-1.....		70
CET112-3.....		70
CET112T2.....		70
CET525.....		70
CEV24F.....		62
CEV24G.....		62
CEV39L.....		72,78
CEV56.....		85
CEV57.....		85
CEV60-120.....		62
CEV60-150.....		62
CEV65.....		62
CEV81.....		60
CEV100-3.....		85
CEV101.....		61
CEV103-5.....		61
CEV103-10.....		61
CEV103AH.....		61
CEV103G3.....		75,77
CEV103G5.....		61
CEV103GH.....		61,75
CEV103P.....		61
CEV103P6.....		63,69,76
CEV104I.....		26
CEV104M.....		28
CEV105.....		67
CEV105I.....		67
CEV105P6.....		64,67,68,69
CEV107I.....		70
CEV110I.....		67
CEV110P6.....		65,67,68,69
CEV112I.....		67
CEV112P6.....		66,68,69
CEV114M.....		28
CEV122I.....		67
CEV125P6.....		64,76
CEV125Q6.....		65,66
CEV127T.....		64
CEV127T3.....		76
CEV127T10.....		65
CEV127T12.....		66
CEV130.....		55
CEV132.....		53
CEV133.....		55
CEV133-1.....		55
CEV134.....		53
CEV134BG.....		53
CEV135-5.....		57
CEV136.....		53
CEV143A3.....		76
CEV143A4.....		76
CEV143A7.....		76
CEV143B7.....		76
CEV145-3.....		68
CEV145-4.....		68
CEV145-6.....		68
CEV145-7.....		68
CEV145A3.....		64
CEV145A5.....		64
CEV145A6.....		64
CEV145A7.....		64
CEV145B6.....		64
CEV145B7.....		64
CEV150-3.....		65,68
CEV150-4.....		65,68
CEV150-6.....		68
CEV150-7.....		68
CEV150A6.....		65
CEV150A7.....		65
CEV150B6.....		65
CEV150B7.....		65
CEV153.....		69
CEV155.....		69
CEV155-1.....		69
CEV155A.....		69
CEV155B.....		69
CEV155P.....		69
CEV156.....		69
CEV156A.....		69
CEV172-1.....		68
CEV172-3.....		68
CEV172A3.....		66
CEV172A4.....		66
CEV172A6.....		66
CEV172A7.....		66
CEV172B6.....		66
CEV172B7.....		66
CEV177-5.....		68
CEV177-12.....		68
CEV177A3.....		76
CEV177A5.....		64
CEV177A10.....		65
CEV177A12.....		66
CEV178-5.....		68
CEV179A5.....		64,76
CEV179A10.....		65
CEV179A12.....		66
CEV181.....		69
CEV181-1.....		69
CEV181A.....		69
CEV203.....		75
CEV204.....		58
CEV204G45.....		58
CEV205.....		58
CEV207.....		58
CEV211.....		15
CEV220.....		13
CEV225.....		56
CEV225A.....		56
CEV229-1A.....		56
CEV229-2A.....		56
CEV229-3A.....		56
CEV229-4A.....		56
CEV229B.....		56
CEV230-1.....		56
CEV278-10.....		68
CEV278-12.....		68
CEV310.....		39
CEV321.....		15
CEV342.....		75
CEV342E.....		75
CEV350.....		75
CEV390B.....		75
CEV391B.....		73
CEV391I.....		72
CEV392D20.....		71
CEV392D30.....		71
CEV392G24.....		71
CEV392P20.....		71
CEV392R15.....		71
CEV392S20.....		71
CEV393I.....		72
CEV394B.....		73
CEV394I.....		72
CEV395B.....		73
CEV395I.....		72
CEV396I.....		72
CEV397.....		74
CEV397BG.....		74
CEV397C.....		74
CEV398I.....		72
CEV404.....		29
CEV405.....		3
CEV405BG.....		3
CEV406.....		5
CEV406I.....		4
CEV407.....		5
CEV407I.....		4

CEV410	9	CEV515M	13	CEV614-1	12
CEV413	58	CEV521	11	CEV615-1	12
CEV416	58	CEV525	3	CEV615H	50
CEV417	69	CEV525C	3	CEV616I	16
CEV418	69	CEV525M	3	CEV617I	12
CEV419	69	CEV526	31	CEV618-5N	48
CEV420	11	CEV530	9	CEV620-1	12
CEV421	11	CEV531-2	19	CEV621-1	10
CEV423	58	CEV531-3	19	CEV622I	39
CEV425	58	CEV532	19	CEV623-1	26
CEV430	21	CEV532I	18	CEV625-1	2
CEV431-65	35	CEV533-2	9	CEV625C	2
CEV432	9	CEV535	33	CEV625D	2
CEV433	9	CEV538-10	24	CEV625E	6
CEV433-1	8	CEV538T	22	CEV625H	50
CEV433D	6	CEV538T5	23	CEV625P	6
CEV433H	8	CEV540-5	23	CEV625P1	2
CEV435-1	33	CEV540-10	24	CEV626-1	30
CEV435-2	33	CEV544	23	CEV627I	30
CEV438	62	CEV544I	22	CEV628-5N	48
CEV438B	62	CEV560	38	CEV628-10	46
CEV440	63	CEV560-1	38	CEV628P5	42
CEV442	25	CEV565	38	CEV630-1A	54
CEV442E	25	CEV566	38	CEV630C	54
CEV443	57	CEV567	38	CEV632-1A	52
CEV444	70	CEV568	38	CEV633-1A	54
CEV445-5	25	CEV570	31	CEV633H	50
CEV445I	25	CEV580-2	53	CEV634-1A	52
CEV448	63	CEV580-3	55	CEV635-1	32
CEV449G	63	CEV600H	50	CEV636-1A	52
CEV453	58	CEV601-1	20	CEV637I	52
CEV453A	58	CEV601P	20	CEV638-5N	48
CEV455	63	CEV602-1	10	CEV641-1	18
CEV456	62	CEV603I	39	CEV641D	17
CEV460	5	CEV604-1	26	CEV641F	17
CEV470	3	CEV605-1	26	CEV642-1	18
CEV470-1	5	CEV605-3	26	CEV642H	50
CEV500-A10	24	CEV605BGP	27	CEV642HD	50
CEV500T	22	CEV605G	26	CEV642TF	17
CEV500T5	23	CEV605GP	27	CEV643I	17
CEV501	21	CEV605GT	26	CEV646B5	41,43,45,51
CEV501-2	21	CEV606-1	14	CEV647B5	41,43,45,51
CEV502	11	CEV607	27	CEV648B5	41,43,45,51
CEV503	29	CEV607-1	27	CEV649-5B	41,43,45,51
CEV504PS	84	CEV607C	27	CEV649-5N	48
CEV504S3	84	CEV607D	27	CEV649-10B	47
CEV504S5	84	CEV608-1	34	CEV649B5	41,43,45,51
CEV505	33	CEV608D	34	CEV649BGP	27
CEV506	15	CEV610-1	8	CEV649C5	27,41,43,45,51
CEV508	35	CEV611-1	14	CEV649G5	7
CEV511	15	CEV611H	50	CEV649M5	7
CEV511M	15	CEV612-1	14	CEV649P5	49
CEV512	15	CEV612D	16	CEV649X5	47
CEV514M	28	CEV612LH	50	CEV650T	22
CEV515	13	CEV613-1	14	CEV650T5D	23

CEV651-1	20	CEV980-3.....	59	CEV8605GT	26
CEV652I.....	20	CEV980A.....	59	CEV8607-1	27
CEV653-1.....	29	CEV980B.....	59	CEV8611H.....	50
CEV658-10.....	46	CEV984.....	59	CEV8612D.....	16
CEV660-1.....	4	CEV985.....	59	CEV8612LH.....	50
CEV661I.....	4	CEV992.....	81	CEV8615-1.....	12
CEV669B.....	52,54	CEV992B.....	73	CEV8615H.....	50
CEV669BG.....	52,54	CEV992C.....	73	CEV8616I.....	16
CEV669E.....	50	CEV992D.....	81	CEV8621-1.....	10
CEV670-1.....	2	CEV995M.....	81	CEV8625-1.....	2
CEV670D.....	6	CEV996M.....	81	CEV8625C.....	2
CEV671-1.....	4	CEV997M.....	81	CEV8625D.....	2
CEV700T5.....	23	CEV998M.....	81	CEV8625E.....	6
CEV700T5D.....	23	CEV1019-5B.....	40	CEV8625H.....	50
CEV720.....	49	CEV1019-5N.....	48	CEV8627I.....	30
CEV720B.....	49	CEV1039-5-B.....	40	CEV8633H.....	50
CEV720BR.....	44	CEV1039-5N.....	48	CEV8641-1.....	18
CEV720R.....	44	CEV1039G5.....	40	CEV8642-1.....	18
CEV721.....	11	CEV1059-10B.....	46	CEV8642H.....	50
CEV725.....	3	CEV2295A.....	56	CEV8642HD.....	50
CEV728BR.....	46	CEV6185R.....	42	CEV8646B5.....	41,43,45,51
CEV728BR5.....	44	CEV6285R.....	42	CEV8647B5.....	41,43,45,51
CEV728T5.....	42	CEV6385R.....	42	CEV8648B5.....	41,43,45,51
CEV728T10.....	46	CEV6795B.....	52,54	CEV8649-5-B.....	41,43,45,51
CEV738-5.....	23	CEV7145A3.....	64	CEV8649-5N.....	48
CEV738T5.....	23	CEV7145A5.....	64	CEV8649B5.....	41,43,45,51
CEV740-5.....	23	CEV7145A6.....	64	CEV8649BCP.....	27
CEV744T5.....	23	CEV7145A7.....	64	CEV8649C5.....	27,41,43,45,51
CEV760.....	5	CEV7145B6.....	64	CEV8649GP.....	27
CEV804.....	28	CEV7145B7.....	64	CEV8649P5.....	49
CEV805I.....	67	CEV7285R.....	44	CEV8649X5.....	47
CEV821G.....	11	CEV7625I.....	7	CEV8650T.....	22
CEV825.....	56	CEV7660I.....	7	CEV8660-1.....	4
CEV830.....	55	CEV8103D5.....	61	CEV8670D.....	6
CEV834.....	53	CEV8103G5.....	61	CEV9103D5.....	61
CEV844.....	77	CEV8104I.....	26	CEV9103G5.....	61
CEV890.....	77	CEV8107I.....	70	CEV9104I.....	26
CEV890C.....	77	CEV8143A3.....	76	CEV9107I.....	70
CEV891I.....	78	CEV8143A4.....	76	CEV9134BG.....	53
CEV891M.....	79	CEV8143A7.....	76	CEV9145A3.....	64
CEV892.....	81	CEV8143B7.....	76	CEV9145A5.....	64
CEV893I.....	78	CEV8145A3.....	64	CEV9145A6.....	64
CEV894I.....	78	CEV8145A5.....	64	CEV9145A7.....	64
CEV894M.....	79	CEV8145A6.....	64	CEV9145B6.....	64
CEV895I.....	78	CEV8145A7.....	64	CEV9145B7.....	64
CEV895M.....	79	CEV8145B6.....	64	CEV9150A3.....	65
CEV896.....	77	CEV8145B7.....	64	CEV9150A4.....	65
CEV896I.....	78	CEV8407I.....	4	CEV9150A6.....	65
CEV897.....	80	CEV8434I.....	8	CEV9150A7.....	65
CEV897C.....	80	CEV8445I.....	25	CEV9150B6.....	65
CEV898.....	77	CEV8538T.....	22	CEV9150B7.....	65
CEV898I.....	78	CEV8544I.....	22	CEV9172A3.....	66
CEV898L.....	77	CEV8604-1.....	26	CEV9172A4.....	66
CEV904.....	28	CEV8605-1.....	26	CEV9172A6.....	66
CEV965.....	23	CEV8605G.....	26	CEV9172A7.....	66

CEV9172B6.....	66	CEV86341A.....	52	MCLP70.....	82
CEV9172B7.....	66	CEV86361A.....	52	MCLT10.....	83
CEV9229-1.....	56	CEV86795B.....	52,54	MCLT11.....	83
CEV9407I.....	4	CEV96321A.....	52	MCLT20.....	83
CEV9445I.....	25	CEV96331A.....	54	MCLT21.....	83
CEV9500T.....	22	CEV96361A.....	52	MCLT30.....	83
CEV9544I.....	22	CEV96795B.....	52,54	MCLT40.....	83
CEV9605-1.....	26	CEV102910R.....	46	MCLT50.....	83
CEV9605G.....	26	CEV864910B.....	8,47	MCO731 ^(D)	85
CEV9605GT.....	26	CEVA101.....	36	MCOS10SC.....	85
CEV9607-1.....	27	CEVA102.....	36		
CEV9611-1.....	14	CEVA103.....	36	S	
CEV9611H.....	50	CEVA104.....	36	SM63 ^(D)	57
CEV9612D.....	16	CEVA110.....	37	SM65 ^(D)	57
CEV9612LH.....	50	CEVA112.....	33		
CEV9614-1.....	12	CEVA528.....	37		
CEV9615-1.....	12	CEVA529.....	37		
CEV9615H.....	50	CM104.....	28		
CEV9617I.....	12	CM106.....	28		
CEV9620-1.....	12	CM107.....	70		
CEV9625-1.....	2	CM109.....	57		
CEV9625C.....	2	CM110.....	13		
CEV9625D.....	2	CM111.....	15		
CEV9625H.....	50	CM112.....	15		
CEV9627I.....	30	CM113.....	15		
CEV9630-1A.....	54	CM115.....	13		
CEV9633H.....	50	CM120.....	13		
CEV9634-1A.....	52	CM131.....	55		
CEV9641-1.....	18	CM170-1.....	62		
CEV9642H.....	50	CM170-2.....	62		
CEV9642HD.....	50	CM170-3.....	62		
CEV9646B5.....	41,43,45,51	CM201.....	60		
CEV9647B5.....	41,43,45,51	CM201-9.....	60		
CEV9648B5.....	41,43,45,51	CM217S ^(D)	60		
CEV9649-5-B.....	41,43,45,51	CM20128H-5.....	60		
CEV9649-5N.....	48	CM20128H-10.....	60		
CEV9649B5.....	41,43,45,51	CM20132H-5.....	60		
CEV9649BGP.....	27	CM20132H10.....	60		
CEV9649C5.....	27,41,43,45,51	CM20132L-5.....	60		
CEV9649P5.....	49	CM20140H5.....	60		
CEV9649X5.....	47	CP391E.....	27,41,43,45,57		
CEV9650T.....	22	CP391M.....	27,41,43,45,57		
CEV9660-1.....	4	CP391V.....	27,41,43,45,57		
CEV9670D.....	6	CP393-2.....	51,71		
CEV9700T5D.....	23	CP393-3.....	51,71		
CEV9738T5.....	23	CP393-4.....	51,71		
CEV9744T5.....	23			M	
CEV10195C.....	40	MCEN54P.....	52,54		
CEV10195D.....	40	MCLP10.....	82		
CEV10195NA.....	48	MCLP20.....	82		
CEV10195R.....	40	MCLP25.....	82		
CEV10395D.....	40	MCLP30.....	82		
CEV10395R.....	40	MCLP40.....	82		
CEV72810R.....	46	MCLP50.....	82		
CEV82291.....	56				

^(D) Distributed by Integra MicroFrance SAS (France)

Integra® MicroFrance®

Integra Laparoscopy catalogue

Manufacturers :



Integra MicroFrance SAS (France)
Le Pavillon ■ 03160 Saint Aubin le Monial ■ France



Medtronic Xomed
6743 Southpoint Drive N.
Jacksonville 32216-0980 ■ Florida USA



Medtronic B.V.
Earl Bakkenstraat 10
6422 PJ Heerlen
The Netherlands



HEMODIA SAS
85, rue du Chêne Vert
31670 LABEGE ■ France



Swann Morton Ltd
Penn Works ■ Owlerton Green ■ Sheffield
South Yorkshire S6 2BJ ■ Royaume Uni

Manufacturers with CE class I products only :



Chemische Fabrik Dr. Weigert GmbH & Co.
Mühlenhagen
D-20539 Hamburg ■ Allemagne

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Please read carefully the instructions for use.
- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.

Additional information for EMEA Customers only :

Products mentioned in this document are CE class I or IIa devices. Contact Integra should you need any additional information on devices classification. All the medical devices mentioned on this document are CE marked according to European council directive 93/42/EEC on medical devices and its relatives.

References that are not medical devices are identified with the indication "NOT A MEDICAL DEVICE".

Products marked with the symbol [®] are distributed by Integra MicroFrance SAS (France)

Date of last révision : 06/2015

For more information, please contact customer service :

Europe, Middle-East, Africa

France / International : +33 (0) 470 67 98 00 ■ +33 (0) 470 67 98 20 (Fax) ■ microfrance@integralife.com
United Kingdom : +44 (0) 1264 312 725 ■ +44 (0) 1264 312 821 (Fax) ■ custsvcuk@integralife.com
Germany : +49 (0) 2102 5535 6200 ■ +49 (0) 2102 5536 636 (Fax) ■ custsvcgermany@integralife.com
Benelux: +32 (0)2 257 4130 ■ +32 (0)2 253 2466 (Fax) ■ custsvcbenelux@integralife.com
integralife.eu/contact

Integra and the Integra logo are registered trademarks of Integra LifeSciences Corporation or its subsidiaries in the United States and/or other countries.

MicroFrance is a trademark of Integra LifeSciences Corporation or its subsidiaries.

All other trademarks and product names are the property of their respective owners.

©2015 Integra LifeSciences Corporation. All rights reserved.

Printed in France. 0335697-1-EN

DOCUMENT FOR USE IN EUROPE, MIDDLE-EAST AND AFRICA ONLY.