

Integra[®] MicroFrance[®]

Laparoscopic Instruments for Bariatric Procedures*

Designed with world-renowned surgeons, MicroFrance Laparoscopic Instruments are known for quality and precision. These instruments are designed to perform a variety of laparoscopic procedures.




*Instruments under the medical grouping can be used in other general surgical procedures.

Bariatric Procedures

Instrument Trays

Reference	Description
7874MF-61	Stainless Steel Case; 610 x 180 x 80 mm
CEV100-3	Perforated Trays (3) for 7874MF-61




Modular Instruments

5mm Diameter

Bipolar Forceps, 410mm Working Length

Reference	Description
CEV9625H	Bipolar Insert Fenestrated Jaws for Tube, 20 x 5mm Jaws
CEV9649-5-B	Tube, Insulated Black
CEV9649B5	Tube, Non-Insulated*
CEV9649C5	Tube, Disposable Sheath, Black, 3-pack*
CEV669E	Rotatable Handle



*Jaw insert, tubes, and handle in this section are interchangeable.
Non-insulated tube CEV9649B5 must be used with disposable sheath CEV9649C5.

Modular Instruments

5mm Diameter, 410mm Working Length

Monopolar Instruments – Jaw Inserts

Reference	Description
CEV9641-1	Medium Curved Dissector Jaws
CEV9605-1	Very Curved Scissor Jaws
CEV9605G	Curved Insulated Scissor Jaws
CEV9625D	Double-Action Fenestrated Forceps
CEV9605GT	Very Curved, Titanium; Insulated Scissors
CEV607	Very Curved Insulated Scissors Insert, 6-pack**
CEV9607-1	Wire Insert**



CEV9641-1



CEV9605-1



CEV9605G



CEV9625D



CEV9605GT



CEV607
CEV9607-1

Monopolar Instruments – Tubes

Reference	Description
CEV9649-5-B	Insulated Black
CEV9646B5	Insulated Yellow
CEV9649B5	Non-Insulated**
CEV9649C5	Disposable Sheath, Black, 3-pack**



CEV9649-5B



CEV9646B5



CEV9649B5



CEV9648B5

Monopolar Instruments – Handles


Reference	Description
CEV10395D	Rotatable Angled Handle
CEV10395R	Rotatable Angled Handle
CEV10195R	Rotatable Handle with Ratchet



CEV10395R



CEV10395D



CEV10195R

**Jaw inserts, tubes, and handles in this section are interchangeable. Insulated scissors tip CEV607 must be used with wire insert CEV9607-1. Non-insulated tube CEV9649B5 must be used with disposable sheath CEV9649C5.

Modular Instruments

5mm Diameter

Forceps

Reference	Description
CEV7625I	Jaw insert; Fenestrated Forceps, 20 x 5 mm Jaws
CEV649M5	Non-Insulated Tube; Graduated Metallic Tube w/ Ruler
CEV628-10	Handle; Rotatable Axial Handle w/ Ratchet



CEV7625I



CEV649M5



CEV628-10

Forceps, 410mm Working Length

Reference	Description
CEV9615-1	Jaw Insert; Alligator Forceps, 15 mm Jaws
CEV9649X5	Non-Insulated Tube; Metallic
CEV102910R	Handle; Rotatable Angled Handle with Ratchet



CEV9615-1



CEV9649X5



CEV102910R

Monopolar Hook, 450mm Working Length

Reference	Description
CEV2295A	Monopolar Hook w/ Size 5 Jaws, 3-Pack
CEV9229-1	Handle; Round Hook Handle



CEV2295A



CEV9229-1

Cables

Bipolar

Reference	Description	Qty
CP393-3	For Erbe, Storz, and Wolf Generators, 4.5 M	1
CP393-4	For Martin, Berchtold, and Aesculap Generators, 4.5 M	1
CP393-2	For Valleylab and Lamidey Generators, 4.5 M	1

Monopolar

Reference	Description	Qty
CP391E	For Erbe Generators, 3.5 M, with 5mm connection tip	1
CP391M	For Martin – Lamidey Generators, 3.5 M, with 4mm connection tip	1
CP391V	For Valleylab Generators, 3.5 M, with 8mm connection tip	1

Standard Instruments

5mm Diameter


Retractors, 160mm Working Length, 5mm Width

Reference	Description
CEV419	Mosnier Double Retractor, Blunt (pair)



Needle Holders, 5mm Width

Reference	Description
CEV9738T5	Needle Holder with Curved Tungsten Carbide Jaws and Large Handle, Ratchet Down, 455mm Working Length
CEV9700T5D	Needle Holder with Straight Tungsten Carbide Jaws and Large Handle, Ratchet Down, 450mm Working Length



Standard Instruments

10mm Diameter


Retractor, 330mm Working Length, 10mm Width

Reference	Description
CEV453A	Articulated Retractor with 5 Blades



Wire Passer, 460mm Working Length, 10mm Width

Reference	Description
CEV438	Bertrand Articulated Wire Pass



Trocars

5mm Diameter

Valve Heads

Reference	Description
CEV177A5	Valve Head, Magnetic
CEV127T	Valve Head, Silicone



CEV177A5



CEV127T

Tips

Reference	Description
CEV9145A3	Sharp Tip, 175mm Long
CEV9145A5	Smooth Conical Tip, 175mm Long



CEV9145A3



CEV9145A5

Cannulae

Reference	Description
CEV9145B6	Threaded Cannula, 120mm Long
CEV9145A6	Threaded Cannula with Stopcock, 120mm Long
CEV9145B7	Smooth Cannula, 120mm Long
CEV9145A7	Smooth Cannula with Stopcock, 120mm Long



CEV9145B6



CEV9145A6



CEV9145B7



CEV9145A7

10mm Diameter

Valve Heads

Reference	Description
CEV177A10	Valve Head, Magnetic
CEV127T10	Valve Head, Silicone



CEV177A10



CEV127T10

Tips

Reference	Description
CEV9150A3	Sharp Tip, 175mm Long
CEV9150A4	Blunt Conical Tip, 175mm Long



CEV9150A3



CEV9150A4

Cannulae

Reference	Description
CEV9150B6	Threaded Cannula, 130mm Long
CEV9150A6	Threaded Cannula with Stopcock, 130mm Long
CEV9150B7	Smooth Cannula, 130mm Long
CEV9150A7	Smooth Cannula with Stopcock, 130mm Long



CEV9150B6



CEV9150A6



CEV9150B7



CEV9150A7

5.5mm Diameter Accessories

Reference	Description	Qty
CEV179A5	Replacement Valves for Magnetic Head CEV177A5	6
CEV125P6	White Seal for Trocar Valve, 3.5mm & 5mm Diameter	6
CEV105P6	Sealing Caps for 5mm Diameter Trocar Heads, Blue	6

Miscellaneous Accessories

CEV60-150	Insufflation Needle with Luer-Lock, 150mm L	1
CEV155	Reduction Sleeve from 10mm to 5mm Diameter, 135mm Length	1


10.5mm Diameter Accessories

Reference	Description	Qty
CEV179A10	Replacement Valves for Magnetic Head CEV177A10	6
CEV125Q6	White Seal for Trocar Valve, 10mm and 12mm Diameter	6
CEV110P6	Sealing Caps for 10mm Diameter Trocar Heads, Red	6

Accessories


Suction Irrigation, 400mm Working Length, 5mm Diameter

Reference	Description
CEV9103D5	Tube
CEV103AH	Handle



Flexible Tension Arm

Reference	Description
CEV984	Operating Table Clamp
CEV980A	Mosnier Flexible Tension Arm, 500 mm
CEV980B	Mosnier Flexible Tension Arm, 600 mm
CEV980-3	Mosnier Instrument Holder



Mosnier Flexible Tension Arm with Instrument Holder, parts sold separately
CEV980A or CEV980B, shown here with CEV980-3

The MicroFrance® brand is a premier line of ENT instruments, designed with world-renowned ENT surgeons to facilitate new techniques and improve patient outcomes.

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.
- Consult product labels and inserts for any indication, contraindications, hazards, warnings, precautions, and instructions for use.

For more information or to place an order, please contact:

United States, Canada, Asia, Pacific, Latin America

USA 800-654-2873 ▪ 888-980-7742 fax

International +1 609-936-5400 ▪ +1 609-750-4259 fax

integralife.com/contact

