

Integra[®] MicroFrance[®] Modular Bipolar Forceps

MicroFrance[®] Modular Bipolar Forceps help provide versatility for laparoscopic surgery, with multiple jaw tips and tube styles available.

Designed with world-renowned surgeons, MicroFrance Laparoscopic Instruments are known for quality and precision. These instruments are designed for a variety of laparoscopic procedures.



Integra® MicroFrance® Modular Bipolar Forceps

Laparoscopic Bipolar Jaw Inserts

Graspers

For 5mm diameter tubes

310 mm length
CEV611H

210 mm length
CEV8611H

410 mm length
CEV9611H



Johann Fenestrated Jaws

For 5mm diameter tubes

310 mm length
CEV625H

210 mm length
CEV8625H

410 mm length
CEV9625H



Alligator Jaws

15mm, for 5mm diameter tubes

310 mm length
CEV615H

210 mm length
CEV8615H

410 mm length
CEV9615H



Dissectors

18mm jaws, for 5mm diameter tubes

Left-curved, 310 mm length
CEV642H

Left-curved, 210 mm length
CEV8642H

Left-curved, 410 mm length
CEV9642H



Triangle Fenestrated Jaws

For 5mm diameter tubes

310 mm length
CEV633H

210 mm length
CEV8633H

410 mm length
CEV9633H



Right-curved, 310 mm length
CEV642HD

Right-curved, 210 mm length
CEV8642HD

Right-curved, 410 mm length
CEV9642HD



Rotatable Insulated Tubes

Insulated Tubes

5mm diameter ■ 310mm length

Yellow

CEV646B5



Blue

CEV647B5



Green

CEV648B5



Black

CEV649-5B



Bipolar Cables

4.5 m length for Valleylab
and Lamidey generators

CP393-2



4.5 m length for Erbe,
Storz, and Wolf generators

CP393-3



4.5 m length for
Martin, Berchtold, and
Aesculap generators

CP393-4



Bipolar Handle

Large Handle

Bipolar ■ For 5 mm diameter, rotatable
insulated tubes and jaw inserts





CEV669E



Bipolar Jaw Inserts

Description	Working Length			Tip
	210mm	310mm	410mm	
Graspers	CEV8611H	CEV611H	CEV9611H	
Alligator Jaws	CEV8615H	CEV615H	CEV9615H	
Triangle Fenestrated Jaws	CEV8633H	CEV633H	CEV9633H	
Johann Fenestrated Jaws	CEV8625H	CEV625H	CEV9625H	
Dissectors, Left-Curved	CEV8642H	CEV642H	CEV9642H	
Dissectors, Right-Curved	CEV8642HD	CEV642HD	CEV9642HD	

Insulated Tubes

Description	Working Length	Tube
		310mm
Insulated Tubes, Yellow	CEV646B5	
Insulated Tubes, Blue	CEV647B5	
Insulated Tubes, Green	CEV648B5	
Insulated Tubes, Black	CEV649-5B	

The MicroFrance® brand is a premier line of ENT instruments, designed with world-renowned ENT surgeons to facilitate new techniques and improve patient outcomes.

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.
- Consult product labels and inserts for any indication, contraindications, hazards, warnings, precautions, and instructions for use.

For more information or to place an order, please contact:

United States, Canada, Asia, Pacific, Latin America

USA 800-654-2873 ■ 888-980-7742 fax

International +1 609-936-5400 ■ +1 609-750-4259 fax

integralife.com/contact

