


Integra®

Neurosurgery IVD Rongeurs

Limit uncertainty with
durability and quality



INTEGRA®
LIMIT UNCERTAINTY

Integra Ruggles Redmond offers a wide selection of inter-vertebral disc (IVD) rongeurs. With the multiple pattern, tip and coating variations, Integra® Ruggles™ Redmond™ IVD Rongeurs are sure to meet the needs of the surgeon.

Wave IVD Rongeurs

Wave design allows the user to identify impurities, remove and clean the shaft of the instrument

Available in cup sizes 2 x 10mm
to 4 x 10mm

Shaft lengths in 7" (180mm) and 9" (230mm)

Straight, angled-up and angled-down jaws



With teeth

Without teeth

Bayonet IVD Rongeurs

Bayonet design allows the user to work down a channel during microsurgery while keeping their hand out of the field of view

Available in Ultra-Coat™ and Diamond™ coating

4 ¾" (120mm) or 8" (200mm) working length

Available in cup sizes 1.5 x 10mm to 4 x 10mm

Straight, angled-up and angled-down jaws



Standard IVD Pituitary Rongeurs

Available in cup sizes 1.5 x 10mm to 5 x 10mm

Shaft lengths from 6" (153mm) to 12" (305mm)

Straight, angled-up and angled-down jaws allow access to multiple anatomic structures

Delicate patterns available in Cushing, Love-Gruenwald and Spurling styles



Decker Forceps

Curved left and right for unique reaches

6" (153mm) shaft length



Williams Style Micro Disc IVD with 1 x 2 Teeth

1 x 2 teeth on end of cup for more secure grasping power

Shaft lengths in 7" (180mm) to 10" (254mm)



Wave IVD Rongeurs - Stainless Steel

Catalog Number	Description
RR5712	Cushing, 7" (180mm), straight, 2x10mm bite
RR5722	Cushing, 7" (180mm), angled up, 2x10mm bite
RR5732	Cushing, 7" (180mm), angled down, 2x10mm bite
RR5912	Cushing, 9" (230mm), straight, 2x10mm bite
RR5922	Cushing, 9" (230mm), angled up, 2x10mm bite
RR5932	Cushing, 9" (230mm), angled down, 2x10mm bite
RR5743	Love-Gruenwald, 7" (180mm), straight, 3x10mm bite
RR5753	Love-Gruenwald, 7" (180mm), angled, up 3x10mm bite
RR5763	Love-Gruenwald, 7" (180mm), angled down, 3x10mm bite
RR5943	Love-Gruenwald, 9" (230mm), straight, 3x10mm bite
RR5953	Love-Gruenwald, 9" (230mm), angled up, 3x10mm bite
RR5963	Love-Gruenwald, 9" (230mm), angled down, 3x10mm bite
RR5744	Spurling, 7" (180mm), straight, 4x10mm bite
RR5754	Spurling, 7" (180mm), angled up, 4x10mm bite
RR5764	Spurling, 7" (180mm), angled down, 4x10mm bite
RR5944	Spurling, 9" (230mm), straight, 4x10mm bite
RR5954	Spurling, 9" (230mm), angled up, 4x10mm bite
RR5964	Spurling, 9" (230mm), angled down, 4x10mm bite



Angled-Up



Straight



Angled-Down

IVD Sizing Chart

IVD Rongeurs Cup Size					
1.5x10mm	2x10mm	2x12mm	2.5x7mm	3x7mm	3x10mm
3x12mm	4x7mm	4x8mm	4x10mm	4x14mm	5x10mm
5x14mm	6x10mm	6x12mm			

Delicate Pattern IVD Rongeurs

Catalog Number	Description
RN7127	Decker, 6" (153mm), right, 2x6mm bite
RN7128	Decker, 6" (153mm), left, 2x6mm bite
RN7129	Decker, 6" (153mm), down, 2x6mm bite
RN7130	Decker, 6" (153mm), up, 2x6mm bite
RN7131	Decker, 6" (153mm), straight, 2x6mm bite
RN7132	Micro Disc Rongeur, 5-1/2" (140mm), straight, 1x2 teeth
RB6154	Williams Style Micro Disc, 7" (180mm), straight, 1x2 teeth
RB6155	Williams Style Micro Disc, 10" (254mm), straight, 1x2 teeth

Bayonet IVD Rongeurs

Catalog Number	Description
RD6811	Diamond, Cushing, 200mm, straight, 1.5x10mm bite
RD6812	Diamond, Cushing, 200mm, angled up, 1.5x10mm bite
RD6813	Diamond, Cushing, 200mm, angled down, 1.5x10mm bite
RD6814	Diamond, Cushing, 200mm, straight, 2x10mm bite
RD6815	Diamond, Cushing, 200mm, angled up, 2x10mm bite
RD6816	Diamond, Cushing, 200mm, angled down, 2x10mm bite
RD6831	Diamond, Love-Gruenwald, 200mm, straight, 3x10mm bite
RD6832	Diamond, Love-Gruenwald, 200mm, angled up, 3x10mm bite
RD6833	Diamond, Love-Gruenwald, 200mm, angled down, 3x10mm bite
RD6851	Diamond, Spurling, 200mm, straight, 4x10mm bite
RD6852	Diamond, Spurling, 200mm, angled up, 4x10mm bite
RD6853	Diamond, Spurling, 200mm, angled down, 4x10mm bite
RB6801	Ultra-Coat, Cushing, 120mm, straight, 1.5x10mm bite
RB6802	Ultra-Coat, Cushing, 120mm, angled up, 1.5x10mm bite
RB6803	Ultra-Coat, Cushing, 120mm, angled down, 1.5x10mm bite
RB6804	Ultra-Coat, Cushing, 120mm, straight, 2x10mm bite
RB6805	Ultra-Coat, Cushing, 120mm, angled up, 2x10mm bite
RB6806	Ultra-Coat, Cushing, 120mm, angled down, 2x10mm bite
RB6821	Ultra-Coat, Love-Gruenwald, 120mm, straight, 3x10mm bite
RB6822	Ultra-Coat, Love-Gruenwald, 120mm, angled up, 3x10mm bite
RB6823	Ultra-Coat, Love-Gruenwald, 120mm, angled down, 3x10mm bite
RB6841	Ultra-Coat, Spurling, 120mm, straight, 4x10mm bite
RB6842	Ultra-Coat, Spurling, 120mm, angled up, 4x10mm bite
RB6843	Ultra-Coat, Spurling, 120mm, angled down, 4x10mm bite

For a complete listing of IVD rongeurs refer to the Integra® Ruggles™ Redmond™ Catalog

For Professional Use only.

Integra®

Neurosurgical IVD Rongeurs

For more information or to place an order, please contact:

Integra ■ 589 Davies Drive, York PA 17042

877-468-5572 USA ■ 800-275-7192 fax

914-789-7094 outside USA ■ 609-750-4218 outside USA fax

integralife.com

Integra and the Integra logo are registered trademarks of Integra LifeSciences Corporation and its subsidiaries in the United States and/or other countries. Ruggles, Redmond, Ultra-Coat, and Diamond are trademarks of Integra LifeSciences Corporation and its subsidiaries.

©2012 Integra LifeSciences Corporation. All rights reserved. Printed in the USA. 5K PM12019-02/12