


# Integra®

Ruggles®-Redmond™ Flexi-Spine™  
Retraction System


Limit uncertainty with  
a robust retraction system  
for spinal procedures



INTEGRA®  
LIMIT UNCERTAINTY

### Versatility... Strength... Ease of use...

With multiple blade and retractor options, the Ruggles®-Redmond™ Flexi-Spine™ Retraction System is used for single-level and multi-level procedures in cervical fusion and microlumbar discectomy cases. The system remains secure during the entire procedure and is designed for easy identification, insertion, and retraction. The Flexi-Spine™ Retraction System meets the needs of today's spine surgeon.




Hooks  
Reference RU2407




Blunt Blades  
Reference RU2405




Short Teeth Blades  
Reference RU2404



Long Teeth Blades  
Reference RU2406



Single-Level Configuration




Multi-Level Configuration

# Flexi-Spine™ Retraction System

## Transverse & Longitudinal Retractors

- Dual blade attachment points allow for expanded retraction to perform single-level and multi-level surgeries.
- Transverse and Longitudinal Retractors for 4-sided retraction.
- Specially-designed curve in retractor arms for easier blade attachment.
- Blades attach to retractor arm along entire width for a more secure fit.



Transverse Retractor  
Reference RU2408



Longitudinal Retractor  
Reference RU2409

## Retractor Handle

Spring-loaded handle for easy insertion and release.



Reference RU2403

## Short-teeth Blades

- Color-coded blades and hooks for easy identification.
- 25mm to 75mm lengths for multiple anatomies and procedures.
- Aluminum construction for high strength with radiolucency



Reference RU2404

## Depth Gauge

Color-coded depth gauge to eliminate guesswork.



Reference RU2422

## Hook and Blade Assortment

- 15mm and 20mm widths
- Hooks, blunt blades, short teeth blades and long teeth blades



### Flexi-Spine™ Retraction System

Reference	Description
RU2400	Complete set including 2 retractor handles and sterilization case
RU2408	Transverse retractor body, double level
RU2409	Longitudinal retractor body, double level
RU2422	Flexi-Spine™ Depth Gauge (color-coded)
RU2403	Retractor handle (2 per set)
RU2404	Short teeth blade set, 12 pair
RU2405	Blunt blade set, 12 pair
RU2406	Long teeth blade set, 12 pair
RU2407	Hook set, 6 hooks
RS2400	Flexi-Spine™ Instrument Tray

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.
- Consult product labels and inserts for any indication, contraindications, hazards, warnings, precautions, and instructions for use.

**For more information or to place an order, please contact:**

**United States, Canada, Asia, Pacific, Latin America**

USA 800-654-2873 ▪ 888-980-7742 fax

International +1 609-936-5400 ▪ +1 609-750-4259 fax

[integralife.com/contact](http://integralife.com/contact)