


Integra®

Ruggles®-Redmond™ Brain Spatulas

Limit uncertainty through
atraumatic retraction




INTEGRA®
LIMIT UNCERTAINTY

- Available in copper-coated black or clear silicone
- Malleable
- Rounded edges help to reduce the risk of lacerations
- Additional silicone at the tips and edges aid in atraumatic use of the instrument
- Available with or without posts

Hand-Held Silicone Brain Spatulas

Black


Reference	Description
RNo711	Small, 6.5 x 10mm
RNo712	Medium, 9 x 15mm
RNo713	Large, 11.5 x 20mm
RNo714	X-Large, 14 x 25mm
RNo710	Hand-Held Set (4) Includes one (1) each of items above



Silicone Brain Spatulas

Black

Reference	Description
RNo711N	Narrow, 6.5mm tip, 6" (152mm) length
RNo711W	Wide, 10mm tip, 6" (152mm) length
RNo712N	Narrow, 9mm tip, 6" (152mm) length
RNo712W	Wide, 15mm tip, 6" (152mm) length
RNo713N	Narrow, 11.5mm tip, 6" (152mm) length
RNo713W	Wide, 20mm tip, 6" (152mm) length
RNo714N	Narrow, 14mm tip, 6" (152mm) length
RNo714W	Wide, 25.5mm tip, 6" (152mm) length



Hand-Held Silicone Brain Spatulas

Clear Reinforced

Reference	Description
NL870-0204	1/4" x 8"
NL870-0206	3/8" x 8"
NL870-0208	1/2" x 8"
NL870-0210	5/8" x 8"
NL870-0214	7/8" x 8"
NL870-0200	Hand-Held Set (5), Includes one (1) each of items above



NL870-0204



NL870-0208



NL870-0214

Silicone Brain Spatulas

Clear Reinforced

Reference	Description
NL7702701	Post: 5cm x 3mm, Body: 7.6 x 0.9cm
NL7702702	Post: 5cm x 3mm, Body: 7.6 x 1.3cm
NL7702703	Post: 5cm x 3mm, Body: 7.6 x 1.1cm
NL7702704	Post: 5cm x 3mm, Body: 7.6 x 2.2cm
NL7702705	Post: 5cm x 3mm, Body: 7.0 x 1.5cm
NL7702707	Post: 5cm x 3mm, Body: 9.5 x 2.0cm
NL7702708	Post: 5cm x 3mm, Body: 7.5 x 2.5cm



NL7702701



NL7702705



NL7702708

T. Fukushima Brain Spatulas

Copper, Ebenol-C Coated

Reference	Description
R8959	Small 2.2mm tip, 87/8" (225mm) length
R8960	Medium 3.8mm tip, 87/8" (225mm) length
R8961	Large 5.8mm tip, 87/8" (225mm) length
R8998	X-Large 7.0mm tip, 87/8" (225mm) length

- Malleable
- Tapered tips
- Non-reflecting black



T. Fukushima Brain Spatulas

Rigid, Stainless Steel

Reference	Description
R8962	Curved 60° tip, 8" (205mm) length
R8963	Small 2.2mm tip, 8" (205mm) length
R8964	Medium 3.8mm tip, 8" (205mm) length
R8999	Large 5.8mm tip, 8" (205mm) length
R8979	X-Large 7.0mm tip, 8" (205mm) length

- Rigid
- Tapered tips
- Curved or straight styles available



Davis Brain Spatulas

Reference	Description
R0924	7" (178mm) long, 1/4" (6mm) wide
R0925	7" (178mm) long, 3/8" (10mm) wide
R0926	7" (178mm) long, 1/2" (13mm) wide
R0927	7" (178mm) long, 5/8" (16mm) wide
R0928	7" (178mm) long, 3/4" (19mm) wide
R0929	7" (178mm) long, 7/8" (22mm) wide
R0930	7" (178mm) long, 1" (25mm) wide



R0925



R0928



R0930

Davis Brain Spatulas

Insulated, Blue Silicone Coating


Reference	Description
R0924A	7" (178mm) long, 1/4" (6mm) wide
R0925A	7" (178mm) long, 3/8" (10mm) wide
R0926A	7" (178mm) long, 1/2" (13mm) wide
R0927A	7" (178mm) long, 5/8" (16mm) wide
R0928A	7" (178mm) long, 3/4" (19mm) wide
R0929A	7" (178mm) long, 7/8" (22mm) wide
R0930A	7" (178mm) long, 1" (25mm) wide



R0924A



R0927A



R0930A

Integra®

Ruggles®-Redmond™ Brain Spatulas

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.
- Consult product labels and inserts for any indication, contraindications, hazards, warnings, precautions, and instructions for use.

For more information or to place an order, please contact:

United States, Canada, Asia, Pacific, Latin America

USA 800-654-2873 ▪ 888-980-7742 fax

International +1 609-936-5400 ▪ +1 609-750-4259 fax

integralife.com/contact

Integra York PA, Inc.
589 Davies Drive
York, PA 17402 ▪ USA