



# DP11A

FLEXIBLE VIDEO RHINOLARYNGOSCOPE



# New Disposable Solution For Nasal Endoscopy



## Optical system

■ Field of view	100°
■ Depth of field	3-150mm
■ Direction of view	0°

## Insertion portion

■ Insertion tube OD	3.2mm
■ Angulation range	U 130° / D 130°
■ Working length	330mm±10%
■ Extension cable	USB、REMO
■ Max. resolution	1280×720



Avoid the risk of cross infection.



Multiple connection modes are not subject to the environment and a single device.



The ultra-fine smooth insertion tube makes the patient more "insensitive".

ZHONGSHAN WESEE MEDITECH CO.,LTD.

Address: 6F, Block C, 2Jianshuo Str. HoChung Village, Nanlang, Zhongshan 528451 P.R. China

Email:info@weseemed.com



WESEE MEDITECH