

MLX 300

We Are Passionate About Putting Surgery in the Best Light

Sunlight is the ideal light for humans to work. It gives you the color distinction you need in surgery. The MLX Xenon light source gives you true white light, not harsh actinic white light with blue shading.

Increased Lamp Warranty - 1000 hours

- True warranty, not prorated
- Better cooling leads to longer lamp life

Membrane Touchpad Controls

- Precise light attenuation in 5% increments
- Lamp and System Age meters

New Floorstand

- Dedicated to MLX, part number: 001337
- New large casters - 75 mm vs. 35 mm
- 10° angle for better viewing of LCD display
- Lower center of gravity

Stand-by Switch

- System starts up in Stand-by Mode
- Allows electronic self-check
- Put unit into Stand-by to block light when removing the cable
- Starts up at the same light output setting



MLX Light Source And Accessories

Catalog Number	Description
00MLX	MLX 300W Xenon light source with 20' power cord
00MLXEU	MLX 300W Xenon light source with 6.1m EU power cord
00MLXAU	MLX 300W Xenon light source with 6.1m Australian power cord
00MLXUK	MLX 300W Xenon light source with 6.1m UK power cord
001337	Floorstand for MLX 300W Xenon light source

