



INSTRUMENTOS Y DISPOSITIVOS DE EXCELENCIA PARA OTORRINOLARINGOLOGÍA

CATÁLOGO DE PRODUCTOS



Tita - Prosthesis type Partial Vario Campana
(Replacement of Malleus & Incus)

Material : CP Titanium Medical Grade ASTM F 67

REF Number	Shaft Diameter	Disk Diameter	Functional Length	Total Length	Units/Box Gamma Sterile
E1101	0.30	3.00	5.00-1.75	5.00	1 Unit



Tita - Prosthesis type Total Vario Cask
(Replacement of Complete Ossicular chain)

Material : CP Titanium Medical Grade ASTM F 67


REF Number	Shaft Diameter	Disk Diameter	Functional Length	Total Length	Units/Box Gamma Sterile
E1102	0.30	3.00	7.00-3.00	7.00	1 Unit



Titanium Universal Middle Ear Prosthesis

Material : CP Titanium Medical Grade ASTM F 67 & Fluoroplastic

REF Number	Shaft Diameter	Disk Diameter	Connector Diameter	Total Length	Units/Box Gamma Sterile
E1138	0.30	3.00	1.00	7.40	1 Unit



TitaHap prosthesis type PORP & TORP

Material : CP Titanium Medical Grade ASTM F 67 & HAp

REF Number	Description	Units/Box Gamma Sterile
E1207	TitaHap prosthesis type Vario Partial (PORP)	1 Unit
E1208	TitaHap prosthesis type Vario Total (TORP)	1 Unit



Clip Stapes Piston

Material : CP Titanium Medical Grade ASTM F 67

REF Number	Functional Length	Shiaft Diameter	Units/Box Gamma Sterile
E1121	7.00 to 3.50	0.40	1 Unit



Tita - prosthesis Piston type Loop

Material : CP Titanium Medical Grade ASTM F 67 & Fluoroplastic

REF Number	Functional Length	Shiaft Diameter	Units/Box Gamma Sterile
E1103	7.00 to 3.50	0.40	1 Unit



Tita - PTFE Vario Piston type Loop

Material : CP Titanium Medical Grade ASTM F 67 & Fluoroplastic (PTFE)

REF Number	Functional Length	Shiaft Diameter	Units/Box Gamma Sterile
E4119	7.00 to 4.00	0.40	1 Unit



Super Titanium Vario Piston type Loop			
Material : CP Titanium Medical Grade ASTM F 67 & Fluoroplastic			
REF Number	Functional Length	Shiaft Diameter	Units/Box Gamma Sterile
E9131	7.00 to 3.50	0.40	1 Unit



Super Tita - PTFE Vario Piston type Loop			
Material : CP Titanium Medical Grade ASTM F 67 & Fluoroplastic (PTFE)			
REF Number	Functional Length	Shiaft Diameter	Units/Box Gamma Sterile
E9132	7.00 to 4.00	0.40	1 Unit



Super Elastic NiTiNOL Vario Piston type Loop			
Material : NiTiNOL & Fluoroplastic			
REF Number	Functional Length	Shiaft Diameter	Units/Box Gamma Sterile
E5120	7.00 to 3.50	0.40	1 Unit



Super Elastic NiTiNOL - PTFE Vario Piston type Loop			
Material : NiTiNOL & Fluoroplastic (PTFE)			
REF Number	Functional Length	Shiaft Diameter	Units/Box Gamma Sterile
E8130	7.00 to 4.00	0.40	1 Unit



Titanium Stapes Piston with Integrated PTFE Shoe			
Material : CP Titanium Medical Grade ASTM F 67 & Fluoroplastic (PTFE)			
REF Number	Functional Length	Total Length	Units/Box Gamma Sterile
E1122	7.00 to 3.50	7.00	1 Unit



Fluoroplastic Partial Prosthesis (PORP)					
Material : Fluoroplastic					
REF Number	Shaft Diameter	Disk Diameter	Functional Length	Total Length	Units/Box Gamma Sterile
E2133	1.8/1.20	3.00	5.00-2.00	5.00	1 Unit



Fluoroplastic Total Prosthesis (TORP)					
Material : Fluoroplastic					
REF Number	Shaft Diameter	Disk Diameter	Functional Length	Total Length	Units/Box Gamma Sterile
E2134	0.8	3.00	7.00-3.00	7.00	1 Unit



Fluoroplastic Offset Partial Prosthesis (PORP)

Material : Fluoroplastic

REF Number	Shaft Diameter OD/ID	Functional Length	Total Length	Units/Box ETO Sterile
E2141	1.8/1.20	5.00-2.00	5.00	1 Unit



Fluoroplastic Offset Total Prosthesis (TORP)

Material : Fluoroplastic

REF Number	Shaft Diameter	Functional Length	Total Length	Units/Box ETO Sterile
E2142	0.80	7.00-3.00	7.00	1 Unit



Fluoroplastic Universal Middle Ear Prosthesis Offset

Material : Fluoroplastic

In case of TORP-Shaft OD 0.8mm | In case of PORP-Shaft ID 1.20mm

REF Number	Functional Length	Units/Box ETO Sterile
E2139	7.00 to 2.00	1 Unit



Fluoroplastic Universal Middle Ear Prosthesis Centered

Material : Fluoroplastic

In case of TORP-Shaft OD 0.8mm | In case of PORP-Shaft ID 1.20mm

REF Number	Functional Length	Units/Box ETO Sterile
E2140	7 .00 to 2.00	1 Unit



Fluoroplastic Stapes Piston

Material : Fluoroplastic

REF Number	Shaft Diameter	Loop Inner Diameter	Functional Length	Total Length	Units/Box ETO Sterile
E2127	0.5	0.60	6.00-3.50	6.00	1 Unit
E2128	0.4	0.60	6.00-3.50	6.00	1 Unit
E2129	0.6	0.60	6.00-3.50	6.00	1 Unit



Middle Ear Implant Cutter

Material : Stainless Steel


REF Number	Description	Units/Box
EI01	Middle Ear Implant Cutter	1 Unit, Non-Sterile



Middle Ear Implant Sizer

Material : Stainless Steel


REF Number	Description	Units/Box Non-Sterile
EI02	Middle Ear Implant Sizer	1 Unit



Tita - prosthesis vent tube type collar button

Material : CP Titanium Medical Grade ASTM F 67

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box Gamma Sterile
E1104	1.40	2.55	2.55	1.20	10 Units
E1105	1.30	2.55	2.55	1.20	10 Units
E1106	1.50	2.75	2.75	1.20	10 Units



Shepard Titanium

Material : CP Titanium Medical Grade ASTM F 67

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box Gamma Sterile
EI 151	1.14	2.40	2.40	2.20	10 Units



Armstrong I Armstrong Blue I Armstrong with Tail

Material : Fluoroplastic


REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box Gamma Sterile
E2109	1.14	3.54	2.40	3.85	10 Units
E2109B	1.14	3.54	2.40	3.85	10 Units
E2109T	1.14	3.54	2.40	3.85	10 Units



Modified Armstrong I Modified Armstrong Blue Modified Armstrong with Tail

Material : Fluoroplastic

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box Gamma Sterile
E2110	1.14	3.46	2.50	3.04	10 Units
E2110B	1.14	3.46	2.50	3.04	10 Units
E2110T	1.14	3.46	2.50	3.04	10 Units



Donaldson

Material : Fluoroplastic


REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box ETO Sterile
E2124	1.14	2.3	2.3	2.2	10 Units
E2124B	1.14	2.3	2.3	2.2	10 Units



Reuter Bobbin

Material : Fluoroplastic


REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box ETO Sterile
E2137	1.14	2.3	2.3	1.6	10 Units
E2137B	1.14	2.3	2.3	1.6	10 Units



Shah and Mini Shah

Material : Fluoroplastic


REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box ETO Sterile
E2137	1.14	3.4	2.3	2.4	10 Units
E2112B	1.14	3.4	2.3	2.4	10 Units
E2112T	1.14	3.4	2.3	2.4	10 Units
E2112W	1.14	3.4	2.3	2.4	10 Units
E2112BW	1.14	3.4	2.3	2.4	10 Units
E2112B	1.14	3.4	2.3	2.4	10 Units
E2113	0.76	2.2	1.5	1.5	10 Units
E2113B	0.76	2.2	1.5	1.5	10 Units
E2113T	0.76	2.2	1.5	1.5	10 Units
E2113W	0.76	2.2	1.5	1.5	10 Units



Shepard & Mini Shepard

Material : Fluoroplastic

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box ETO Sterile
E2114	1.14	2.4	2.4	2.2	10 Units
E2114B	1.14	2.4	2.4	2.2	10 Units
E2114T	1.14	2.4	2.4	2.2	10 Units
E2114W	1.14	2.4	2.4	2.2	10 Units
E2114BW	1.14	2.4	2.4	2.2	10 Units
E2114BT	1.14	2.4	2.4	2.2	10 Units
E2138	1.00	2.4	2.4	2.2	10 Units
E2138W	1.00	2.4	2.4	2.2	10 Units
E2138T	1.00	2.4	2.4	2.2	10 Units
E2138B	1.00	2.4	2.4	2.2	10 Units
E2138BW	1.00	2.4	2.4	2.2	10 Units
E2138BT	1.00	2.4	2.4	2.2	10 Units



Collar Button (ID 1.27 mm) & (ID 1.14 mm)

Material : Fluoroplastic


REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box ETO Sterile
E2123	1.27	3.0	3.0	2.15	10 Units
E2123B	1.27	3.0	3.0	2.15	10 Units
E2123T	1.27	3.0	3.0	2.15	10 Units
E2125	1.14	2.50	2.50	2.0	10 Units
E2125B	1.14	2.50	2.50	2.0	10 Units
E2125T	1.14	2.50	2.50	2.0	10 Units



Baxter Beveled Bobbin

Material : Fluoroplastic

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box ETO Sterile
E2136	1.27	2.8	2.8	2.2	10 Units
E2136B	1.27	2.8	2.8	2.2	10 Units



Bevel Bobbin

Material : Fluoroplastic


REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box ETO Sterile
E2117	1.14	2.8	2.8	1.7	10 Units
E2117B	1.14	2.8	2.8	1.7	10 Units



Straight Tube

Material : Fluoroplastic

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box ETO Sterile
E2116	1.14	2.7	1.54	7.25	10 Units



Donaldson Silicone

Material : Medical Grade Silicone


REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box ETO Sterile
E3111	1.10	2.3	2.3	2.2	10 Units



T - Tube

Material : Medical Grade Silicone

REF Number	Inner Diameter	Inner Flange Diameter	Total Length	Units/Box ETO Sterile
E3115S	1.14	8.0	6.0	10 Units
E3115M	1.14	8.0	9.0	10 Units
E3115L	1.14	8.0	12.0	10 Units



Paparella

Material : Medical Grade Silicone

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box ETO Sterile
E3126	1.10	3.5	2.5	1.65	10 Units



Gold vent Tube Type Collar Button With Wire

Material : Gold Plated Silver

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box ETO Sterile
E1170	1.25	2.55	2.55	1.20	10 Units



Nasal Septal Button

Material : Medical grade Silicone

REF Number	Description	Size (cm)	Units/Box ETO Sterile
E3145S	Nasal Septal Button - Standard	3.0	1 Unit
E3145L	Nasal Septal Button - Large	5.0	1 Unit
E3145E	Nasal Septal Button - Eliptical	3.0	1 Unit
E3145EL	Nasal Septal Button - Eliptical Large	5.0	1 Unit



Internal Nasal Airway Splint

Material : Medical grade Silicone

REF Number	Description	No. of ETO Sterile Pairs/Packages
E3147S	Internal Nasal Airway Splint - Standard	10 Pairs
E3147L	Internal Nasal Airway Splint - Large	10 Pairs



Internal Nasal Pre-Cut Splint

Material : Medical grade Silicone

REF Number	Description	No. of ETO Sterile Pairs/Packages
E3147P	Internal Nasal Pre-cut Splint	10 Pairs



Internal Bi-Valve Splints

Material : Medical grade Silicone & Fluoroplastic

REF Number	Description	Thickness (mm)	Units/Box ETO Sterile
E3148SS	Silicone-standard	0.50	10 Pairs
E3148ST	Silicone-Thick	1.00	10 Pairs
E3148FS	Fluoroplastic-Standard	0.50	10 Pairs
E3148FT	Fluoroplastic-Thick	0.25	10 Pairs
E3148FSP	Pediatric-Standard	0.50	10 Pairs
E3148FTP	Pediatric-Thin	0.25	10 Pairs




External Thermoplastic Splints

REF Number	Description	Units/Box ETO Sterile/Non Sterile
E5150S	External Thermoplastic Splint Small	10 Units
E5150S	External Thermoplastic Splint Medium	10 Units
E5150S	External Thermoplastic Splint Large	10 Units



External Aluminium Splints

REF Number	Description	Units/Box ETO Sterile/Non Sterile
E6149S	External Aluminium Splints Small	10 Units
E6149M	External Aluminium Splints Medium	10 Units
E6149L	External Aluminium Splints Large	10 Units



Silicone Sheets

Material : Medical grade Silicone

REF Number	Size (mm)	Thickness (mm)	Units/Box ETO Sterile
E3146S	55x75	3.00	10 Units
E3146T	55x75	5.00	10 Units
E3146NS	55x75	1.00	10 Units
E3146MS	55x75	0.50	10 Units
E3146US	55x75	0.25	10 Units
E3146UUS	38x50	0.13	10 Units



Silicone Blocks

Material: Medical Grade Silicone

REF Number	Size (mm)	Thickness (mm)	Units/Box ETO Sterile
E3146B5	50x30x10	10.00	1 Unit
E3146BM	60x40x15	15.00	1 Unit
E3146BL	70x50x20	20.00	1 Unit



Epistaxis Balloon Catheter

Material: Medical Grade Silicone

REF Number	Product Name	Specification	Units/Box ETO Sterile
E31030	Intra Pack	30 CC I ml	10 Units
E31010	Post pack	10 CC/ ml	10 Units
E31031	Double Balloon	30 & 10 CC/ ml	10 Units



External Splint Kits


Aluminium Splint kit	E6149K
Aluminium Splint (Small) OR	1 Unit
Aluminium Splint (Medium) OR	1 Unit
Aluminium Splint (Large)	1 Unit
Alcohol Swab	1 Unit
Skin Preparation Swab	1 Unit
Skin Closure Tape	1 Unit
Dorsum Pad	1 Unit




Thermoplastic Splint kit	E5150K
Thermoplastic Splint (Small) OR	1 Unit
Thermoplastic Splint (Medium) OR	1 Unit
Thermoplastic Splint (Large)	1 Unit
Alcohol Swab	1 Unit
Skin Preparation Swab	1 Unit
Skin Closure Tape	1 Unit
Dorsum Pad	1 Unit



Ear Wick				
REF Number	Product Name	Specification	Product Description	Units/Box Gamma Sterile
E7161	Ear Wick	Dia.9.0mmX15.0mm	No String	20 Units
E7162	Ear Wick	Dia.9.0mmX24.0mm	No String	20 Units
E7163	Ear Wick	Dia.12.0mmX15.0mm	No String	20 Units
E7164	Ear Wick	Dia.12.0mmX24.0mm	No String	20 Units
E7169	Ear Wick	Dia.9.0mmX15.0mm	Fenestrated	20 Units
E7170	Ear Wick	Dia.9.0mmX24.0mm	Fenestrated	20 Units




Ear Packs			
REF Number	Specification	Product Description	Units/Box Gamma Sterile
E7165	Dia.9.0mmX15.0mm	With String	20 Units
E7166	Dia.9.0mmX24.0mm	With String	20 Units
E7167	Dia.12.0mmX15.0mm	With String	20 Units
E7168	Dia.12.0mmX24.0mm	With String	20 Units





Anatomical Nasal Dressings			
REF Number	Specification	Product Description	Units/Box Gamma Sterile
E7188	8.0cmX1.5cmX3.0cm	With String	10 Units
E7189	8.0cmX1.5cmX3.0cm	With String & Tube	10 Units
E7190	8.0cmX1.5cmX2.0cm	Thin	10 Units



Sinus packs			
REF Number	Specification	Product Description	Units/Box Gamma Sterile
E7191	3.5cmX1.2cmX0.6cm	With String	20 Units
E7192	3.5cmX1.2cmX0.9cm	With String	20 Units





Standard Nasal Dressings			
REF Number	Specification	Product Description	Units/Box Gamma Sterile
E7176	4.5cmX1.5cmX2.0cm	No String	10 Units
E7177	4.5cmX1.5cmX2.0cm	With String	10 Units
E7178	4.5cmX1.5cmX2.0cm	With String & Tube	10 Units
E7179	4.5cmX1.0cmX2.0cm	Thin	10 Units
E7180	8.0cmX1.5cmX2.0cm	No String	20 Units
E7181	8.0cmX1.5cmX2.0cm	With String	20 Units
E7182	8.0cmX1.5cmX2.0cm	With String & Tube	10 Units
E7183	8.0cmX1.0cmX2.0cm	Thin	10 Units





Epistaxis Nasal Dressings			
REF Number	Specification	Product Description	Units/Box Gamma Sterile
E7184	10.0cmX1.5cmX2.5cm	No String	10 Units
E7185	10.0cmX1.5cmX2.5cm	With String	10 Units
E7186	10.0cmX1.5cmX2.5cm	With String & Tube	10 Units
E7187	10.0cmX1.0cmX2.5cm	Thin	10 Units
E7193	5.5cmX1.0cmX2.5cm	Pediatric	10 Units



Coatad With Proprietary Surface Treatment Technology			
REF Number	Specification	Product Description	Units/Box Gamma Sterile
C7181	8.0cmX1.5cmx2.0cm	With String	20 Units
C7185	10.0cmX1.5cmx2.5cm	With String	20 Units



Bacilk Bacteriostatic Treatment Technology			
REF Number	Specification	Product Description	Units/Box Gamma Sterile
EB7181	8.0cmX1.5cmx2.0cm	With String	20 Units
EB7185	10.0cmX1.5cmx2.5cm	With String	20 Units



Lamina Laminated Structured Surface Technology			
REF Number	Specification	Product Description	Units/Box Gamma Sterile
L7181	8.0cmX1.5cmx2.0cm	With String	20 Units
L7185	10.0cmX1.5cmx2.5cm	With String	20 Units



Microspon Antimicrobial Treatment Technology			
REF Number	Specification	Product Description	Units/Box Gamma Sterile
EA7181	8.0cmX1.5cmx2.0cm	With String	20 Units
EA7185	10.0cmX1.5cmx2.5cm	With String	20 Units

Reusable Suction Tube for Nose

Material : Stainless Steel | 10 Units/Box ETO Sterile

REF Number	Specification	Specification	Specification
EI1320120R	Dia.2.0mmX120.0mm	EI1335120L	Dia.3.5mmX120.0mm
EI1320120L	Dia.2.0mmX120.0mm	EI1335140R	Dia.3.SmmX140.0mm
EI1320140R	Dia.2.0mmX140.0mm	EI1335140L	Dia.3.SmmX140.0mm
EI1320140L	Dia.2.0mmX140.0mm	EI1335160R	Dia.3.SmmX160.0mm
EI1320160R	Dia.2.0mmX160.0mm	EI1335160L	Dia.3.SmmX160.0mm
EI1320160L	Dia.2.0mmX160.0mm	EI1340120R	Dia.4.0mmX120.0mm
EI1325120R	Dia.2.SmmX120.0mm	EI1340120L	Dia.4.0mmX120.0mm
EI1325120L	Dia.2.5mmX120.0mm	EI1340140R	Dia.4.0mmX140.0mm
EI1325140R	Dia.2.5mmX140.0mm	EI1340140L	Dia.4.0mmX140.0mm
EI1325140L	Dia.2.5mmX140.0mm	EI1340160R	Dia.4.0mmX160.0mm
EI1325160R	Dia.2.5mmX160.0mm	EI1340160L	Dia.4.0mmX160.0mm
EI1325160L	Dia.2.SmmX160.0mm	EI1350120R	Dia.5.0mmX120.0mm
EI1330120R	Dia.3.0mmX120.0mm	EI1350120L	Dia.5.0mmX120.0mm
EI1330120L	Dia.3.0mmX120.0mm	EI1350140R	Dia.5.0mmX140.0mm
EI1330140R	Dia.3.0mmX140.0mm	EI1350140L	Dia.5.0mmX140.0mm
EI1330140L	Dia.3.0mmX140.0mm	EI1350160R	Dia.5.0mmX160.0mm
EI1330160R	Dia.3.0mmX160.0mm	EI1350160L	Dia.5.0mmX160.0mm
EI1330160L	Dia.3.0mmX160.0mm	EI131RS	Nose Reusable Set
EI1335120R	Dia.3.5mmX120.0mm		

Disposable Throat Suction Tube

Material : Stainless Steel & PE

REF Number	Specification	Units/Box ETO Sterile
EI1620250	Dia.2.0mmX250.0mm	10 Units
EI1630250	Dia.3.0mmX250.0mm	10 Units

Disposable Sinus Suction Tube


Material : Stainless Steel & PE

REF Number	Specification	Units/Box ETO Sterile
EI1825150	Dia.2.SmmXIS0.0mm	10 Units
EI1830150	Dia.3.0mmXIS0.0mm	10 Units
EI1840150	Dia.4.0mmXIS0.0mm	10 Units

Ear Specula

Material : Polyethylene

REF Number	Tip Diameter	Units/Box Non-Sterile
EI0951	3mm	100 Units
EI0952	4mm	100 Units
EI0953	5mm	100 Units
EI0954	6mm	100 Units



Myringotomy Knife

Material : SS & Polyethylene

REF Number	Length	Units/Box Non-Sterile
EI0966	157mm	1 Unit




Rosen Inserter

Material : Polyethylene

REF Number	Units/Box Non-Sterile
EI0959	1 Unit

ENT Burs

10 Units/Box, Non-Sterile



Stainless Steel	Tungsten Carbide	Fine Diamond	Course Diamond	Length (mm)	Head Diameter (mm)
EI2105SS	EI2105TC	EI2105FD	EI2105CD	70.0	0.5
EI2110SS	EI2110TC	EI2110FD	EI2110CD	70.0	1.0
EI2115SS	EI2115TC	EI2115FD	EI2115CD	70.0	1.5
EI2120SS	EI2120TC	EI2120FD	EI2120CD	70.0	2.0
EI2125SS	EI2125TC	EI2125FD	EI2125CD	70.0	2.5
EI2130SS	EI2130TC	EI2130FD	EI2130CD	70.0	3.0
EI2135SS	EI2135TC	EI2135FD	EI2135CD	70.0	3.5
EI2140SS	EI2140TC	EI2140FD	EI2140CD	70.0	4.0
EI2145SS	EI2145TC	EI2145FD	EI2145CD	70.0	4.5
EI2150SS	EI2150TC	EI2150FD	EI2150CD	70.0	5.0
EI2155SS	EI2155TC			70.0	5.5
EI2160SS	EI2160TC	EI2160FD	EI2160CD	70.0	6.0
EI2165SS	EI2165TC			70.0	6.5
EI2170SS	EI2170TC		EI2170CD	70.0	7.0



Straight Sinus Blades

Material : Stainless Steel

REF Number	Size	Angle	Speed	Units/Box ETO Sterile / Non-Sterile
EI2201	4mm	0°	5000RPM	1 Unit



Curved Sinus Blades

Material : Stainless Steel

REF Number	Size	Angle	Speed	Units/Box ETO Sterile / Non-Sterile
EI2202	4mm	40°	5000RPM	1 Unit
EI2203	4mm	60°	5000RPM	1 Unit



Tonsillectomy and Adenoidectomy Blades

Material : Stainless Steel


REF Number	Size	Angle	Speed	Units/Box ETO Sterile / Non-Sterile
EI2204	4mm	40°	5000RPM	1 Unit

INSTRUMENTOS ENT



Myringotomy Kit		
EI0972 Pediatric		EI0973 Standard
Jobson Horne Probe		Jobson Horne Probe
Ear Speculas- 3,4,5 MM		Ear Speculas- 4,5,6 MM
Myringotomy Knife		Myringotomy Knife
Suction Tube - 16,18 Gauge		Suction Tube - 18,20 Gauge
Micro Alligator Ear Forcep		Micro Alligator Ear Forcep
Suction Controller		Suction Controller
Rosen Inserter		Rosen Inserter
Grommet		Grommet

ESPONJAS DE GELATINA



Absorbable Haemostatic Gelatin Sponge			
REF Number	Product	Specification	Units/Box Gamma Sterile
E91018050	Gelatin Sponge	80mm x 50mm x 1mm	20 Units
E91017050	Gelatin Sponge	70mm x 50mm x 1mm	20 Units
E91048050	Gelatin Sponge	Bilateral 80mm x 50mm x 4mm	10 Units
E91048025	Gelatin Sponge	Unilateral 80mm x 25mm x 4mm	10 Units
E91078020	Gelatin Sponge	60mm x 20mm x 7mm	10 Units
E91078015	Gelatin Sponge	Nasal 80mm x 15mm x 7mm	10 Units
E91078025	Gelatin Sponge	Gynec 80mm x 25mm x 7mm	10 Units
E91072020	Gelatin Sponge	20mm x 20mm x 7mm	10 Units
E91108050	Gelatin Sponge	80mm x 50mm x 10mm	10 Units
E91107050	Gelatin Sponge	70mm x 50mm x 10mm	10 Units
E91108040	Gelatin Sponge	Strifen 80mm x 40mm x 10mm	10 Units
E91108020	Gelatin Sponge	Strip 80mm x 20mm x 10mm	10 Units
E9115200	Gelatin Sponge	200mm x 150mm x 15mm	10 Units
E911010	Gelatin Sponge	Dental 10mm x 10mm x 10mm	32 Units
E911030	Gelatin Sponge	Dial 30mm x 30mm x 10mm	10 Units
E9110120	Gelatin Sponge	120 120mm x 100mm x 10mm	02 Units
E9110125	Gelatin Sponge	125mm x 80mm x 10mm	02 Units
E9110250	Gelatin Sponge	250 250mm x 80mm x 10mm	02 Units
E9105200	Gelatin Sponge	Film 200mm x 70mm x 0.5mm	10 Units
E910203	Gelatin Sponge	MicroSponge 3mm x 2mm x 2mm	200 Units
E910404	Gelatin Sponge	MicroSponge 4mm x 4mm x 4mm	100 Units
E918005	Gelatin Sponge	Nasal Tampon 5mm Dia x 80mm	10 Units
E918030	Gelatin Sponge	Anal Tampon 30mm Dia x 80mm	06 Units
E911208	Gelatin Sponge	Dental Tampon 8mm Dia x 12mm	32 Units
E9101	Gelatin Powder	1 Gram	06 Units



www.eonmeditech.com

Para mayor información o contactos de venta, comuníquese con:

Thiemed®
 TECNOLOGÍAS MÉDICAS

Sujeto a modificaciones técnicas.
 Copyright® EON Meditech PVT. LTD.



+56 2 2245 3428



ventas@thiemed.cl



www.thiemed.cl



Sistema de
 Gestión
 ISO 9001:2015

www.tuv.com
 ID 9000030212