

BESDATA



E-pocket

The world's first integrated ENT
endoscope with flexible and
rigid quick switch



BESDATA E-pocket A/0

E-pocket

BESDATA E-pocket: The world's first integrated ENT endoscope with screen, compatible with both flexible and rigid scopes. 4.0-inch IPS HD display, WIFI real-time transmission, 2.0–4.0mm full-size scopes for all clinical scenarios.

HD **HD Imaging**
4.0" IPS HD screen for sharp visualization

WIFI Transmission
Real-time wireless sync to monitors & computers

Full-Size Scopes
2.0–4.0mm flexible&rigid scopes for all scenarios



4.0" Touch Screen
4.0-inch touch screen, responsive operation for quick and convenient control.



120° Up/Down
Flexible 120° up/down tilt for comfortable viewing in various scenarios.




WIFI Wireless Sync Transmission

Real-time wireless transmission of images to external medical monitors / computers for convenient synchronous viewing.



Full-Size Scope Compatibility

Integrated screen ENT endoscope, compatible with 2.0–4.0mm full-size flexible and rigid scopes, covering nasal, throat examinations, foreign body removal, biopsy sampling and other routine scenarios, suitable for outpatient clinics and primary medical instruments.



Hoster Machine Specifications

Model	Display	Resolution	Memory Type	Output
VOT-1000	3.95 inch LCD	720*720(RGB)	128GB	USB, HDMI

Flexible and Rigid Scopes Specifications

Model	Product Name	Working Length	Max OD	Min ID	Field of View	Direction of View (DOV)
VRN-1010	Rigid Scope	100mm	2.0mm	/	120°	0°
VRN-1020	Rigid Scope	100mm	3.1mm	/	130°	0°
VRN-1030	Rigid Scope	175mm	3.1mm	/	130°	0°
VRN-1040	Mellable Scope	175mm	4.2mm	/	130°	0°
VRN-1050	Rigid Scope 30°	175mm	4.2mm	/	130°	30°
VRN-1060	Rigid Scope 70°	175mm	5.0mm	/	130°	70°
OTT-1010	Otoscope	20mm	4.2mm	/	130°	0°
VFL-1010	Skin Scope	/	/	/	16°	0°
VRN-101S	Disposable Flexible Scope	330mm	2.8mm	1.2mm	120°	0°
VRN-102S	Disposable Flexible Scope	330mm	3.5mm	/	120°	0°
VRN-103S	Disposable Flexible Scope	330mm	3.8mm	2.0mm	130°	0°
VRN-101R	Reusable Flexible Scope	330mm	3.0mm	/	130°	0°
VRN-102R	Reusable Flexible Scope	330mm	2.8mm	1.2mm	120°	0°
VRN-103R	Reusable Flexible Scope	330mm	4.0mm	2.0mm	120°	0°
VRN-104R	Reusable Flexible Scope	330mm	3.5mm	/	120°	0°
OTT-1020	Disposable Otoscope	30mm	/	/	130°	0°
NA	Disposable Flexible Scope	100mm	4.0mm	/	130°	0°

SHENZHEN BESDATA TECHNOLOGY CO.,LTD

503, Building #2, Xinnantian Industrial Park, No. 3179, Danzi Avenue, Laokeng Community, Longtian Street, Pingshan District, Shenzhen, China.

☎ +86-755-84416070

✉ support@besdatatech.com

🌐 www.besdatatech.com