



CATÁLOGO DE PRODUCTOS

Instrumentos y
dispositivos de
excelencia para
otorrinolaringología





TITA - PROSTHESES TYPE PARTIAL VARIO CAMPANA

(Replacement of Malleus & Incus)
Material: CP Titanium Medical Grade ASTM F 67

REF Number	Shaft Diameter	Disk Diameter	Functional Length	Total Length	Units/Box Gamma Sterile
F1101	0.30	3.00	5.00 - 1.75	5.00	1 Unit



TITA - PROSTHESES TYPE TOTAL VARIO CASK

(Replacement of Complete Ossicular chain)
Material: CP Titanium Medical Grade ASTM F 67

REF Number	Shaft Diameter	Disk Diameter	Functional Length	Total Length	Units/Box Gamma Sterile
F1102	0.30	3.00	7.00 - 3.00	7.00	1 Unit



TITAHAP PROSTHESIS TYPE PORP & TORP

Material: CP Titanium Medical Grade ASTM F 67 & HAp

REF Number	Shaft Diameter	Disk Diameter	Functional Length	Total Length	Description	Units/Box Gamma Sterile
F1207	0.30	3.00	5.00 - 1.75	5.00	TitaHAp prosthesis type Vario Partial (PORP)	1 Unit
F1208	0.30	3.00	7.00 - 3.00	7.00	TitaHAp prosthesis type Vario Total (TORP)	1 Unit



TITA - PROSTHESIS PISTON TYPE LOOP

Material : CP Titanium Medical Grade ASTM F 67 & Fluoroplastic

REF Number	Functional Length	Shaft Diameter	Units/Box Gamma Sterile
F1103	7.00 to 3.50	0.40	1 Unit



TITA - PTFE VARIO PISTON TYPE LOOP

Material : CP Titanium Medical Grade ASTM F 67 & Fluoroplastic (PTFE)

REF Number	Functional Length	Shaft Diameter	Units/Box Gamma Sterile
F4119	7.00 to 4.00	0.40	1 Unit



SUPER TITANIUM VARIO PISTON TYPE LOOP

Material : CP Titanium Medical Grade ASTM F 67 & Fluoroplastic

REF Number	Functional Length	Shaft Diameter	Units/Box Gamma Sterile
F9131	7.00 to 3.50	0.40	1 Unit



SUPER TITA - PTFE VARIO PISTON TYPE LOOP

Material : CP Titanium Medical Grade ASTM F 67 & Fluoroplastic (PTFE)

REF Number	Functional Length	Shaft Diameter	Units/Box Gamma Sterile
F9132	7.00 to 4.00	0.40	1 Unit



SUPER ELASTIC NITINOL VARIO PISTON TYPE LOOP

Material : NiTiNOL & Fluoroplastic

REF Number	Functional Length	Shaft Diameter	Units/Box Gamma Sterile
F5120	7.00 to 3.50	0.40	1 Unit



SUPER ELASTIC NITINOL - PTFE VARIO PISTON TYPE LOOP

Material : NiTiNOL & Fluoroplastic (PTFE)

REF Number	Functional Length	Shaft Diameter	Units/Box Gamma Sterile
F8130	7.00 to 4.00	0.40	1 Unit



FLUOROPLASTIC PARTIAL PROSTHESIS (PORP)

Material : Fluoroplastic

REF Number	Shaft Diameter	Disk Diameter	Functional Length	Total Length	Units/Box Gamma Sterile
F2133	1.8/1.20	3.00	5.00 - 2.00	5.00	1 Unit



FLUOROPLASTIC TOTAL PROSTHESIS (TORP)

Material : Fluoroplastic

REF Number	Shaft Diameter	Disk Diameter	Functional Length	Total Length	Units/Box Gamma Sterile
F2134	0.8	3.00	7.00 - 3.00	7.00	1 Unit



FLUOROPLASTIC OFFSET PARTIAL PROSTHESIS (PORP)

Material : Fluoroplastic

REF Number	Shaft Diameter OD/ID	Functional Length	Total Length	Units/Box ETO Sterile
F2141	1.8/1.20	5.00 - 2.00	5.00	1 Unit



FLUOROPLASTIC OFFSET TOTAL PROSTHESIS (TORP)

Material : Fluoroplastic

REF Number	Shaft Diameter	Functional Length	Total Length	Units/Box ETO Sterile
F2142	0.80	7.00 - 3.00	7.00	1 Unit



FLUOROPLASTIC UNIVERSAL MIDDLE EAR PROSTHESIS OFFSET

Material : Fluoroplastic

In case of TORP-Shaft OD 0.8mm | In case of PORP-Shaft ID 1.20mm

REF Number	Functional Length	Units/Box ETO Sterile
F2139	7.00 to 2.00	1 Unit



FLUOROPLASTIC UNIVERSAL MIDDLE EAR PROSTHESIS CENTERED

Material : Fluoroplastic

In case of TORP-Shaft OD 0.8mm | In case of PORP-Shaft ID 1.20mm


REF Number	Functional Length	Units/Box ETO Sterile
F2140	7 .00 to 2.00	1 Unit



FLUOROPLASTIC STAPES PISTON

Material : Fluoroplastic

REF Number	Shaft Diameter	Loop Inner Diameter	Functional Length	Total Length	Units/Box ETO Sterile
F2127	0.5	0.60	6.00 - 3.50	6.00	1 Unit
F2128	0.4	0.60	6.00 - 3.50	6.00	1 Unit
F2129	0.6	0.60	6.00 - 3.50	6.00	1 Unit



MIDDLE EAR IMPLANT CUTTER

Material : Stainless Steel

REF Number	Description	Units/Box
FI01	Middle Ear Implant Cutter	1 Unit, Non-Sterile



MIDDLE EAR IMPLANT SIZER

Material : Stainless Steel

REF Number	Description	Units/Box Non-Sterile
FI02	Middle Ear Implant Sizer	1 Unit



TITA - PROSTHESIS VENT TUBE TYPE COLLAR BUTTON

Material : CP Titanium Medical Grade ASTM F 67


REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box Gamma Sterile
F1104	1.40	2.55	2.55	1.20	10 Units
F1105	1.30	2.55	2.55	1.20	10 Units
F1106	1.50	2.75	2.75	1.20	10 Units



SHEPARD TITANIUM

Material : CP Titanium Medical Grade ASTM F 67


REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box Gamma Sterile
FI 151	1.14	2.40	2.40	2.20	10 Units



ARMSTRONG I ARMSTRONG BLUE I ARMSTRONG WITH TAIL

Material : Fluoroplastic


REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box Gamma Sterile
F2109	1.14	3.54	2.40	3.85	10 Units
F2109B	1.14	3.54	2.40	3.85	10 Units
F2109T	1.14	3.54	2.40	3.85	10 Units



MODIFIED ARMSTRONG I MODIFIED ARMSTRONG BLUE
MODIFIED ARMSTRONG WITH TAIL

Material : Fluoroplastic


REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box Gamma Sterile
F2110	1.14	3.46	2.50	3.04	10 Units
F2110B	1.14	3.46	2.50	3.04	10 Units
F2110T	1.14	3.46	2.50	3.04	10 Units



DONALDSON

Material : Fluoroplastic

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box ETO Sterile
F2124	1.14	2.3	2.3	2.2	10 Units
F2124B	1.14	2.3	2.3	2.2	10 Units



REUTER BOBBIN

Material : Fluoroplastic

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box ETO Sterile
F2137	1.14	2.3	2.3	1.6	10 Units
F2137B	1.14	2.3	2.3	1.6	10 Units



SHAH AND MINI SHAH

Material : Fluoroplastic


REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box ETO Sterile
F2137	1.14	3.4	2.3	2.4	10 Units
F2112B	1.14	3.4	2.3	2.4	10 Units
F2112T	1.14	3.4	2.3	2.4	10 Units
F2112W	1.14	3.4	2.3	2.4	10 Units
F2112BW	1.14	3.4	2.3	2.4	10 Units
F2112B	1.14	3.4	2.3	2.4	10 Units
F2113	0.76	2.2	1.5	1.5	10 Units
F2113B	0.76	2.2	1.5	1.5	10 Units
F2113T	0.76	2.2	1.5	1.5	10 Units
F2113W	0.76	2.2	1.5	1.5	10 Units

SHEPARD & MINI SHEPARD

Material : Fluoroplastic




REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box ETO Sterile
F2114	1.14	2.4	2.4	2.2	10 Units
F2114B	1.14	2.4	2.4	2.2	10 Units
F2114T	1.14	2.4	2.4	2.2	10 Units
F2114W	1.14	2.4	2.4	2.2	10 Units
F2114BW	1.14	2.4	2.4	2.2	10 Units
F2114BT	1.14	2.4	2.4	2.2	10 Units
F2138	1.00	2.4	2.4	2.2	10 Units
F2138W	1.00	2.4	2.4	2.2	10 Units
F2138T	1.00	2.4	2.4	2.2	10 Units
F2138B	1.00	2.4	2.4	2.2	10 Units
F2138BW	1.00	2.4	2.4	2.2	10 Units
F2138BT	1.00	2.4	2.4	2.2	10 Units



COLLAR BUTTON (ID 1.27 MM) & (ID 1.14 MM)

Material : Fluoroplastic


REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box ETO Sterile
F2123	1.27	3.0	3.0	2.15	10 Units
F2123B	1.27	3.0	3.0	2.15	10 Units
F2123T	1.27	3.0	3.0	2.15	10 Units
F2125	1.14	2.50	2.50	2.0	10 Units
F2125B	1.14	2.50	2.50	2.0	10 Units
F2125T	1.14	2.50	2.50	2.0	10 Units



BAXTER BEVELED BOBBIN

Material : Fluoroplastic

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box ETO Sterile
F2136	1.27	2.8	2.8	2.2	10 Units
F2136B	1.27	2.8	2.8	2.2	10 Units



BEVEL BOBBIN

Material : Fluoroplastic

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box ETO Sterile
F2117	1.14	2.8	2.8	1.7	10 Units
F2117B	1.14	2.8	2.8	1.7	10 Units



STRAIGHT TUBE

Material : Fluoroplastic

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box ETO Sterile
F2116	1.14	2.7	1.54	7.25	10 Units



DONALDSON SILICONE

Material : Medical Grade Silicone

REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box ETO Sterile
F3111	1.10	2.3	2.3	2.2	10 Units

T - TUBE


Material : Medical Grade Silicone



REF Number	Inner Diameter	Inner Flange Diameter	Total Length	Units/Box ETO Sterile
F3115S	1.14	8.0	6.0	10 Units
F3115M	1.14	8.0	9.0	10 Units
F3115L	1.14	8.0	12.0	10 Units

PAPARELLA


Material : Medical Grade Silicone



REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box ETO Sterile
F3126	1.10	3.5	2.5	1.65	10 Units

GOLD VENT TUBE TYPE COLLAR BUTTON WITH WIRE

Material : Gold Plated Silver



REF Number	Inner Diameter	Inner Flange Diameter	Outer Flange Diameter	Total Length	Units/Box ETO Sterile
F1170	1.25	2.55	2.55	1.20	10 Units

MYRINGOTOMY KIT



EI0972 Pediatric	EI0973 Standard
Jobson Horne Probe	Jobson Horne Probe
Ear Speculas - 3,4,5 MM	Ear Speculas - 4,5,6 MM
Myringotomy Knife	Myringotomy Knife
Suction Tube - 16,18 Gauge	Suction Tube - 18,20 Gauge
Micro Alligator Ear Forcep	Micro Alligator Ear Forcep
Suction Controller	Suction Controller
Rosen Inserter	Rosen Inserter
Grommet	Grommet

MICRO EAR ALLIGATOR FORCECP

Material : Stainless Steel



REF Number	Jaw Length	Length	Units/Box, Non-Sterile
FI0969	3.5mm	75mm	1



EAR SPECULA

Material : Polythylene

REF Number	Tip Diameter	Units/Box, Non-Sterile
FI0951	3mm	100 Units
FI0952	4mm	100 Units
FI0953	5mm	100 Units
FI0954	6mm	100 Units



OTOSCOPE SPECULA - WELCH ALLYN TYPE I

Material : Polyethylene


Ref. Number	Tip Diameter	Units / Box Non-Sterile
FI0981WA	2.5mm	100 Units
FI0982WA	3mm	100 Units
FI0983WA	4mm	100 Units
FI0984WA	5mm	100 Units



SUCTION CONTROLLER

Material : Polyethylene

Ref. Number	Length	Units / Box Non Sterile
FI0968	39mm	10 Units



EAR CURETTES

Material : Polyethylene

Ref. Number	Description	Size	Units / Box Non-Sterile
FI0962	Jobson Home Probe-Round	3mm	100 Units
FI0963	Jobson Home Probe-Oval	2.5mm	100 Units
FI0964	Cerumen Hook	Standard	100 Units

NASAL SEPTAL BUTTON

Material : Medical grade Silicone



REF Number	Description	Size (cm)	Units/Box ETO Sterile
F3145S	Nasal Septal Button - Standard	3.0	1 Unit
F3145L	Nasal Septal Button - Large	5.0	1 Unit
F3145F	Nasal Septal Button - Eliptical	3.0	1 Unit
F3145FL	Nasal Septal Button - Eliptical Large	5.0	1 Unit



INTERNAL NASAL AIRWAY SPLINT

Material : Medical grade Silicone

REF Number	Description	No. of ETO Sterile Pairs/Packages
F3147S	Internal Nasal Airway Splint - Standard	10 Pairs
F3147L	Internal Nasal Airway Splint - Large	10 Pairs



INTERNAL NASAL PRE-CUT SPLINT

Material : Medical grade Silicone

REF Number	Description	No. of ETO Sterile Pairs/Packages
F3147P	Internal Nasal Pre-cut Splint	10 Pairs



INTERNAL BI-VALVE SPLINTS

Material : Medical grade Silicone & Fluoroplastic

REF Number	Description	Thickness (mm)	Units/Box ETO Sterile
F3148SS	Silicone-standard	0.50	10 Pairs
F3148ST	Silicone-Thick	1.00	10 Pairs
F3148FS	Fluoroplastic-Standard	0.50	10 Pairs
F3148FT	Fluoroplastic-Thick	0.25	10 Pairs
F3148FSP	Pediatric-Standard	0.50	10 Pairs
F3148FTP	Pediatric-Thin	0.25	10 Pairs



EXTERNAL THERMOPLASTIC SPLINTS

REF Number	Description	Units/Box ETO Sterile/Non Sterile
F5150S	External Thermoplastic Splint Small	10 Units
F5150S	External Thermoplastic Splint Medium	10 Units
F5150S	External Thermoplastic Splint Large	10 Units



EXTERNAL ALUMINIUM SPLINTS

REF Number	Description	Units/Box ETO Sterile/Non Sterile
F6149S	External Aluminium Splints Small	10 Units
F6149M	External Aluminium Splints Medium	10 Units
F6149L	External Aluminium Splints Large	10 Units

SILICONE SHEETS


Material : Medical grade Silicone



REF Number	Size (mm)	Thickness (mm)	Units/Box ETO Sterile
F3146S	55x75	3.00	10 Units
F3146T	55x75	5.00	10 Units
F3146NS	55x75	1.00	10 Units
F3146MS	55x75	0.50	10 Units
F3146US	55x75	0.25	10 Units
F3146UUS	38x50	0.13	10 Units

SILICONE BLOCKS

Material: Medical Grade Silicone



REF Number	Size (mm)	Thickness (mm)	Units/Box ETO Sterile
F3146B5	50x30x10	10.00	1 Unit
F3146BM	60x40x15	15.00	1 Unit
F3146BL	70x50x20	20.00	1 Unit

EPISTAXIS BALLOON CATHETER

Material: Medical Grade Silicone




REF Number	Product Name	Specification	Units/Box ETO Sterile
F31030	Intra Pack	30 CC / ml	10 Units
F31010	Post pack	10 CC / ml	10 Units
F31031	Double Balloon	30 /10 CC / ml	10 Units

INTERNAL NASAL AIRWAY BACTERIOSTATIC SPONGE SPLINT (COMBO)

Material : Medical Grade Silicone



Ref. Number	Description	No. of ETO Sterile Pairs / Packages
F3147SP	Standard	10 Pairs
F3147LP	Large	10 Pairs



EAR WICK

REF Number	Product Name	Specification	Product Description	Units/Box Gamma Sterile
F7161	Ear Wick	Dia. 9.0mmX15.0mm	No String	20 Units
F7162	Ear Wick	Dia. 9.0mmX24.0mm	No String	20 Units
F7163	Ear Pack	Dia. 12.0mmX15.0mm	No String	20 Units
F7164	Ear Pack	Dia. 12.0mmX24.0mm	No String	20 Units
F7169	Ear Wick	Dia. 9.0mmX15.0mm	Fenestrated	20 Units
F7170	Ear Wick	Dia. 9.0mmX24.0mm	Fenestrated	20 Units



EAR PACKS

REF Number	Specification	Product Description	Units/Box Gamma Sterile
F7165	Dia. 9.0mmX15.0mm	With String	20 Units
F7166	Dia. 9.0mmX24.0mm	With String	20 Units
F7167	Dia. 12.0mmX15.0mm	With String	20 Units
F7168	Dia. 12.0mmX24.0mm	With String	20 Units



ANATOMICAL NASAL DRESSINGS

REF Number	Specification	Product Description	Units/Box Gamma Sterile
F7188	8.0cmX1.5cmX3.0cm	With String	10 Units
F7189	8.0cmX1.5cmX3.0cm	With String & Tube	10 Units
F7190	8.0cmX1.0cmX2.0cm	Thin	10 Units




SINUS PACKS

REF Number	Specification	Product Description	Units/Box Gamma Sterile
F7191	3.5cmX1.2cmX0.6cm	With String	20 Units
F7192	3.5cmX1.2cmX0.9cm	With String	20 Units




STANDARD NASAL DRESSINGS

REF Number	Specification	Product Description	Units/Box Gamma Sterile
F7176	4.5cmX1.5cmX2.0cm	No String	10 Units
F7177	4.5cmX1.5cmX2.0cm	With String	10 Units
F7178	4.5cmX1.5cmX2.0cm	With String & Tube	10 Units
F7179	4.5cmX1.0cmX2.0cm	Thin	10 Units
F7180	8.0cmX1.5cmX2.0cm	No String	10 Units
F7181	8.0cmX1.5cmX2.0cm	With String	10 Units
F7182	8.0cmX1.5cmX2.0cm	With String & Tube	10 Units
F7183	8.0cmX1.0cmX2.0cm	Thin	10 Units



EPISTAXIS NASAL DRESSINGS

REF Number	Specification	Product Description	Units/Box Gamma Sterile
E7184	10.0cmX1.5cmX2.5cm	No String	10 Units
E7185	10.0cmX1.5cmX2.5cm	With String	10 Units
E7186	10.0cmX1.5cmX2.5cm	With String & Tube	10 Units
E7187	10.0cmX1.0cmX2.5cm	Thin	10 Units
E7193	5.5cmX1.0cmX2.5cm	Pediatric	10 Units



COAT AID

With Proprietary Surface Treatment Technology

REF Number	Specification	Product Description	Units/Box Gamma Sterile
F7181-C	8.0cmX1.5cmX2.0cm	With String	20 Units
F7185-C	10.0cmX1.5cmX2.5cm	With String	20 Units



BACILK

Bacteriostatic Treatment Technology

REF Number	Specification	Product Description	Units/Box Gamma Sterile
F7181-B	8.0cmX1.5cmX2.0cm	With String	20 Units
F7185-B	10.0cmX1.5cmX2.5cm	With String	20 Units



LAMINA

Laminated Structured Surface Technology

REF Number	Specification	Product Description	Units/Box Gamma Sterile
F7181-L	8.0cmX1.5cmX2.0cm	With String	20 Units
F7185-L	10.0cmX1.5cmX2.5cm	With String	20 Units



MICROSPON

Antimicrobial Treatment Technology

REF Number	Specification	Product Description	Units/Box Gamma Sterile
F7181-A	8.0cmX1.5cmX2.0cm	With String	20 Units
F7185-A	10.0cmX1.5cmX2.5cm	With String	20 Units

ABSORBABLE HAEMOSTATIC GELATIN SPONGE



REF Number	Product	Specification	Units/Box Gamma Sterile
E91018050	Gelatin Sponge	80mm x 50mm x 1mm	20 Units
E91017050	Gelatin Sponge	70mm x 50mm x 1mm	20 Units
E91048050	Gelatin Sponge	Bilateral 80mm x 50mm x 4mm	10 Units
E91048025	Gelatin Sponge	Unilateral 80mm x 25mm x 4mm	10 Units
E91078020	Gelatin Sponge	60mm x 20mm x 7mm	10 Units
E91078015	Gelatin Sponge	Nasal 80mm x 15mm x 7mm	10 Units
E91078025	Gelatin Sponge	Gynec 80mm x 25mm x 7mm	10 Units
E91072020	Gelatin Sponge	20mm x 20mm x 7mm	10 Units
E91108050	Gelatin Sponge	80mm x 50mm x 10mm	10 Units
E91107050	Gelatin Sponge	70mm x 50mm x 10mm	10 Units
E91108040	Gelatin Sponge	Strifen 80mm x 40mm x 10mm	10 Units
E91108020	Gelatin Sponge	Strip 80mm x 20mm x 10mm	10 Units
E9115200	Gelatin Sponge	200mm x 150mm x 15mm	10 Units
E911010	Gelatin Sponge	Dental 10mm x 10mm x 10mm	32 Units
E911030	Gelatin Sponge	Dial 30mm x 30mm x 10mm	10 Units
E9110120	Gelatin Sponge	120 120mm x 100mm x 10mm	02 Units
E9110125	Gelatin Sponge	125mm x 80mm x 10mm	02 Units
E9110250	Gelatin Sponge	250 250mm x 80mm x 10mm	02 Units
E9105200	Gelatin Sponge	Film 200mm x 70mm x 0.5mm	10 Units
E910203	Gelatin Sponge	MicroSponge 3mm x 2mm x 2mm	200 Units
E910404	Gelatin Sponge	MicroSponge 4mm x 4mm x 4mm	100 Units
E918005	Gelatin Sponge	Nasal Tampon 5mm Dia x 80mm	10 Units
E918030	Gelatin Sponge	Anal Tampon 30mm Dia x 80mm	06 Units
E911208	Gelatin Sponge	Dental Tampon 8mm Dia x 12mm	32 Units
E9101	Gelatin Powder	1 Gram	06 Units



DISPOSABLE THROAT SUCTION TUBE

Material : Stainless Steel & PE


REF Number	Specification	Units/Box ETO Sterile
FI1620250	Dia. 2.0mmX250.0mm	10 Units
FI1630250	Dia. 3.0mmX 250.0mm	10 Units



DISPOSABLE SINUS SUCTION TUBE

Material : Stainless Steel & PE

REF Number	Specification	Units/Box ETO Sterile
FI1825150	Dia. 2.5mmX150.0mm	10 Units
FI1830150	Dia. 3.0mmX150.0mm	10 Units
FI1840150	Dia. 4.0mmX150.0mm	10 Units



MYRINGOTOMY KNIFE

Material : SS & Polyethylene

REF Number	Length	Units/Box Non-Sterile
FI0966	157mm	1 Unit



ROSEN INSERTER

Material : Polyethylene

REF Number	Units/Box Non-Sterile
FI0959	1 Unit



STRAIGHT SINUS BLADES

Material : Stainless Steel

REF Number	Size	Angle	Speed	Units/Box ETO Sterile / Non-Sterile
FI2201	4mm	0°	5000RPM	1 Unit



CURVED SINUS BLADES

Material : Stainless Steel

REF Number	Size	Angle	Speed	Units/Box ETO Sterile / Non-Sterile
FI2202	4mm	40°	5000RPM	1 Unit
FI2203	4mm	60°	5000RPM	1 Unit
FI2205	4mm	12°	5000RPM	1 Unit



TONSILLECTOMY AND ADENOIDECTOMY BLADES

Material : Stainless Steel

REF Number	Size	Angle	Speed	Units/Box ETO Sterile / Non-Sterile
FI2204	4mm	40°	5000RPM	1 Unit



ENDROSCRUB

Material : SS & ABS

Ref. Number	Size	Angle	Units / Box ETO-Sterile
FI2901	4mm	0°	5 Unit
FI2902	4mm	30°	5 Unit

DISPOSABLE SUCTION TUBE FOR EAR

Material : Stainless Steel & PE (10 Units / Box ETO Sterile)



Ref. Number	Specification (mm)	Ref. Number	Specification (mm)
FI152580	Dia. 2.5 X 80.0	FI151080	Dia. 1.0 X 80.0
FI153080	Dia. 3.0 X 80.0	FI1510100	Dia. 1.0 X 100.0
FI152060	Dia. 2.0 X 60.0	FI150870	Dia. 0.8 x 70.0
FI152070	Dia. 2.0 x 70.0	FI1508100	Dia. 0.8 x 100.0
FI152080	Dia. 2.0 X 80.0	FI150760	Dia. 0.7 X 60.0
FI1520100	Dia. 2.0 X 100.0	FI150770	Dia. 0.7 X 70.0
FI1516100	Dia. 1.6 X 100.0	FI150780	Dia. 0.7 x 80.0
FI151460	Dia. 1.4 X 60.0	FI1507100	Dia. 0.7 X 100.0
FI151470	Dia. 1.4 X 70.0	FI150670	Dia. 0.6 X 70.0
FI151480	Dia. 1.4 X 80.0	FI150460	Dia. 0.4 X 60.0
FI1514100	Dia. 1.4 X 100.0	FI150470	Dia. 0.4 X 70.0
FI151270	Dia. 1.2 X 70.0	FI150480	Dia. 0.4 X 80.0
FI1512100	Dia. 1.2 X 100.0	FI1504100	Dia. 0.4 X 100.0
FI151060	Dia. 1.0 X 60.0	F10120	Ear Disposable Set
FI151070	Dia. 1.0 x 70.0		

DISPOSABLE SUCTION TUBE FOR NOSE

Material : Stainless Steel & PMMA (1 O Units / Box ETO Sterile)



Ref. Number	Specification (mm)	Ref. Number	Specification (mm)
FI1220100	Dia. 2.0 x 100.0	FI1235140	Dia. 3.5 X 140.0
FI1220120	Dia. 2.0 X 120.0	FI1235160	Dia. 3.5 X 160.0
FI1220140	Dia. 2.0 X 140.0	FI1240100	Dia. 4.0 x 100.0
FI1220150	Dia. 2.0 X 150.0	FI1240120	Dia. 4.0 X 120.0
FI1220160	Dia. 2.0 X 160.0	FI1240140	Dia. 4.0 X 140.0
FI1230100	Dia. 3.0 x 100.0	FI1240150	Dia. 4.0 X 150.0
FI1230120	Dia. 3.0 X 120.0	FI1240160	Dia. 4.0 X 160.0
FI1230140	Dia. 3.0 X 140.0	FI1250120	Dia. 5.0 X 120.0
FI1230150	Dia. 3.0 X 150.0	FI1250140	Dia. 5.0 X 140.0
FI1230160	Dia. 3.0 X 160.0	FI1250160	Dia. 5.0 X 160.0
FI1235120	Dia. 3.5 x 120.0	FI121DS	Nose Disposable Set

LARYNGECTOMY TUBE


Material : Medical Grade Silicone (1 Unit / Box ETO Sterile)



Ref. Number Standard	Ref. Number Fenestrated	Size	Inner Diameter (A)mm	Outer Diameter (B) mm	Length Diameter (L) mm
F3110836S	F3110836F	8/36	9.5	12	36
F3110855S	F3110855F	8/55	9.5	12	55
F3110936S	F3110936F	9/36	10.5	13.5	36
F3110955S	F3110955F	9/55	10.5	13.5	55
F3111036S	F3111036F	10/36	12	15	36
F3111055S	F3111055F	10/55	12	15	55
F3111236S	F3111236F	12/36	13.5	17	36
F3111255S	F3111255F	12/55	13.5	17	55

LARYNGECTOMY BUTTON


Material : Medical Grade Silicone (1 Unit / Box, ETO Sterile / Non-Sterile)



Ref. Number	Size	Outer Diameter (A)mm	Length (B) mm
F3131208	12/8	12	8
F3131408	14/8	14	8
F3131608	16/8	16	8
F3131808	18/8	18	8
F3131218	12/18	12	18
F3131418	14/18	14	18
F3131618	16/18	16	18
F3131818	18/18	18	18

LARYNGEAL UMBRELLA KEEL

Material : Medical Grade Silicone (1 Unit / Box, ETO Sterile)



Ref. Number Clear	Ref. Number White	Size (A) mm	Size (B) mm
F31212C	F31212W	12	16
F31214C	F31214W	14	16
F31216C	F31216W	16	16

HEAT / MOISTURE EXCHANGER (HME)



Ref. Number	Description	Units / Box ETO Sterile / Non-Sterile
FA91001	Air Heat / Moisture Exchanger (HME)	5 Units
FA91002	Moist Heat / Moisture Exchanger (HME)	5 Units
FA91006	Slim HME	30 Units

TOUCH FREE HEAT / MOISTURE EXCHANGER (HME)



Ref. Number	Description
FA910011	TouchFree HME Series - Standard (Valve Cap + 2 Cassettes Air & Moist)
FA910012	TouchFree HME Series - Medium (Valve Cap + 2 Cassettes Air & Moist)
FA910013	TouchFree HME Series - Heavy (Valve Cap + 2 Cassettes Air & Moist)
FA910014	TouchFree HME Series - Light (Valve Cap + 2 Cassettes Air & Moist)
FA910015	TouchFree HME Cassettes - Air (Without Valve Cap)
FA910016	TouchFree HME Cassettes - Moist (Without Valve Cap)
FA910017	TouchFree HME Series - Set (4 Different Valve Cap + 2 Cassettes Air & Moist)




SALIVARY BY-PASS TUBE

Material : Medical Grade Silicone (1 Unit / Box ETO Sterile)




Ref. Number Clear	Ref. Number White	Size (A) mm	Size (B) mm	Length (C) mm
F31410C	F31410W	10	24	188
F31412C	F31412W	12	24	188
F31414C	F31414W	14	27	188



STANDARD TRACHEAL T-TUBE


Material : Medical Grade Silicone (1 Unit / Box ETO Sterile)



Ref. Number Clear	Ref. Number White	Outer Dia. mm	(A) mm	(B) mm	(C) mm	(D) mm
F31506CS	F31506WS	6	10	6	48	5
F31508CS	F31508WS	8	16	12	36	6
F31510CS	F31510WS	10	23	17	44	8
F31512CS	F31512WS	12	31	21	55	11
F31514CS	F31514WS	14	32	23	64	11
F31516CS	F31516WS	16	31	23	63	12

LONG & EXTRA LONG TRACHEAL T-TUBE

Material: Medical Grade Silicone (1 Unit / Box, ETO Sterile)



Ref. No. Clear	Ref. No. White	Outer Oia. mm	(A)mm	(B) mm	(C) mm	(O) mm	Size
F31506CL	F31506WL	6	34	7	48	5	Long
F31508CL	F31508WL	8	58	12	39	7	Long
F31510CL	F31510WL	10	63	17	50	8	Long
F31512CL	F31512WL	12	72	20	59	11	Long
F31514CL	F31514WL	14	72	23	66	11	Long
F31516CL	F31516WL	16	81	24	70	12	Long
F31506CE	F31506WE	6	40	30	48	5	Extra Long
F31508CE	F31508WE	8	60	49	40	7	Extra Long
F31510CE	F31510WE	10	68	48	41	8	Extra Long
F31512CE	F31512WE	12	98	39	64	11	Extra Long
F31514CE	F31514WE	14	97	48	74	11	Extra Long




SPEAKVALVE

Ref. Number	Product Name	Units / Box Non-Sterile
FA91021	Speakvalve	1 Unit

LARYNGOLOGY ACCESSORIES





Ref. Number	Product Name	Units Box Sterile / Non-Sterile
FA91051	Shower Protect	1 Unit
FA91004	Laryngeal Tube Brush	1 Unit
FA91005	Laryngeal Tube Holder	1 Unit
FA91071	Tube Cleaning Brush 8mm	10 Units
FA91072	Tube Cleaning Brush 12mm	10 Units
FA91031	StomaGrip Adhesives Round (Disposable)	30 Units
FA91032	StomaGrip Adhesives Round Max (Disposable)	30 Units
FA91033	StomaGrip Adhesives Oval (Disposable)	30 Units
FA91034	StomaGrip Adhesive Butterfly (Disposable)	30 Units



EON Inspiring for Splendid Healthcare®
MEDITECH
 PVT. LTD.

www.eonmeditech.com



Para mayor información o contactos de venta, comuníquese con:

Thiemed®
 TECNOLOGÍAS MÉDICAS

Sujeto a modificaciones técnicas.
 Copyright® EON Meditech PVT. LTD.



+56 2 2245 3428



ventas@thiemed.cl



www.thiemed.cl



Sistema de
 Gestión
 ISO 9001:2015

www.tuv.com
 ID 9000030212